



## Devon™ Surgical Vessel Loops

Vessel Loops are used in surgical situations requiring retraction identification or occlusion of arteries, veins, nerves, ureters or tendons. They are oval in cross-section, allowing for even pressure distribution, thus minimizing damage to tissue

Product Code	Colour	Size	Loop Dimensions, mm	Loop Length, mm	GTIN Number
<b>MINI</b>					
31145728	White	Mini	1.3 +/- 0.3 x 0.9 +/- 0.3	435 mm	10884521152120
31145710	Red	Mini	1.3 +/- 0.3 x 0.9 +/- 0.3	435 mm	10884521152113
31145702	Blue	Mini	1.3 +/- 0.3 x 0.9 +/- 0.3	435 mm	10884521152106
31145660	Yellow	Mini	1.3 +/- 0.3 x 0.9 +/- 0.3	435 mm	10884521152069
<b>MAXI</b>					
31145686	White	Maxi	1.15 +/- 0.20 x 2.4 +/- 0.3	435 mm	10884521152083
31145744	Red	Maxi	1.15 +/- 0.20 x 2.4 +/- 0.3	435 mm	10884521152137
31145678	Blue	Maxi	1.15 +/- 0.20 x 2.4 +/- 0.3	435 mm	10884521152076
31145694	Yellow	Maxi	1.15 +/- 0.20 x 2.4 +/- 0.3	435 mm	10884521152090

### Material

Platinum-cured, radiopaque biocompatible silicone

### Packaging

Two Vessel loops are attached to a tray and placed into a printed pouch measuring 5.5" x 2.375".

Each pouch contains x2 vessel loops

There are 10 pouches per box

### Shipping information *(applies to all codes)*

Case Quantity	20 Loops (10 pouches containing x 2 loop per pouch)
Case Dimension (mm)	142 x 66 x 62 mm
Weight	0.078Kgs

Shelf Life 5 years

Sterile Sterile by irradiation

Latex Free

CE Marked Class II in Accordance with Annex II