



# The Essential **Medical Products** product guide



**CardinalHealth**

*Essential to care™*

## Introduction

At Cardinal Health we provide and deliver a diverse and comprehensive range of medical products focusing on patient care while reducing costs, enhancing efficiency and improving quality. Manufacturing a quality portfolio of leading brands from Compression to Surgical Gloves our range of products provide clinicians and carers with a choice of products needed to meet the daily demand for complete healthcare solutions.

---

### Fluid Management

We offer a complete portfolio of fluid management products, all in one place.

Our Flex Advantage™ single-use systems provide hygienic, safety features and efficient disposal of liquid medical waste. Our Argyle™ yankauers and tubing offer you a wide choice of design for various needs.

### Surgical gloves

Our Protexis™ surgical gloves are designed to protect the hands that heal. That's why we build quality, clinical best practice, and trust directly into every Protexis™ Surgical Glove, to provide you with the Perfect Fit for every type of procedure.

### Compression

Cardinal Health is dedicated to working together with clinicians to reduce the risk of Venous Thromboembolism (VTE), a potentially lethal complication. Our solutions are clinically proven to reduce the incidence of Deep Vein Thrombosis (DVT) and Pulmonary Embolism (PE) in a variety of Peer Reviewed Clinical Studies without increasing the risk of bleeding.

### Electrocardiology

Most acute hospitals require a variety of medical devices for electrocardiography. Available in a multiple choice of size, material and application, the Cardinal Health solutions help HCPs to optimize work flow and cost efficiency.

### Thermometry

The Genius™ 3 Tympanic Thermometer and Base provides enhanced features and can be calibrated in the field using the Genius™ Thermometer Checker / Calibrator.

### Operating Room

From protection against needle sticks to helping prevent patient cross contamination Cardinal Health offers a complete line of products and accessories to help keep clinicians and patients safe in the operating room. Our full range of OR essentials provide you with a variety of sizes and options to meet clinical needs and ensure flexibility for care setting.

### Cardiothoracic

Offers a complete line of chest drainage products and accessories, thoracic catheters and procedural trays, including perfusion products to meet the needs of clinicians in many situations.

### Nutritional Delivery

Kangaroo™ has been a trusted name in Enteral Feeding for over 30 years. At Cardinal Health we offer The latest Kangaroo™ ePump and Joey feeding pumps with the option of continuous and intermittent feed and flush in one compact, easy-to-use device. Complimented by our range of universal feeding sets and short and long term enteral access devices.

# Contents

Cardiothoracic	4
Compression	9
Dental	14
Electrodes	16
Nutritional Insufficiency	24
Obstetrics and Perinatal	37
Operating Room	40
Safety Needles & Syringes and Blood Collection	50
Sharps Disposal	53
Surgical Gloves	60
Surgical Suction	65
Thermometry	76
Urology	78
Wound Care	86
Index	96

# Cardiothoracic



Sentinel Seal™ Chest Drainage Units

Altitude™ Chest Drainage Units

Aqua-Seal™ Chest Drainage Units

Thora-Seal™ Chest Drainage Units

Thoracic Catheters and Trocars





## Sentinel Seal™ Chest Drainage Units

### Four Chamber, Dry Chest Drainage Unit

The Sentinel Seal™ four chamber chest drainage unit with it's patient assesment manometer allows continuous monitoring of the pressure inside the patients pleural space.

Code	Description	Collection Chamber	Sales Unit
8888571562	Sentinel Seal™ Chest Drainage Unit	2500 ml	5
8888571513	Sentinel Seal™ Chest Drainage Unit with dual collection chamber	1370 ml and 1030 ml	5
8888571489	Sentinel Seal™ Chest Drainage Unit with tube clamp & connector	2500 ml	5



## Altitude™ Chest Drainage Units

### Altitude™ Three Chamber, Dry Chest Drainage Unit

Single-use, thoracic drainage dry system with three chambers for use under gravity or suction with mechanical suction regulator.

Code	Description	Collection Chamber	Sales Unit
8888571370	Altitude™ Chest Drainage Unit	2500 ml	5
8888571371	Altitude™ Chest Drainage Unit with dual collection chamber	1370 ml and 1030 ml	5



## Aqua-Seal™ Chest Drainage Units

### Three Chamber System

Single-use, thoracic drainage system with three chambers for use under gravity or suction.

Code	Description	Collection Chamber	Sales Unit
8888571299	Aqua-Seal™ Chest Drainage Unit	2300 ml	5
8888571406	Aqua-Seal™ Chest Drainage Unit with dual collection chamber	1265 ml and 1050 ml	5
8888571315	Aqua-Seal™ Wet Suction Easy Change connector	2300 ml	5



## Thora-Seal™ Chest Drainage Units

### Three Chamber System

Three chamber chest drainage unit, with water seal chamber and separate collection chamber.

Code	Description	Collection Chamber	Sales Unit
8884713308	Thora-Seal™ three chamber system	2500 ml	4
8884713900	Thora-Seal™ replacement collection chamber	2500 ml	6



### Two Chamber System

Two chamber chest drainage unit with a water seal chamber and a separate collection chamber.

Code	Description	Collection Chamber	Sales Unit
1180571570	Thora-Seal™ two chamber system	2600 ml	5



### One Chamber System

Single chamber chest drainage unit.

Code	Description	Collection Chamber	Sales Unit
8884713100	Thora-Seal™ complete system	2000 ml	6
1814713105	Thora-Seal™ bottle only	2000 ml	6

## Thoracic Catheters and Trocars

Used to drain air and fluid from the chest or abdominal cavity.



### Turkel™ Safety Thoracentesis and Paracentesis System

The Turkel™ safety system is a blunt safety cannula, housed within a sharp, beveled hollow needle.

Code	Description	Sales Unit
8888566067	Turkel™ Safety Thoracentesis and Paracentesis System 8 Fr/Ch Safety Needle and Catheter – 12.5 cm 50 ml Monoject™ Luer Lock Syringe, 2000 ml Drainage Bag	10



### Argyle™ Thoracic Catheters

Single-use, thin-walled catheters for maximum drainage capacity per outside diameter (50 cm) with radio-opaque control line, latex free.

Code	Description	Diameter Fr/Ch	N. of Bores	Sales Unit
<b>Thoracic Catheters</b>				
8888570507	Straight in PVC	12	4	10
8888570515	Straight in PVC	16	4	10
8888570523	Straight in PVC	20	5	10
8888570531	Straight in PVC	24	6	10
8888570549	Straight in PVC	28	6	10
8888570556	Straight in PVC	32	6	10
8888570564	Straight in PVC	36	6	10
8888570572	Straight in PVC	40	6	10
8888571018	Angled in PVC	16	4	10
8888571026	Angled in PVC	20	5	10
8888571034	Angled in PVC	24	6	10
8888571042	Angled in PVC	28	6	10
8888571059	Angled in PVC	32	6	10
8888571067	Angled in PVC	36	6	10
<b>Y Connector for Catheters</b>				
8884714200	"Y" extension for two thoracic catheters			10



## Argyle™ Trocar Catheters

Thermosensitive PVC catheter, pre-mounted on the trocar for immediate use with smooth elliptical holes and with x-ray line, latex free.

Code	Description	Diameter Fr/Ch	Length	Sales Unit
<b>Trocar Catheters</b>				
8888560805	Sharp tip	8	23 cm	10
8888561019	Sharp tip	10	23 cm	10
8888561027	Sharp tip	12	23 cm	10
8888561035	Sharp tip	16	25 cm	10
8888561043	Sharp tip	20	40 cm	10
8888561050	Sharp tip	24	40 cm	10
8888561068	Sharp tip	28	40 cm	10
8888561076	Sharp tip	32	40 cm	10

# Compression



Kendall SCD™ Comfort Sleeves  
& 700 Series Controller

Kendall SCD™ Express sleeves

A-V Impulse™ System

T.E.D.™ Anti-embolism Stockings

T.E.D.™ Anti-Embolism Stockings  
for Continuing Care



## Kendall SCD™ Comfort Sleeves & 700 Series Controller

### Kendall SCD™ Sequential Compression Comfort Sleeves – Thigh Length

Code	Description	Thigh Circumference	Size	Sales Unit
74010	Thigh length	Up to 40.6 cm	X-Small	5 pairs
74011	Thigh length	40.6 to 55.9 cm	Small	5 pairs
74012	Thigh length	55.9 to 71.1 cm	Medium	5 pairs
74013	Thigh length	71.1 to 91.4 cm	Large	3 pairs

### Kendall SCD™ Sequential Compression Comfort Sleeves – Knee Length

Code	Description	Calf Circumference	Size	Sales Unit
74021	Knee length	Up to 35.5 cm	Small	5 pairs
74022	Knee length	35.5 to 53.3 cm	Medium	5 pairs
74023	Knee length	53.3 to 66.0 cm	Large	5 pairs

### Kendall SCD™ Sequential Compression Tear-Away – Thigh Length

Code	Description	Thigh Circumference	Size	Sales Unit
74041	Thigh length	40.6 to 55.9 cm	Small	5 pairs
74042	Thigh length	55.9 to 71.1 cm	Medium	5 pairs
74043	Thigh length	71.1 to 91.4 cm	Large	3 pairs



### Kendall SCD™ 700 Series Controller

Code	Description	Sales Unit
295250	SCD™ 700 series controller – EU	1
295251	SCD™ 700 series controller – UK	1

## Kendall SCD™ Express Sleeves

### Kendall SCD™ Express Sleeves



Code	Description	Calf/Thigh Circumference	Size	Sales Unit
73022	Knee length	34.3 to 53.3 cm	Medium	5 pairs
73023	Knee length	53.3 to 66.0 cm	Large	5 pairs
9790	Knee length	55.9 to 81.2 cm	X-Large	5 pairs
73011	Thigh length	31.1 to 55.9 cm	Small	5 pairs
73012	Thigh length	55.9 to 71.1 cm	Medium	5 pairs
73013	Thigh length	71.1 to 91.4 cm	Large	3 pairs



## Kendall SCD™ Express Sleeves

### Kendall SCD™ Express Tear-Away Sleeves

Code	Description	Thigh Circumference	Size	Sales Unit
73041	Thigh length	31.1 to 55.9 cm	Small	5 pairs
73042	Thigh length	55.9 to 71.1 cm	Medium	5 pairs
73043	Thigh length	71.1 to 91.4 cm	Large	3 pairs

### Kendall SCD™ Express Sterile Sleeve

Code	Description	Thigh Circumference	Size	Sales Unit
9736	Thigh length	55.9 to 71.1 cm	Medium	5 singles

### Kendall SCD™ Express Foot Cuff

Code	Description	Size	Sales Unit
73032	Foot cuff	Regular	10 singles
73033	Foot cuff	Large	10 singles

### Kendall SCD™ Express Tubing

Code	Description	Sales Unit
9528	Compression system tubing set 2.1 m	1
9595	Compression system tubing set (ambulatory extension) 1.2 m	1



## A-V Impulse™ System

### A-V Impulse™ 6000 Series Controller

Code	Description	Sales Unit
AV6000-EUR	Dual channel controller	1 each
AV6000-UK	Dual Channel controller UK version	1 each

### ImPad™ Rigid Sole Foot Cover

Code	Description	Sales Unit
AV740-5	Medium	5 pairs
AV750-5	Large	5 pairs

### A-V Impulse™ Tubing Assembly

Code	Description	Sales Unit
AV810-01	Air hose blue, 3 m	1 each
AV820-01	Air hose red, 3 m	1 each
AV830-00	Air hose grey, 3 m	1 each



## T.E.D.™ Anti-embolism Stockings

### Thigh Length Style

Thigh Circumference	Calf Girth	Code	Leg Length	Size	Length Colour	Sales Unit
Less than 63.5 cm	Small	3071LF	Short	A	Yellow	6 pairs
		3130LF	Regular	B	White	6 pairs
		3222LF	Long	C	Yellow	6 pairs
	Medium	3310LF	Short	D	White	6 pairs
		3416LF	Regular	E	Yellow	6 pairs
		3549LF	Long	F	White	6 pairs
	Large	3634LF	Short	G	Yellow	6 pairs
		3728LF	Regular	H	White	6 pairs
		3856LF	Long	J	Yellow	6 pairs
63.5 to 81.3 cm	X-Large	4010LF	Short	K	White	6 pairs
		4114LF	Regular	L	Yellow	6 pairs
		4216LF	Long	M	White	6 pairs
		3180LF	Short	N	Yellow	6 pairs
		3181LF	Regular	P	White	6 pairs
		3182LF	Long	Q	Yellow	6 pairs
81.3 to 91.4 cm	XX-Large	3183LF	Short	R	White	6 pairs
		3184LF	Regular	S	Yellow	6 pairs
		3185LF	Long	T	White	6 pairs



### Thigh Length with Belt Style

Thigh Circumference	Calf Girth	Code	Leg Length	Size	Length Colour	Sales Unit	
Less than 63.5 cm	X-Small	3306	Regular	AA+	White	6 pairs	
		3320	Long	BB+	Blue	6 pairs	
	Small	3039	Regular	A+	White	6 pairs	
		3364	Long	B+	Blue	6 pairs	
	Medium	3144	Regular	C+	White	6 pairs	
		3449	Long	D+	Blue	6 pairs	
	Large	3221	Regular	E+	White	6 pairs	
		3523	Long	F+	Blue	6 pairs	
	63.5 to 81.3 cm	X-Large	3922	Regular	G+	White	6 pairs
			3995	Long	H+	Blue	6 pairs



### Knee Length Style

Thigh Circumference	Calf Girth	Code	Leg Length	Size	Length Colour	Sales Unit
Less than 30.5 cm	Small	7071	Regular	A-	White	12 pairs
		7339	Long	B-	Blue	12 pairs
30.5 to 38 cm	Medium	7115	Regular	C-	White	12 pairs
		7480	Long	D-	Blue	12 pairs
38 to 44.5 cm	Large	7203	Regular	E-	White	12 pairs
		7594	Long	F-	Blue	12 pairs
44.5 to 51 cm	X-Large	7604	Regular	G-	White	12 pairs
		7802	Long	H-	Blue	12 pairs
51 to 58.4 cm	XX-Large	7470LF	Regular	J-	White	12 pairs
		7471LF	Long	K-	Blue	12 pairs
58.4 to 66 cm	XXX-Large	7472LF	Regular	L-	White	12 pairs
		7473LF	Long	M-	Blue	12 pairs



## T.E.D.™ Anti-embolism Stockings for Continuing Care

### Knee Length Style

Calf Circumference	Length	Description	Item Code			Sales Unit
			White	Beige	Black	
<30.5 cm	<41 cm	Small Regular	4278-	4265-	4434	1 pair
	>41 cm	Small Long	4282	4317-	4572-	1 pair
30.5 to 38 cm	<43 cm	Medium Regular	4279-	4271	4435	1 pair
	>43 cm	Medium Long	4283	4323-	4573	1 pair
38 to 44.5 cm	<46 cm	Large Regular	4280	4289	4436	1 pair
	>46 cm	Large Long	4284	4335	4574	1 pair
44.5 to 51 cm	<46 cm	X-Large Regular	4281	4296-	4437-	1 pair
	>46 cm	X-Large Long	4285	4344	4575	1 pair



### Thigh Length Style

Thigh Circumference	Calf Circumference		Length			Sales Unit
			Short <74 cm	Regular 74-84 cm	Long >84 cm	
<63.5 cm	Small	30.5 cm	4297-	4300-	4303	1 pair
	Medium	30.5 to 38 cm	4298-	4301-	4304	1 pair
	Large	38 to 44.5 cm	4299	4302	4305-	1 pair

Available in white only

# Dental



Endodontic Needles



## Endodontic Needles

### Monoject™ Endodontic Needles – Sterile

Code	Needle Gauge and Length	Colour	Sales Unit
8881471232	23G x 32 mm	Orange	4 x (25 x 1 pc)
8881471273	27G x 32 mm	Yellow	4 x (25 x 1 pc)



### Monoject™ Endodontic Irrigation Syringes with Needle – Sterile

Code	Size Syringe	Needle Gauge and Length	Sales Unit
1180323114E	3 ml	23G x 32 mm (1.25 inches)	800 (8x100)
1180327114E	3 ml	27G x 32 mm (1.25 inches)	800 (8x100)

# Electrodes



ECG Electrodes for Adults

Accessories

Neonatal ECG Electrodes

Infant ECG electrodes

Electrodes for Polysomnography


Multifunction Defibrillation Electrodes

Kendall DL™ Disposable Lead-Wires



## Kendall™ ECG Electrodes for Adults

### General monitoring Electrodes - Long Term monitoring



Code	Match Code	Description	Sales Unit
31.1946.21	H34SG	Radiolucent, suitable for X-ray, CT and MRI with solid gel, strong adhesive (72 hour wearable, shower and wash resistant), with integrated abrader for skin preparation, 50 x 45 mm Ø, with stud adaptor	30/300
31.1925.21	H92SG	With solid gel for long-term use, 57 x 34 mm, with stud adaptor	50/500
31.1355.21	H135SG	With solid gel for long-term use, also suitable for neurology and paediatric applications, 43 x 35 mm Ø, stud adaptor	50/500
31.1355821	H135TSG	As H135SG, but radiolucent, suitable for CT and MRI	50/500
31.1245.21	H124SG	With special hydrogel for long-term use, also suitable for neurology and paediatric applications, 30 x 24 mm Ø, stud adaptor	50/500
31.1917.21	H91SSG	As H92SG, but equipped with strong adhesive (recommended for long-term ECG monitoring and stress test), 57 x 34 mm, with stud adaptor	50/500
31.1915821	H91TSG	As H91SSG, but radiolucent, suitable for X-ray, CT and MRI, 57 x 34 mm, with stud adaptor	50/500
31.1935.21	H93SG	As H92SG, however only 42 x 24 mm	50/500
31439725	Medi-Trace™ 133	With conductive adhesive hydrogel, packed in perforated strips of 3, 30 mm Ø with stud adaptor	3/600
31050522	Medi-Trace™ 200	With conductive adhesive hydrogel, packed in perforated strips of 5, 36,5 mm tear drop with stud adaptor	100/1000
31118733	Medi-Trace™ 100	With conductive adhesive hydrogel, Packed in perforated strips of 5, 30 mm Ø with stud adaptor	100/1000
31013926	Medi-Trace™ 530	With conductive adhesive hydrogel, pack of 30 loose electrodes, 45 mm tear drop with stud adaptor	30/600
31.1914.21	H914SG	With solid gel for long term use , 57 x 34 mm, with 4 mm banana adaptor	30/300
31.1914821	H914TSG	As H914SG but radiolucent, suitable for CT and MRI, with stud adaptor	30/300
RR3305		With conductive adhesive hydrogel, radiolucent, perforated strips of 5, 31.8x36.8 mm with stud adaptor, repositionable 2 times within 72 hours	5/600
ES40076	Medi-Trace™ H59P	Cloth backing with conductive adhesive hydrogel, repositionable, 38 x 22 mm with stud adaptor	30/450
ES40281	Medi-Trace™ H69P	Cloth backing with conductive adhesive hydrogel, repositionable, 22 mm squared with stud adaptor	30/450
22750	Medi-Trace™ 750	Clear Tape Electrodes, Conductive Adhesive Hydrogel, Radiolucent	50/1000

### General monitoring Electrodes - Short Term monitoring

Code	Match Code	Description	Sales Unit
31.1986.21	H98SG	Microporous non-woven backing, with solid gel, 55 mm Ø, with stud adaptor	30/300
31.1996.21	H99SG	Non-woven backing, with solid gel, 55 mm Ø, with stud adaptor	30/300



## Kendall™ ECG Electrodes for Adults

### Diagnostic Monitoring Electrodes

Code	Match Code	Description	Sales Unit
31.1985.21	H98LG	Microporous non-woven backing, with liquid gel, 55 mm Ø, with stud adaptor	30/300
31.1345.21	H34LG	With liquid gel, wearable for 72 hours, shower and wash resistant, particularly suitable for diagnostic and stress test, 50 x45 mm Ø, with stud adaptor	30/300
31.1663.21	H66LG	With liquid gel, 55 mm Ø, with stud adaptor	30/300
31.1995.21	H99LG	Non-woven backing, with liquid gel, 55 mm Ø, with stud adaptor	30/300



### Disposable ECG TAB Electrodes

With conductive and adhesive hydrogel

Code	Match Code	Description	Sales Unit
31433538	Kendall™ 5400	Tab Electrodes 23 x 33 mm	100/4000
30807732	Kendall™ 5500	Tab Electrodes 23 x 33 mm	100/2000



## Accessories




Code	Match Code	Description	Sales Unit
22660A	A34	Universal adaptor for stud and TAB electrodes, adapts to 4 mm banana plug	10


## Kendall™ Neonatal ECG Electrodes

### Disposable ECG Electrodes – Transparent backing

With conductive and adhesive hydrogel and transparent tape backing.




Code	Match Code/ Lead Length	Colour Code	Connector Type	Shape	Size (mm)	Special Features	Sales Unit
31424780	1041PTS/60 cm	AAMI (white, black, green)	1.5 mm DIN safety plug	Round	25 Ø	Suitable for incubator, X-ray radiolucent, repositionable	3/300
31424785	1042PTS/60 cm	AAMI (white, black, green)	1.5 mm DIN safety plug	Rectangular	25 x 30	Suitable for incubator, X-ray radiolucent, repositionable	3/300



## Disposable ECG Electrodes – Cloth backing

With conductive and adhesive hydrogel and cloth backing.

Code	Match Code/ Lead Length	Colour Code	Connector Type	Shape	Size (mm)	Special Features	Sales Unit
31424743	1051NPSM/60 cm	AAMI (white, black, green)	1.5 mm, DIN safety plug	Round	25 Ø	Suitable for incubator, repositionable, radiolucent	3/300
31424735	1052NPSM/60 cm	AAMI (white, black, green)	1.5 mm, DIN safety plug	Round	28 Ø	Suitable for incubator, repositionable	3/300
31424719	1053NPSM/60 cm	AAMI (white, black, green)	1.5 mm, DIN safety plug	Round	28 Ø	Suitable for incubator, repositionable, radiolucent	3/300
31424768	1050NPSM/60 cm	AAMI (white, black, green)	1.5 mm, DIN safety plug	Round	0 Ø	Suitable for incubator, 3 /300 repositionable	



## Disposable ECG Electrodes – Foam backing

With conductive and adhesive hydrogel and foam backing.

Code	Match Code/ Lead Length	Colour Code	Connector Type	Shape	Size (mm)	Special Features	Sales Unit
30710001	4103/45.7 cm	AAMI (white, black, green)	1.5 mm, DIN safety plug	Round	23 Ø	Suitable for incubator, repositionable	3/300
31368072	4203/45.7 cm	IEC (red, yellow, black)	1.5 mm, DIN safety plug	Round	23 Ø	Suitable for incubator, repositionable	3/300
31086575	4105/45.7 cm	AAMI (white, black, green)	1.5 mm, DIN safety plug	Oval	20 x 30	Suitable for incubator, repositionable	3/300


## Kendall™ Neonatal ECG Electrodes

### Disposable ECG Electrodes – non-woven tape backing

With conductive and adhesive hydrogel and non-woven backing.



Code	Match Code/ Lead Length	Colour Code	Connector Type	Shape	Size (mm)	Special Features	Sales Unit
31.5872.04	H87PG/F/60 cm	IEC (red, yellow, black)	1.5 mm DIN safety plug	Round	30 Ø	Suitable for incubator, repositionable	3/120
31.5072.03	H207PG/F/60 cm	IEC (red, yellow, black)	1.5 mm DIN safety plug	Round	22 x 30	Suitable for incubator, repositionable	3/120
31.5032.26	H203PG/60 cm	White	4 mm banana socket	Oval	22 x 30	Suitable for incubator, repositionable	15/120
31.5033.26	H203PGK/60 cm	White	4 mm banana, socket, extra short lead	Oval	22 x 30	Suitable for incubator, repositionable	15/120
31.5072824	H207PGT/90 cm	Carbon leads, transparent leads jacketing	1.5 mm DIN safety plug	Oval	22 x 30	Suitable for incubator, radiolucent, repositionable	3/90
31.5039124	H203PGT/90 cm	Carbon leads, transparent leads jacketing	4 mm banana socket	Oval	22 x 30	Suitable for incubator, radiolucent, repositionable	3/90
31.3876.04	H87H/F/60 cm	IEC (red, yellow, black)	1.5 mm DIN safety plug	Round	30 Ø	Suitable for standard care, repositionable	3/120
31.3076.26	H207H/60 cm	White	1.5 mm DIN safety plug	Oval	22 x 30	Suitable for standard care, repositionable	15/120
31.3276.03	H27HF/60 cm	IEC (red, yellow, black)	1.5 mm DIN safety plug	Round	24 Ø	Suitable for standard care, repositionable	3/120



## Kendall™ Infant ECG Electrodes

### Infant ECG Electrode for Prematures & Speciality Application

With conductive and adhesive hydrogel and cloth backing.

Code	Match Code/ Lead Length	Colour Code	Connector Type	Shape	Size (mm)	Special Features	Sales Unit
EP30016	LB36-LB/45.7 cm	AAMI (white, black, green)	1.5 mm, DIN safety plug	Rectangular	17 x 22	Suitable for incubator, with integrated fastening band, repositionable	3/300



## Electrodes for Polysomnography

Suitable for adult and children the Cardinal Health PSG electrodes meet the specific requirements of sleep lab polysomnography applications.

Code	Match Code	Description	Sales Unit
31.6071.26	H207LG/F60	Electrode for PSG applications with integrated 60 cm lead wire red, yellow, black	12/180
31.6072.26	H207LGN	Electrode for PSG applications with integrated 150 cm lead wire white	10/140
31.6073.26	H207LG/F150	Electrode for PSG applications with integrated 150 cm lead wire red, yellow, black	12/180
31.6074.26	H207LG/F200	Electrode for PSG applications with integrated 150 cm lead wire red, yellow, black	6/90



## Multifunction Defibrillation Electrodes

### Medi-Trace™ Classic

Code	Match Code	Can be used with compatible defibrillator modes	Quantity/ Pouch	Pairs/ Case
31319281	1310P	LIFEPAK (LP 9, LP 12, LP 15, LP 20, LP 20E, LP 500, LP 1000), LIFEPAK CR Plus, LIFEPAK EXPRESS	1 Pair	10
31177721	1210H	CodeMaster XL+ (M1772A, M1722B, M1723A, M1723B), Heartstream XLT M3500B <sup>1</sup> , HeartStart MRx (M3535A <sup>1</sup> , M3536A <sup>1</sup> ), HeartStart XL M4735A <sup>1</sup> , HeartStart 4000 M5500B <sup>1</sup>	1 Pair	10
31469219	1410Z	Zoll PD1400, PD2000, R Series BLS, R Series Plus, R Series ALS, X Series, Propaq MD	1 Pair	10
40000006	1710H	Heartstream XLT M3500B <sup>2</sup> , HeartStart MRx (M3535A <sup>2</sup> , M3536A <sup>2</sup> ), HeartStart XL M4735A <sup>2</sup> , HeartStart 4000 M5500B <sup>2</sup> , HeartStart FR2+ (M3840A, M3841A, M3860A, M3861A <sup>3</sup> )	1 Pair	10

1 Requires Philips M3508A hands-free cable

2 AED self-test will fail if electrodes are pre-connected

3 Requires Philips M3507A (barrel style) hands-free cable



## Multifunction Defibrillation Electrodes

### Medi-Trace™ Cadence PC (pre-connect)

Code	Can be used with compatible defibrillator modes	Quantity/ Pouch	Pairs/ Case
22550PC	LIFEPAK (LP 9, LP 12, LP 15, LP 20, LP 20E, LP 500, LP 1000), LIFEPAK CR Plus, LIFEPAK EXPRESS	1 Pair	10
22660PC	Heartstream XLT M3500B <sup>1</sup> , HeartStart MRx (M3535A, M3536A <sup>1</sup> ), HeartStart XL M4735A, HeartStart 4000 M5500B, HeartStart FR2+ (M3840A, M3841A, M3860A, M3861A <sup>2</sup> )	1 Pair	10
22661PC	CodeMaster XL+ (M1722A, M1722B, M1723A, M1723B), Heartstream XLT M3500B <sup>3</sup> , HeartStart MRx (M3535A, M3536A <sup>3</sup> ), HeartStart XL M4735A <sup>3</sup> , HeartStart 4000 M5500B <sup>3</sup>	1 Pair	5
22770PC	Zoll PD1400, PD2000, R Series BLS, R Series Plus, R Series ALS, X Series, Propaq MD	1 Pair	10

All above codes are suitable for patients  $\geq 8$  years or  $\geq 25$  kg

1 Requires Philips M3508A hands-free cable

2 AED self-test will fail if electrodes are pre-connected

3 Requires Philips M3507A (barrel style) hands-free cable



### Medi-Trace™ Cadence RTS (radio-transparent)

Code	Can be used with compatible defibrillator modes	Quantity/ Pouch	Pairs/ Case
22550A	LIFEPAK (LP 9, LP 12, LP 15, LP 20, LP 20E, LP 500, LP 1000), LIFEPAK CR Plus, LIFEPAK EXPRESS	1 Pair	10
22550R	LIFEPAK (LP 9, LP 12, LP 15, LP 20, LP 20E, LP 500, LP 1000), LIFEPAK CR Plus, LIFEPAK EXPRESS	1 Pair	10
22550P*	LIFEPAK (LP 9, LP 12, LP 15, LP 20, LP 20E, LP 500, LP 1000)	1 Pair	5
22770R	Zoll PD1400, PD2000, R Series, R Series Plus, R Series ALS, Propaq MD	1 Pair	10
22770P*	Zoll PD1400, PD2000, R Series BLS, R Series Plus, R Series ALS, X Series, Propaq MD	1 Pair 1 Pair	5 5
22660R	Heartstream XLT M3500B <sup>1</sup> , HeartStart MRx (M3535A, M3536A), HeartStart XL M4735A <sup>1</sup> , HeartStart 4000 M5500B <sup>1</sup> , HeartStart FR2+ (M3840A, M3841A, M3860A, M3861A <sup>2</sup> )	1 Pair	10
22660P*	Heartstream XLT M3500B, HeartStart MRx (M3535A, M3536A, M4735A), HeartStart 4000 M5500B	1 Pair	10
22661R	CodeMaster XL+ (M1722A, M1722B, M1723A, M1723B), Heartstream XLT M3500B <sup>3</sup> , HeartStart MRx (M3535A, M3536A <sup>3</sup> ), HeartStart XL M4735A <sup>3</sup> , HeartStart 4000 M5500B <sup>3</sup>	1 Pair	10
22661P	CodeMaster XL+ (M1722A, M1722B, M1723A, M1723B), Heartstream XLT M3500B <sup>3</sup> , HeartStart MRx (M3535A, M3536A <sup>3</sup> ), HeartStart XL M4735A <sup>3</sup> , HeartStart 4000 M5500B <sup>3</sup>	1 Pair	5

\* These pads are designed for patients  $< 8$  years or  $< 25$  kg.

1 Requires Philips M3508A hands-free cable

2 AED self-test will fail if electrodes are pre-connected

3 Requires Philips M3507A (barrel style) hands-free cable



## Kendall™ DL – Disposable Lead-Wires

### Kendall™ DL Bedside Lead Wire System

IEC	AAMI	Description	Quantity (each)
33103EG	-	3 Lead bedside red, green, yellow	10
33103EB	-	3 Lead bedside black, red, yellow	10
33105E	33105	5 Lead bedside	10
33110E	33110	10 Lead bedside	5
33111E	33111	Chest Leads (V-leads)	10

### Kendall™ DL Telemetry Lead Sets Direct Connect



IEC	AAMI	Description	Quantity (each)
33112E	33112	For use with 5 Lead Telemetry System, manufactured by GE Medical Systems Information Technologies, Inc., under its trademark ApexPro™	10
33113E	33113	For use with 5 Lead Telemetry System manufactured by Koninklijke Philips Electronics N.V. under its trademark Veridia™	10
33114E	33114	For use with 5 Lead Telemetry Systems manufactured by Koninklijke Philips Electronics N.V. under its trademark IntelliVue™	10
-	33125	For use with 6 Lead CH Telemetry System manufactured by GE Medical Systems Information Technologies, Inc., under its trademark ApexPro™	10
-	33126	For use with 6 lead FH Telemetry Systems manufactured by GE Medical Systems Information Technologies, Inc., under its trademark ApexPro™	10
33132E		For use with 5 Lead Draeger Medical, Inc., Dual Pin Cables	10
-	33133	For use with 5 Lead Spacelabs Healthcare Telemetry Systems	10
-	33134	For use with 6 Lead Nihon Kohden Telemetry Systems	10
-	33147	5 lead non-SPO2. For use with Philips IntelliVue™ MX40 Telemetry System	100
-	33148	6 lead non-SPO2. For use with Philips IntelliVue™ MX40 Telemetry System – Direct Connect	100





IEC	AAMI	Description	Quantity (each)
-	33150	6 lead SPO2. For use with Philips IntelliVue™ MX40 Telemetry System, compatible with Nellcor SpO <sub>2</sub> Sensors – Direct Connect	100
-	33151	5 lead SPO2. For use with Philips IntelliVue™ MX40 Telemetry System, compatible with Masimo SpO <sub>2</sub> Sensors	100
-	33152	6 lead SPO2. For use with Philips IntelliVue™ MX40 Telemetry System, compatible with Masimo SpO <sub>2</sub> Sensors	100

### Kendall™ DL Dual Connect System

IEC	AAMI	Description	Quantity (each)
33135E	-	Dual connect Lead Set + Cable - 5 lead	5
33136E	-	Dual connect Lead Set only - 5 lead	25
-	33145	Dual connect Cable only	5

### Radiolucent Disposable Lead Sets

IEC	AAMI	Description	Quantity (each)
-	33136R36	Radiolucent Lead Set, 5 Lead 90 cm	10
-	33136R72	Radiolucent Lead Set, 5 Lead 180 cm	10
-	33136RV36	Radiolucent Chest Lead Set, 5 Lead 90 cm	10

# Nutritional Insufficiency



Enteral Feeding Pumps and Sets

Paediatric Feeding Tubes

Nasogastric Feeding Tubes

Short Term Access – Naso-jejunal Tube

Gastrointestinal Access – Dual Lumen Tubes

Replacement Long Term Access – Skin Level Devices

Replacement Long Term Access –  
Gastrostomy Tubes

## Complete Nutritional Delivery System

We are committed to assist in the **prevention of nutritional insufficiency** while promoting patient monitoring, ensuring enteral feeding and/or hydration.



Kangaroo™ enteral feeding pumps and sets



Kangaroo™ enteral access



## KANGAROO™ ePump Enteral Feeding Pump

Code	Description	Sales Unit
482400J	Kangaroo™ ePump including pole clamp, power adaptor and operating manual	1

### Kangaroo™ ePump Enteral Pump Feeding Sets, DEHP-Free PVC, Sterile

Feeding sets with spikes

Code	Description	Sales Unit
<b>FEED SETS ONLY</b>		
777401	ENPlus spike set	30
777402	Dual ENPlus spike set	30
777403	ENPlus 3-in-1 set	30
<b>FEED AND WATER SETS</b>		
777405	ENPlus spike & 1000ml flush set	30
777406	ENPlus 3-in-1 & 1000ml flush set	36
777407	ENPlus spike feed & vented spike flush set	30
777408	ENPlus 3-in-1 feed & 3-in-1 flush set	30

Feeding sets with integrated bags

Code	Description	Sales Unit
<b>FEED SETS ONLY</b>		
777064	500 ml set	30
777106	1000 ml set	30
777163	1600 ml set	30
<b>FEED AND WATER SETS</b>		
778104	1000 ml feed & 1000 ml flush set	30
778161	1600 ml feed & 1000 ml flush set	30
673662	1000 ml feed & 1000 ml flush set, nonsterile	30

Adaptors

Code	Description	Sales Unit
777222	ENFit™ to ENLock Oral-luer medication port adaptor	30
777444	ENFit™ to funnel feeding tube adaptor	30
1814717995	3in1 adaptor	50



ENplus spike



Dual spike



3 in 1 adaptor



Water bag



Water vented spike



Water 3 in 1 adaptor  
Green adaptor



ENFit™ medication port



ENFit™ feeding tube connector



Feed bag



ENFit™ to ENLock Oral-luer medication port adaptor



ENFit™ to funnel feeding tube adaptor

## Kangaroo™ ePump Accessory Products

Code	Description	Sales Unit
382491	Power adaptor	1
382492	Pole clamp	1
1054715*	Battery pack	1
382493	Plug attachment set	1

\* This product is not a medical device and is not CE marked.



## Kangaroo™ Joey Enteral Feeding Pump

Code	Description	Sales Unit
383400	Kangaroo™ Joey Enteral Feeding Pump including pole clamp, power adaptor and operating manual	1

### Kangaroo™ Joey Enteral Pump Feeding Sets, DEHP-Free PVC, Sterile

Feeding sets with spikes

Code	Description	Sales Unit
<b>FEED SETS ONLY</b>		
777501	ENPlus spike set	30
777502	Dual ENPlus spike set	30
777503	ENPlus 3-in-1 set	30
<b>FEED AND WATER SETS</b>		
777505	ENPlus spike & 1000ml flush set	30
777506	ENPlus 3-in-1 & 1000ml flush set	36

Feeding sets with integrated bags

Code	Description	Sales Unit
<b>Feed only sets</b>		
666064	500 ml set	30
666106	1000 ml set	30
<b>FEED AND WATER SETS</b>		
668104	1000 ml feed & 1000 ml flush set	30



ENplus spike



Dual spike



3 in 1 adaptor



Water bag



Water 3 in 1 adaptor  
Green adaptor



Water vented spike



ENFit™ medication port



ENFit™ feeding tube connector



Feed bag



ENFit™ to ENLock Oral-luer medication port adaptor



ENFit™ to funnel feeding tube adaptor

Adaptors

Code	Description	Sales Unit
777222	ENFit™ to ENLock Oral-luer medication port adaptor	30
777444	ENFit™ to funnel feeding tube adaptor	30
1814717995	3in1 adaptor	50

## Kangaroo™ Joey Enteral Pump Accessories and Backpacks



Code	Description	Sales Unit
383491	Power adaptor	1
383493	Pole clamp	1
F010506*	Battery pack	1
382493	Plug attachment set	1
770026	Mini backpack – blue	1
770028	Super mini backpack – blue	1
770034	Mini backpack – pink	1
770029	Super mini backpack – pink	1
770036	Large backpack – black	1

\* This product is not a medical device and is not CE marked.



## Enteral Nutrition Gravity Feeding Sets

Enteral Nutrition Delivery Gravity Sets, DEHP-Free PVC, Sterile

Code	Description	Sales Unit
1180262154	Paediatric 60ml feeding set	36
777901	ENPlus gravity feeding set	30





## Paediatric Feeding Tubes

Kangaroo™ Neonatal and Paediatric Feeding Tubes – Polyurethane (PUR), Sterile

Code	Fr/Ch	Length	Sales Unit
461355E	3.5	31 cm with 1 cm depth markings	10
461412E	5	51 cm with 1 cm depth markings	10
461420E	6.5	51 cm with 1 cm depth markings	10
461438E	6.5	91 cm with 1 cm depth markings	10
461701E	8	51 cm with 1 cm depth markings	10
461503E	5	91 cm with 1 cm depth markings	10
461800E	8	107 cm with 1 cm depth markings	10
461042E	10	107 cm with 1 cm depth markings	10

Argyle™ Neonatal and Paediatric Feeding Tubes – DEHP-Free Polyvinylchloride (PVC), Sterile

Code	Fr/Ch	Length	Sales Unit
461206E	3.5	31 cm	10
460802E	5	41 cm	50
461008E	5	91 cm	50
460711E	6.5	41 cm	50
460604E	8	41 cm	50
460406E	8	107 cm	50
460208E	10	107 cm	50



## KANGAROO™ Nasogastric Feeding Tubes with ENFit™ Connection


Code	Description	Length	Weight	Stylet	Sales Unit
8884720817E	8 Fr/Ch	91 cm	3 g	NO	10
8884720825E	8 Fr/Ch	91 cm	3 g	YES	10
8884720841E	8 Fr/Ch	109 cm	3 g	NO	10
8884720858E	8 Fr/Ch	109 cm	3 g	YES	10
8884710859E	8 Fr/Ch	109 cm	7 g	YES	10
8884710867E	8 Fr/Ch	109 cm	5 g	YES	10
8884710826E	8 Fr/Ch	140 cm	7 g	YES	10
8884721088E	10 Fr/Ch	109 cm	5 g	YES	10
8884721096E	10 Fr/Ch	109 cm	5 g	NO	10
8884711006E	10 Fr/Ch	109 cm	7 g	YES	10
8884721055E	10 Fr/Ch	140 cm	5 g	YES	10
8884721211E	12 Fr/Ch	91 cm	5 g	NO	10
8884721237E	12 Fr/Ch	109 cm	5 g	NO	10
8884721252E	12 Fr/Ch	109 cm	5 g	NO	10
8884711246E	12 Fr/Ch	109 cm	7 g	NO	10
8884711253E	12 Fr/Ch	109 cm	7 g	YES	10
8884721255E	12 Fr/Ch	140 cm	5 g	YES	10



## Short Term Access – Naso-jejunal Feeding Tubes

Kangaroo™ Naso-Jejunal Feeding Tube – Polyurethane (PUR), Non-Sterile

Code	Fr/Ch	Length	Weight	Guide/wire	Sales Unit
8884752505E	12	152 cm	7 g	YES	2



Kangaroo™ Naso-Jejunal Feeding and Gastric Decompression Tube – Polyurethane (PUR), Non-Sterile

Code	Fr/Ch	Length	Weight	Guide/wire	Sales Unit
8884752513E	16	Gastric tube section	-	YES	2



## Gastrointestinal Access – Dual Lumen Tubes

Salem Sump Stomach Tube – Silicone, Sterile

Code	Fr/Ch	Length	Sales Unit
8888265108	10	90 cm	10
8888265116	12	120 cm	10
8888265124	14	120 cm	10
8888265132	16	120 cm	10
8888265140	18	120 cm	10



## Gastrointestinal Access – Dual Lumen Tubes

Salem Sump Stomach Tube with Anti-Reflux Valve –  
Silicone, Sterile

Code	Fr/Ch	Length	Sales Unit
8888266221	12	120 cm	10
8888266247	14	120 cm	10
8888266262	16	120 cm	10
8888266288	18	120 cm	10



Salem Sump Gastroduodenal Tube – Polurethane (PUR), Sterile

Code	Fr/Ch	Length	Sales Unit
1180264523	14	120 cm	25
1180264531	16	120 cm	25
1180264549	18	120 cm	25
1180264556	20	120 cm	25



Salem Sump Gastroduodenal Tube – DEHP-Free Polyvinylchloride  
(PVC), Sterile

Code	Fr/Ch	Length	Sales Unit
8888264804	10	90 cm	25
8888264820	12	108 cm	25
8888264846	14	120 cm	25
8888264861	16	120 cm	25
8888264887	18	120 cm	25
1180264895	20	120 cm	25



## Gastrointestinal Access – Dual Lumen Tubes

Salem Sump Gastroduodenal Tube with Moulded Connector – DEHP-Free Polyvinylchloride (PVC), Sterile

Code	Fr/Ch	Length	Sales Unit
1180264408	10	90 cm	25
1180264416	12	120 cm	25
1180264424	14	120 cm	25
1180264432	16	120 cm	25
1180264440	18	120 cm	25
1180264457	20	120 cm	25



Salem Sump Stomach Tube with Anti-Reflux Valve – DEHP-Free Polyvinylchloride (PVC), Sterile

Code	Fr/Ch	Length	Sales Unit
266122	14	120 cm	50
266130	16	120 cm	50
266148	18	120 cm	50



Salem Sump Dual Lumen Stomach Tube With ENFit Connection

Code	Fr/Ch	Length	Sales Unit
8888268060E	6	61 cm	10
8888268086E	8	61 cm	10
8888264911E	10	91 cm	50
8888264929E	12	122 cm	50
8888264945E	14	122 cm	50
8888264960E	16	122 cm	50
8888264986E	18	122 cm	50
266197	Anti-Reflux Valve		10



### Salem Sump Anti-Reflux Valve, Sterile

Code	Sales Unit
266197	10



### Kangaroo™ P.E.G. Accessory Products, Non-Sterile

Code	Fr/Ch	Length	Sales Unit
8884751622E	16	PEG adaptor – ENFit™	20
8884752042E	20	PEG adaptor – ENFit™	20



## Replacement Long Term Access – Skin Level Devices

### Kangaroo™ Skin Level Balloon Gastrostomy Kit – Silicone, Sterile

Code	Fr/Ch	Length	Sales Unit
712080	12	0.8 cm	1
712100	12	1.0 cm	1
712120	12	1.2 cm	1
712150	12	1.5 cm	1
712170	12	1.7 cm	1
712200	12	2.0 cm	1
712230	12	2.3 cm	1
712250	12	2.5 cm	1
712270	12	2.7 cm	1
712300	12	3.0 cm	1
712350	12	3.5 cm	1
712400	12	4.0 cm	1
712450	12	4.5 cm	1
712500	12	5.0 cm	1
714080	14	0.8 cm	1
714100	14	1.0 cm	1
714120	14	1.2 cm	1
714150	14	1.5 cm	1
714170	14	1.7 cm	1
714200	14	2.0 cm	1
714230	14	2.3 cm	1
714250	14	2.5 cm	1
714270	14	2.7 cm	1
714300	14	3.0 cm	1
714350	14	3.5 cm	1
714400	14	4.0 cm	1
714450	14	4.5 cm	1
714500	14	5.0 cm	1
716080	16	0.8 cm	1
716100	16	1.0 cm	1
716120	16	1.2 cm	1
716150	16	1.5 cm	1
716170	16	1.7 cm	1
716200	16	2.0 cm	1
716230	16	2.3 cm	1
716250	16	2.5 cm	1
716270	16	2.7 cm	1
716300	16	3.0 cm	1
716350	16	3.5 cm	1
716400	16	4.0 cm	1
716450	16	4.5 cm	1
716500	16	5.0 cm	1
718080	18	0.8 cm	1
718100	18	1.0 cm	1
718120	18	1.2 cm	1
718150	18	1.5 cm	1



## Replacement Long Term Access – Skin Level Devices

### Kangaroo™ Skin Level Balloon Gastrostomy Kit – Silicone, Sterile

Code	Fr/Ch	Length	Sales Unit
718170	18	1.7 cm	1
718200	18	2.0 cm	1
718230	18	2.3 cm	1
718250	18	2.5 cm	1
718270	18	2.7 cm	1
718300	18	3.0 cm	1
718350	18	3.5 cm	1
718400	18	4.0 cm	1
718450	18	4.5 cm	1
718500	18	5.0 cm	1
720080	20	0.8 cm	1
720100	20	1.0 cm	1
720120	20	1.2 cm	1
720150	20	1.5 cm	1
720170	20	1.7 cm	1
720200	20	2.0 cm	1
720230	20	2.3 cm	1
720250	20	2.5 cm	1
720270	20	2.7 cm	1
720300	20	3.0 cm	1
720350	20	3.5 cm	1
720400	20	4.0 cm	1
720450	20	4.5 cm	1
720500	20	5.0 cm	1
724080	24	0.8 cm	1
724100	24	1.0 cm	1
724120	24	1.2 cm	1
724150	24	1.5 cm	1
724170	24	1.7 cm	1
724200	24	2.0 cm	1
724230	24	2.3 cm	1
724250	24	2.5 cm	1
724270	24	2.7 cm	1
724300	24	3.0 cm	1
724350	24	3.5 cm	1
724400	24	4.0 cm	1
724450	24	4.5 cm	1
724500	24	5.0 cm	1



### Skin Level Gastrostomy Accessory Products

Code	Fr/Ch	Length	Sales Unit
8884741813		Bolus feeding set, non-sterile	10
8884761126		Bolus feeding set, sterile	5
8884741821		Continuous feeding set, non-sterile	10
8884761100		Continuous feeding set, sterile	5



## KANGAROO™ gastrostomy feeding tube with ENFIT™ Y-Port

### KANGAROO™ gastrostomy feeding tube

#### Features and Benefits

- Manufactured from clear medical grade silicone
- Silicone internal retention balloon seals internal stoma site to reduce the risk of fluid leakage
- External retention disk with “raised feet” design facilitates air circulation
- Short rounded distal tip
- Graduated gastrostomy shaft enables easy measurement of stoma depth
- Replaceable Y-port enables simple exchange when it becomes soiled or stretched.  
Can also be purchased separately
- Medication port, compatible with ENFit™ syringes, enables administration of medicines without the need to disconnect pump giving set

Code	Description	Sales Unit
8884715122E	12 Fr/Ch diameter, 5 ml balloon	5
8884715148E	14 Fr/Ch diameter, 5 ml balloon	5
8884720163E	16 Fr/Ch diameter, 20 ml balloon	5
8884720189E	18 Fr/Ch diameter, 20 ml balloon	5
8884720205E	20 Fr/Ch diameter, 20 ml balloon	5
8884720221E	22 Fr/Ch diameter, 20 ml balloon	5
8884720247E	24 Fr/Ch diameter, 20 ml balloon	5
8884720265E	26Fr/Ch diameter, 20 ml balloon	5
8884720285E	28 Fr/Ch diameter, 20 ml balloon	5

### Accessory

Code	Description	Sales Unit
8884752042E	Replacement ENFit™ Y-port	20

# Obstetrics and Perinatal



Fetal Spiral Electrodes

Belts

Accessories for Neonatology



## Fetal Spiral Electrodes

### Kendall™ Fetal Spiral Electrode System

- One step removal of the drive and guide tube
- Flexible guide tube
- Simple connection to leg plate, no orientation required
- Leg plate for secure connections
- Attaches to the patient via a self-adhesive grounding electrode



Code	Description	Sales Unit
<b>KENDALL™ FETAL SPIRAL ELECTRODES</b>		
31479549	Single helix with safety connector	50
<b>KENDALL™ FETAL SPIRAL ELECTRODE ATTACHMENT PADS</b>		
31481396	Small attachment pad 4.8 cm	50
50000095	Large attachment pad 9.5 cm x 4.8 cm	50



## Belts

### Kendall™ Abdominal Belts

- Soft materials
- Available in pre-cut or roll form
- Elastic to retain shape and hold transducers snugly

Code	Size	Hole frequency	Finished ends	Colour	Sales Unit
<b>KENDALL™ BUTTONHOLE BELTS</b>					
31410270	6 cm x 1.2 m	3.2 cm	Yes	Pink and blue	50 pairs
40000007	6 cm x 1.2 m	3.2 cm	Yes	Pink and blue	100 pairs

## Accessories for Neonatology

### Argyle™ IV Support Boards



Code	Description	Dimension	Sales Unit
MI00880	Argyle™ I.V. Support Board, Premie	2.5 cm x 7.6 cm	50
MI00775	Argyle™ I.V. Support Board, Small	4.1 cm x 9.8 cm	50
MI00850	Argyle™ I.V. Support Board, Large	4.8 cm x 12.4 cm	50



## Argyle™ Replogle Suction Catheters

Specially designed featuring a clear dual lumen tube and clearly marked from 5 cm to 25 cm for measuring the depth of insertion.

- Clear dual lumen tube
- X-ray opaque sentinel line tubing with sentinel eye catheter
- Latex and DEHP free

Code	Fr size	Length	Sales Unit
8888256545	6	61 cm	10/case
8888256529	8	61 cm	10/case
8888256503	10	61 cm	10/case

# Operating Room



Devon™ Instrument Handling Products

Devon™ Hands-Free-Transfer Products

Devon™ Surgical Markers

Devon™ Utility Markers

Devon™ Light Handle Systems

Devon™ Needle Counters

Devon™ OR Essentials

Devon™ Positioning Products

Devon™ Vessel Loops



## Devon™ Instrument Handling Products



### Devon™ Sponge Counting System

- Sponge counting system provides clear, constant visibility for assessing blood loss and counting up to five lab sponges or ten 4x4's per bag
- Single-use product
- Latex free
- Non sterile

Code	Description	Sales Unit
23250-50C	Sponge Counter Bag 5-10 count Devon™ Clear Non-Sterile	1
23250-50W	Sponge Counter Bag 5-10 count Devon™ White Non-Sterile	1
31144119	Dispenser box, clear back	50 box/250 case
31144127	Dispenser box, opaque back	50 box/250 case
31144168	Dispenser rack for I.V. pole	1
31144176	Dispenser rack for wall	1
31144184	Dispenser rack for I.V. pole, triple	1

These codes are not medical devices and are not CE marked.



## Devon™ Hands-Free-Transfer Products

### Devon™ Magnetic Drapes

Magnetic drapes introduce an effective way to minimise percutaneous injuries and accidental sticks during instrument transfer in the OR.

- Disposable magnetic drapes attract and hold surgical instruments on the sterile field
- Flexible drapes conform to any shape with non-slip foam backing for secure positioning
- Available in a variety of styles and two sizes
- Sterile
- The economical No-Skid Pad holds instruments in place without the use of magnets
- Raised magnets provide holding power while allowing instruments to be easily grasped

Code	Description	Size	Sales Unit
31140570	Magnetic drape	25.4 cm x 40.6 cm	30
31140562	Magnetic drape	40.6 cm x 50.8 cm	20
31140596	Magnetic drape with stays	25.4 cm x 40.6 cm	30
31140547	Magnetic drape with stays	40.6 cm x 50.8 cm	20

## Devon™ Hands-Free-Transfer Products



### Devon™ Magnetic Drapes with Neutral Zone

Non-magnetic neutral zone may be used to hold or transfer sharp instruments on the sterile field.

- Sterile

Code	Description	Neutral Zone	Size	Sales Unit
31140612	Magnetic drape with stay	1	25.4 cm x 40.6 cm	30
31140539	Magnetic drape with stays	1	40.6 cm x 50.8 cm	20
31140505	Magnetic drape	2	25.4 cm x 40.6 cm	26
31140430	Magnetic drape	2	40.6 cm x 50.8 cm	16

## Devon™ Surgical Markers

Devon™ surgical markers are sterile and single-use products. They are designed to meet the needs of plastic, orthopaedic, neurological and cardiovascular surgeons.

- Regular, fine or dual tip markers
- Sterile labels included with markers as indicated
- Flexible rulers (Flexi-Ruler) available

### Devon™ Surgical Markers with Ruler Cap

Code	Description	Sales Unit
31145926	Skin marker ruler cap-reg	25 box/100 case
31145918	Skin marker ruler cap-reg-ruler/label	25 box/100 case
31145942	Skin marker ruler cap-fine	25 box/100 case
31145892	Skin marker ruler cap-fine-ruler	25 box/100 case
31145843	Skin marker ruler cap-fine-ruler/label	25 box/100 case





### Devon™ Dual Tip Surgical Markers with Ruler Cap

Code	Description	Sales Unit
31145900	Dual Tip	25 box/100 case
31145868	Dual Tip with 9 labels	25 box/100 case
31160560	Dual tip with Flexi-ruler	25 box/100 case
31145934	Dual tip with 9 labels and flexi ruler	25 box/100 case



### Devon™ Surgical Markers

Code	Description	Sales Unit
19005-01CE	Surgical marking pen traditional regular tip	25 box/100 case
19005-02CE	Surgical marking pen traditional regular tip with ruler	25 box/100 case
19005-03CE	Surgical marking pen traditional regular tip with ruler and labels	25 box/100 case
19004-01CE	Surgical marking pen traditional fine tip	25 box/100 case
19004-02CE	Surgical marking pen traditional fine tip with ruler	25 box/100 case
19004-03CE	Surgical marking pen traditional fine tip with ruler and labels	25 box/100 case



### Devon™ Skin Marker – Non-Sterile

Code	Description	Sales Unit
31146020	Devon™ Skin Marker Non-Sterile	200 /case



### Devon™ Utility Markers

Used to identify instruments, specimens and other objects within the sterile field.

- Sterile, single-use utility marker
- Nine labels included

Code	Description	Sales Unit
31145884	Black ink, labels, sterile	25 box/100 case



## Devon™ Light Handle Systems

The Devon™ disposable light handle system provides a complete solution to your specific surgical lighting needs. The disposable and sterile system eliminates costly reprocessing and simplifies ordering and inventory management.



### Devon™ Light Glove

Flexible light handle cover – convenient and economical.

Code	Description	Sales Unit
LT-F01B	Light glove	1 pack/25 box/100 case
LT-F02	Light glove	2 pack/40 box/160 case
LT-F03	Light glove	3 pack/48 box/192 case



### Devon™ Light Handle

Rigid, one piece light handle – screws securely in place into Cardinal Health light adapters.

Code	Description	Sales Unit
31140133	Light handle	1 pack/16 box/128 case
31140125	Light handle	2 pack/28 box/112 case
31140117	Light handle	3 pack/21 box/84 case



## Devon™ Operating Room, Camera Light Cover

Provides a sterile interface between the non-sterile operating room camera and sterile environment.

- Clear lens allows for visualisation of surgical field
- Durable construction
- Saves time and cost associated with having to sterilise reusable camera covers
- Individual sterile packaging
- Latex free

Code	Description	Sales Unit
15010-1CE	OR Camera Light Cover	10 (10 covers per box)



## Devon™ Light Sleeve

Semi-rigid light handle cover – solid feel.

- Single use, disposable

Code	Description	Sales Unit
31140273	Light sleeve	1 pack/70 case
31140240	Light sleeve	2 pack/120 case



## Devon™ Needle Counters

The Devon™ needle counters provide an effective solution of secure needle placement, retention and count after use in surgery. A fast, easy scalpel blade removal system is also available.

The easy-to-use needle counters have valuable features available in a variety of combinations:

- Magnetic and foam surfaces
- Halves can be detached for use in two separate locations
- Adhesive backing secures foam block to any surface
- Outer adhesive tabs secure container halves to any work surface
- Boxlocks feature allows for secure closure and disposal, yet re-opens for recounts
- Puncture resistant, closed on all sides
- 10-100 needle count options in standard configurations



### Devon™ Foam Block Needle Counters with Boxlocks

Code	Description	Count Capacity	Sales Unit
31142352	Blade removal, double black magnet	100/100	8 box/64 case
31142360	Blade removal, adhesive	20/40	12 box/96 case
31142386	Blade removal, single black magnet	20/40	12 box/96 case
31142394	Blade removal, single black magnet	10/20	12 box/96 case
31142469	Double foam	80/80	8 box/64 case
31142402	Adhesive	40/70	8 box/64 case
31142311	Single black magnet	40/70	8 box/64 case
31142501	Adhesive	20/40	12 box/96 case
31142493	Single black magnet	20/40	12 box/96 case



This range of products are not medical devices and they are not CE marked.

## Devon™ Needle Counters

### Devon™ Double Blade Removers



Code	Description	Count Capacity	Sales Unit
31181459	Single foam block with magnet	30/30	12 box/96 case
31181467	Double magnet	25/25	12 box/96 case
31313425	Single foam strips with single magnet	25/25	12 box/96 case

This range of products are not medical devices and they are not CE marked.





## Devon™ Magnetic Needle Counters

Code	Description	Count Capacity	Sales Unit
31142238	Boxlocks, blade removal, double white magnet	30/30	12 box/96 case
31142113	Boxlocks, blade removal, double black magnet	30/30	12 box/96 case
31142154	Boxlocks, double black magnet	60/60	8 box/64 case
31142147	Boxlocks, double black magnet	30/30	12 box/96 case
31142196	Double black magnet	20/20	21 box/168 case
31142139	Clear top, boxlocks, single black magnet	30/30	8 box/64 case
31142410	Clear top, single black magnet	10/10	21 box/168 case

This range of products are not medical devices and they are not CE marked.



## Devon™ Needle Counters

### Devon™ Foam Strip Needle Counters with Boxlocks


Code	Description	Count Capacity	Sales Unit
31142451	Blade removal, double foam	30/30	12 box/96 case
31142295	Blade removal, black magnet	15/30	12 box/96 case
31142303	Double foam	30/30	12 box/96 case
31142188	Black magnet	30/60	8 box/64 case
31142212	Clear top	15/15	12 box/96 case

This range of products are not medical devices and they are not CE marked.

### Devon™ Needle Counters – Non-Sterile

Code	Description	Count Capacity	Sales Unit
31155370	Needle counter, plastic box, single magnet, clear top	10	200
31155362	Needle counter, plastic box, single foam strips, clear top	15	200
31155487	Needle counter, plastic box, double magnet	20	200
31155263	BladeGuard, plastic box, double white magnet	30	200
31155552	Pop'n Count plastic box, foam block and magnet, medium	10	200
31155503	Pop'n Count plastic box, foam block and adhesive, medium	10	200
31155545	Pop'n Count BladeGuard plastic box, foam block and magnet, medium	10	200
31155651	Pop'n Count BladeGuard plastic box, foam box and magnet, extra large	20	200
31155107	Pop'n Count BladeGuard plastic box, foam block and double magnet, extra large	80	100

This range of products are not medical devices and they are not CE marked.



## Devon™ OR Essentials

### Devon™ Cautery Tipolisher

Effectively cleans electrosurgical pencil tips, strong adhesive backing and X-ray detectable.

Code	Description	Sales Unit
300-50CE	Cleaner cautery tip sterile	25 box/100 case



## Devon™ Positioning Products

The Devon™ positioning products offer a reliable, efficient and economical solution for general use by healthcare professionals, intended to provide comfortable and safe positioning and support in the healthcare and surgical setting. They help to provide stability and elevation during surgeries, imaging and rehabilitation, and correct anatomical patient presentation during scans, surgical procedures, and recovery.

### Head Positioning

- Resilient, firm density foam designed to reduce pressure while maximizing support
- Whole-in-One positioners can also be used for legs, heels, knees or other extremities
- Slotted head positioner and head rest are designed to accommodate an endotracheal tube



Code	Description	Sales Unit
FP-HDCCE	Child head positioner	1 pack/24 case
FP-HDCRCE	Adult head positioner	1 pack/24 case
FP-HDSLCE	Slotted adult head positioner	1 pack/24 case
FP-HEAD9CE	Head Donut 9in head positioner	1 pack/36 case
FP-HEAD7CE	Head Donut 7in Head Positioner	1 pack/72 case
FP-HDMCE	Foam Multi Ring Head Donut	1 pack/36 case
FP-HDFRCE	Frame Head Positioner	1 pack/28 case
FP-HDSFCE	Soft touch head positioner	1 pack/28 case



## Positioning Torso, Upper and Lower Extremities

- Cushioning and protection for arms, shoulders, elbows, ulnar nerves, bony prominences and other pressure points
- Support of legs, knees, heels, ankles and feet
- Helps to reduce pressure while maintaining proper circulation
- Designed to provide adequate exposure to operative site
- Provides protection, correct alignment and comfort for a variety of procedures

Code	Description	Size	Sales Unit
FP-UN1CE	Convoluted ulnar nerve protector	41 x 15 x 5 cm	2 pack/72 case
FP-LARMCCE	Foam Laminectomy Arm Cradle	61 x 13 x 8 cm	2 pack/24 case
FP-ARTKNCE	Arthroscopic knee holder	-	1 pack/10 case
FP-ARTLECE	Well leg holder	-	1 pack/18 case
FP-HEEL1CE	Convoluted foot & heel protector	-	2 pack/72 case
FP-ORTB2CE	Convoluted OR table pad	183 x 51 x 5 cm	1 pack/12 case
FP-CHESTCE	Chest roll	43 x 13 cm	2 pack/12 case
FP-OVER1CE	Utility pad	51 x 30 x 30 cm	1 pack/12 case
FP-ARMB1CE	Convoluted arm board pad	51 x 20 x 5 cm	12 pair
FP-ABDSCE	Small abduction pillow	46 x 30 x 15 cm	6 case
FP-ABDMCE	Medium abduction pillow	56 x 38 x 15 cm	6 case
FP-ABDLCE	Large abduction pillow	63 x 46 x 15 cm	6 case



## Devon™ Positioning Products

### Body Positioners

Code	Description	Size	Sales Unit
FP-BAW2CE	Foam Body Alignment Wedge 18 inch	46 x 20 x 20 cm	8
FP-BAW1CE	Foam Body Alignment Wedge 22 inch	55 x 30 x 22 cm	8



## Devon™ Vessel Loops


Made of Silicone. Radiopaque bands that find practical application in situation where arteries veins, nerves or tendons must be looped for identification, retraction or occlusion. Available in four different colors and 2 sizes.

Code	Description	Sales Unit
2400-1RCE	Silicone Vessel Loop, Red, Mini, Radiopaque, 2/PK, Sterile	20 pack/200 per case
2400-1BCE	Silicone Vessel Loop, Blue, Mini, Radiopaque, 2/PK, Sterile	20 pack/200 per case
2400-1YCE	Silicone Vessel Loop, Yellow, Mini, Radiopaque, 2/PK, Sterile	20 pack/200 per case
2400-1WCE	Silicone Vessel Loop, White, Mini, Radiopaque, 2/PK, Sterile	20 pack/200 per case
2400-2RCE	Silicone Vessel Loop, Red, Maxi, Radiopaque, 2/PK, Sterile	20 pack/200 per case
2400-2BCE	Silicone Vessel Loop, Blue, Maxi, Radiopaque, 2/PK, Sterile	20 pack/200 per case
2400-2YCE	Silicone Vessel Loop, Yellow, Maxi, Radiopaque, 2/PK, Sterile	20 pack/200 per case
2400-2WCE	Silicone Vessel Loop, White, Maxi, Radiopaque, 2/PK, Sterile	20 pack/200 per case

# Safety Needles & Syringes and Blood Collection



Safety Needles



## Safety Needles

### Magellan™ Safety Hypodermic Needles

The Magellan™ safety hypodermic needle is designed to reduce the risk of inadvertent needlestick injuries among healthcare workers.

The device consists of a hypodermic needle with an integral safety shield that securely and permanently locks over the needle after use. With one-handed technique, the needle is easily locked using the thumb, finger or surface.

Code	Needle Gauge and Length	Sales Unit
8881850558	25G x 5/8" (16 mm)	10 x 50 pcs
8881850510	25G x 1" (25 mm)	10 x 50 pcs
8881850358	23G x 5/8" (16 mm)	10 x 50 pcs
8881850310	23G x 1" (25 mm)	10 x 50 pcs
8881850210	22G x 1" (25 mm)	10 x 50 pcs
8881850215	22G x 1-1/2" (38 mm)	10 x 50 pcs
8881850158	21G x 5/8" (16 mm)	10 x 50 pcs
8881850110	21G x 1" (25 mm)	10 x 50 pcs
8881850115	21G x 1-1/2" (38 mm)	10 x 50 pcs
8881850010	20G x 1" (25 mm)	10 x 50 pcs
8881850015	20G x 1-1/2" (38 mm)	10 x 50 pcs
8881850910	19G x 1" (25 mm)	10 x 50 pcs
8881850915	19G x 1-1/2" (38 mm)	10 x 50 pcs
8881850810	18G x 1" (25 mm)	10 x 50 pcs
8881850815	18G x 1-1/2" (38 mm)	10 x 50 pcs



### Magellan™ Insulin and Tuberculin Safety Syringe

Code	Size	Needle Gauge and Length	Colour	Cannula	Sales Unit
<b>INSULIN</b>					
8881892910	1 ml	29G x 1/2" (13 mm)	Orange	Fixed	50 box/500 case
8881892950	0.5 ml	29G x 1/2" (13 mm)	Orange	Fixed	50 box/500 case
8881892930	0.3 ml	29G x 1/2" (13 mm)	Orange	Fixed	50 box/500 case
8881893010	1 ml	30G x 5/16" (8 mm)	Orange	Fixed	50 box/500 case
8881893050	0.5 ml	30G x 5/16" (8 mm)	Orange	Fixed	50 box/500 case
8881893030	0.3 ml	30G x 5/16" (8 mm)	Orange	Fixed	50 box/500 case
<b>TUBERCULIN</b>					
8881882812	1 ml	28G x 1/2" (13 mm)	Blue	Fixed	50 box/500 case

## Safety Needles



### Monoject™ Safety Insulin Syringes – Sterile

Code	Size	Graduation	Needle Gauge and Length	Cannula	Sales Unit
8881511344	0.3 ml	U-100	30G x 11/16" (8 mm)	Fixed	5 x 100 pcs
8881511336	0.5 ml	U-100	30G x 11/16" (8 mm)	Fixed	5 x 100 pcs



### Monoject™ Hypodermic Safety Needle

The Monoject™ Hypodermic Safety needle has a hinged safety device designed to protect healthcare workers from exposure to needle-stick injuries.

- Integrated Safety Sheath
- Double-Lock, Hinged Arm System
- One-Handed Activation
- ISO Color-Coded Safety Shield

Code	Pack Colour	Needle Gauge and Length	Box/Case Qty
1181810		18G x 1" (1.27 mm x 2.5 cm)	100/800
1181815		18G x 1-1/2" (1.27 mm x 3.8 cm)	100/800
1181910		19G x 1" (1.067 mm x 2.5 cm)	100/800
1181915		19G x 1-1/2" (1.067 mm x 3.8 cm)	100/800
1182010		20G x 1" (0.902 mm x 2.5 cm)	100/800
1182015		20G x 1-1/2" (0.902 mm x 3.8 cm)	100/800
1182110		21G x 1" (0.813 mm x 2.5 cm)	100/800
1182115		21G x 1-1/2" (0.813 mm x 3.8 cm)	100/800
1182210		22G x 1" (0.711 mm x 2.5 cm)	100/800
1182215		22G x 1-1/2" (0.711 mm x 3.8 cm)	100/800
1182310		23G x 1" (0.635 mm x 2.5 cm)	100/800
1182558		25G x 5/8" (0.508 mm x 1.6 cm)	100/800
1182510		25G x 1" (0.508 mm x 2.5 cm)	100/800
1182705		27G x 1/2" (0.406 mm x 1.3 cm)	100/800
1182734		27G x 3/4" (0.406 mm x 1.9 cm)	100/800
1183005		30G x 1/2" (0.305 mm x 1.3 cm)	100/800



### Monoject™ Bluntfill Needle

Code	Needle Gauge and Length	Box/Case Qty
11811022F	Bluntfill with filter 18G x 1-1/2" (1.27 mm x 3.8 cm)	100/1000



# Sharps Disposal



Phlebotomy and Multi-Purpose Containers

Multi-Purpose Containers

In-Room System

Large Volume Containers

Chemotherapy Containers

PGII Large Volume Sharps Containers

Transportable Sharps Containers

Brackets and Holders

Wall Enclosures

Foot Pedal Carts and Floor Brackets

## Phlebotomy and Multi-Purpose Containers

### SharpSafety Phlebotomy Sharps Containers Key

● Red ● Yellow ● Transparent Red ● Transparent Yellow ○ Transparent ○ Moulded with recycled medical waste plastic binCode



Code	Description	Litres Capacity	Dimensions HxDxW	Brackets & Holders	Colour	Sales Unit
8900SA	SharpSafety phlebotomy sharps container adapts to present blood drawing trays and technique	1	15.9 x 10.8 x 11.4 cm	1	●	100
8906				11	●	100
1316676202	Small Monoject™ sharps container, with clear top for easier assessment of fill level	1.5	13.3 x 20 x 13.3 cm	-	●	40
1316676228				-	●	40
1522SA	SharpSafety phlebotomy sharps container with rotor opening	2	18.4 x 12 x 16.5 cm	2	●	60

### Brackets & Holders




This range of products are not medical devices and they are not CE marked.

## Monoject™ Multi-Purpose Sharps Containers

Code	Description	Litres Capacity	Dimensions HxDxW	Brackets & Holders	Colour	Sales Unit
1100676376	Monoject™ sharps container with chimney drop top lid. Clear top on translucent container for easier assessment of fill level. Can be decontaminated in an autoclave.	4	13.9 x 17.3 x 28.6 cm	11	●	40
1100676517					●	40
1100676384	Monoject™ sharps container with chimney drop top lid. Clear top on translucent container for easier assessment of fill level. Can be decontaminated in an autoclave.	8	23.7 x 17.3 x 28.6 cm	11	●	20
1100676533					●	20
1100676392	Monoject™ sharps container with chimney drop top lid. Clear top on translucent container for easier assessment of fill level. Can be decontaminated in an autoclave.	11	40.8 x 17.3 x 28.6	11	●	10
1100676558					●	10



### Brackets & Holders



This range of products are not medical devices and they are not CE marked.

## Multi-Purpose Containers

### SharpSafety Multi-Purpose Sharps Containers



Code	Description	Litres Capacity	Dimensions HxDxW	Brackets & Holders	Colour	Sales Unit
8970 8970Y	Multi-purpose container for all sharps waste with rotor opening	7.4	25.4 x 18.4 x 26.6 cm	2 7 8	● ●	20
8990SA	Multi-purpose container with clear hinged lid	7.4	25.4 x 18.4 x 26.6 cm	2 7 8	●	20
8990Y	Multi-purpose container with special gasket hinged lid	7.4	25.4 x 18.4 x 26.6 cm	2 7 8	●	20
8967Y	Multi-purpose container with horizontal drop opening lid	7.4	32.4 x 18.4 x 26.6 cm	2 7 8	● ○	20
8527R	Multi-purpose container with rotor opening and hinged lid to accommodate a variety of sharps sizes and provide temporary and permanent closure lid	11	35 x 15.2 x 35 cm	4 9	● ●	10

### Brackets & Holders



Code 8974



Code 8975

This range of products are not medical devices and they are not CE marked.

## In-Room System

### SharpSafety Counterbalanced Lid Containers



Code	Description	Litres Capacity	Dimensions HxDxW	Brackets & Holders	Colour	Sales Unit
8506Y 8507SA	Sharps container with counterbalanced lid and counter-balanced door	4.6	31.7 x 14 x 27.3 cm	3 6 15	● ●	20
8534SA	Sharps container with counterbalanced lid and counter-balanced door	7.4	30 x 16.5 x 35 cm	9 16	○ ●	10
8536SA 8537SA 8537Y	Sharps container with counterbalanced lid and counter-balanced door	11	42 x 15.2 x 35 cm	4 9 16	○ ●	10

This range of products are not medical devices and they are not CE marked.

### SharpSafety Always-Open Lid Containers

Code	Description	Litres Capacity	Dimensions HxDxW	Brackets & Holders	Colour	Sales Unit
85321R	Sharps container with always-open lid provides convenient, quick and easy point-of-use sharps disposal	7.4	29.2 x 15.2 x 35 cm	9 14	●	10
85221R	Sharps container with always-open lid provides convenient, quick and easy point-of-use sharps disposal	11	41.3 x 15.2 x 35 cm	4 9 15	○ ● ○	10

This range of products are not medical devices and they are not CE marked.

## Large Volume Containers

### SharpSafety Large Volume Containers



Code	Description	Litres Capacity	Dimensions HxDxW	Brackets & Holders	Colour	Sales Unit
8980S 8935,	Large volume sharps container with sliding lid	30	45 x 28 x 39.4 cm	18, 19	●	10
8935Y 8938		45	47.6 x 32.4 x 46.3 cm	17, 18, 22	● ●	10
		67	66 x 32.4 x 46.3 cm	17, 18, 22	● ●	5
8980, 8980Y	Large volume sharps container hinged lid and separate rotor opening for smaller sharps waste	30	45 x 28 x 39.4 cm	18, 20	● ●	10
8991		45	47.6 x 32.4 x 46.3 cm	18, 21		10
		67	66 x 32.4 x 46.3 cm	18, 21		5

This range of products are not medical devices and they are not CE marked.

# Chemotherapy Containers

## SharpSafety Large Volume Chemotherapy Containers



Code	Description	Litres Capacity	Dimensions HxDxW	Brackets & Holders	Colour	Sales Unit
8982	Chemotherapy sharps container with hinged lid	7.4	25.4 x 18.4 x 26.7 cm	2, 7, 8	●	20
8985	Chemotherapy sharps container with hinged lid	30	45 x 28 x 39.4 cm	18, 20	●	10
8931		45	47.6 x 32.4 x 46.3 cm	18, 21	●	10
8989		67	66 x 32.4 x 46.3 cm	18, 21	●	5
8934	Chemotherapy sharps container with sliding lid	45	47.6 x 32.4 x 46.3 cm	17, 18, 22	●	10
8939		65	66 x 32.4 x 46.3 cm	17, 18, 22	●	5

## Brackets & Holders

2



Code 8974

7



Code 8975

17



Code 8938FP

18



Code 8992H

19



Code 8981FP

20



Code 8980FP

21



Code 8991FP

This range of products are not medical devices and they are not CE marked.

## PGII Large Volume Sharps Containers

### SharpSafety Large Volume PGII Containers



Code	Description	Litres Capacity	Dimensions HxDxW	Brackets & Holders	Colour	Sales Unit
8998PG2	PGII container large volume, with hinged and sealing gasket lid	30	45 x 28 x 39.4 cm	18, 20	●	10
		45	47.6 x 32.4 x 46.3 cm	18, 21	●	10
		67	66 x 32.4 x 46.3 cm	18, 21	●	5

This range of products are not medical devices and they are not CE marked.

### SharpSafety Chemotherapy PGII Containers



Code	Description	Litres Capacity	Dimensions HxDxW	Brackets & Holders	Colour	Sales Unit
8985PG2	Chemotherapy sharps container	30	45 x 28 x 39.4 cm	18, 19	●	10
8989PG2	with hinged lid	45	47.6 x 32.4 x 46.3 cm	18, 21	●	10
		67	66 x 32.4 x 46.3 cm	18, 21	●	5

This range of products are not medical devices and they are not CE marked.

## Brackets and Holders

### Brackets and Holders



Code	Description	Dimensions HxDxW	Sales Unit
8974	Holder for 7.4 litre phlebotomy, multi-purpose and ChemoSafety containers. For containers code, 8996, 8970, 8970Y, 89671, 8967Y, 8982, 8990SA, 8990Y	8.2 x 19 x 28 cm	5



8518X	Locking brackets for 2 & 4.6 litre In-Room™ containers. For containers code 850201, 850301, 851201, 85120Y	21 x 16.5 x 28.6 cm	5
-------	--	---------------------	---



8975	Wall/cart bracket for 7.4 litre multi-purpose and ChemoSafety containers. For containers code 8970, 8970Y, 89671, 8967Y, 8990Y, 8996, 8982	19 x 16.5 x 25.4 cm	5
------	--	---------------------	---



8912SA	Wall Bracket for 1 litre phlebotomy container. For containers code 8900SA, 8906	–	5
--------	---	---	---

This range of products are not medical devices and they are not CE marked.

## Wall Enclosures

### Wall Enclosures

	<b>Code</b>	<b>Description</b>	<b>Dimensions HxDxW</b>	<b>Sales Unit</b>
12	8556H	Wall enclosure for 2 & 4.6 litre sharps container with glove dispenser and incl. a 4.6 litre SharpSafety™ horizontal lid container	47 x 14 x 33 cm	1
13	8550LG	Glove dispenser for 85161H and 85301H wall enclosure	17.1 x 11.4 x 30.5 cm	10
16	8550B	Separate glove dispenser with attached wall mounting bracket	18.4 x 11.4 x 31.1 cm	5

This range of products are not medical devices and they are not CE marked.

## Foot Pedal Carts and Floor Brackets

### Foot Pedal Carts and Floor Brackets

	<b>Code</b>	<b>Description</b>	<b>Dimensions HxDxW</b>	<b>Sales Unit</b>
17	8938FP	Sharps cart foot pedal cart for 45 & 67 litre sliding lid large volume and ChemoSafety containers (container not included)	78.7 x 57.8 x 38.7 cm	1
18	8992H	Floor bracket with wheels for 30, 45 & 67 litre large volume and ChemoSafety containers (container not included)	55.9 x 32.4 x 57.1 cm	1
19	8981FP	Sharps cart foot pedal cart for 30 litre sliding lid large volume and ChemoSafety containers (container not included)	55.9 x 54.6 x 40.6 cm	1
20	8980FP	Sharps cart foot pedal cart for 30 litre hinged lid large volume and ChemoSafety containers (container not included)	55.9 x 54.6 x 40.6 cm	1
21	8991FP	Sharps cart foot pedal cart for 45 & 67 litre hinged lid large volume and ChemoSafety containers (container not included)	78.7 x 57.8 x 38.7 cm	1

This range of products are not medical devices and they are not CE marked.



# Surgical Gloves




Protexis™ PI range

Protexis™ Neoprene range

Protexis™ Latex PF range

# Protexis™ Polyisoprene, Powder-Free

## Protexis™ PI Surgical Gloves



Code	Size	Material	Color	Category	Fingertip thickness	Coatings	Sales Unit
2D72PT55X	5½	Synthetic polyisoprene (PI)	Cream	Standard	0.23 mm		50 pairs/box
2D72PT60X	6						4 boxes/case
2D72PT65X	6½						
2D72PT70X	7						
2D72PT75X	7½						
2D72PT80X	8						
2D72PT85X	8½						
2D72PT90X	9						

## Protexis™ PI Micro Surgical Gloves



Code	Size	Material	Color	Category	Fingertip thickness	Coatings	Sales Unit
2D73PM55	5½	Synthetic polyisoprene (PI)	Cream	Thin	0.20 mm		50 pairs/box
2D73PM60	6						4 boxes/case
2D73PM65	6½						
2D73PM70	7						
2D73PM75	7½						
2D73PM80	8						
2D73PM85	8½						
2D73PM90	9						

## Protexis™ PI with Neu-Thera™ Surgical Gloves



Code	Size	Material	Color	Category	Fingertip thickness	Coatings	Sales Unit
2D73TE55	5½	Synthetic polyisoprene (PI)	Cream	Standard	0.23 mm	Neu-Thera™ Emollient Coating	50 pairs/box
2D73TE60	6						4 boxes/case
2D73TE65	6½						
2D73TE70	7						
2D73TE75	7½						
2D73TE80	8						
2D73TE85	8½						
2D73TE90	9						

## Protexis™ PI Blue with Neu-Thera™ Surgical Gloves



Code	Size	Material	Color	Category	Fingertip thickness	Coatings	Sales Unit
2D73EB55	5½	Synthetic polyisoprene (PI)	Royal blue	Underglove	0.20 mm	Neu-Thera™ Emollient Coating	50 pairs/box
2D73EB60	6						4 boxes/case
2D73EB65	6½						
2D73EB70	7						
2D73EB75	7½						
2D73EB80	8						
2D73EB85	8½						
2D73EB90	9						

## Protexis™ PI Ortho Surgical Gloves



Code	Size	Material	Color	Category	Fingertip thickness	Coatings	Sales Unit
2D73ET60	6	Synthetic polyisoprene (PI)	Cream	Thick	0.30 mm		40 pairs/box 4 boxes/case
2D73ET65	6½						
2D73ET70	7						
2D73ET75	7½						
2D73ET80	8						
2D73ET85	8½						
2D73ET90	9						

## Protexis™ PI Orthopaedic Surgical Gloves



Code	Size	Material	Color	Category	Fingertip thickness	Coatings	Sales Unit
2D73HT60	6	Synthetic polyisoprene (PI)	Brown	Thick	0.34 mm	Hydrogel	40 pairs/box 4 boxes/case
2D73HT65	6½						
2D73HT70	7						
2D73HT75	7½						
2D73HT80	8						
2D73HT85	8½						
2D73HT90	9						

## Protexis™ PI Textured Surgical Gloves



Code	Size	Material	Color	Category	Fingertip thickness	Coatings	Sales Unit
2D72TG55	5½	Synthetic polyisoprene (PI)	Cream	Medium	0.27 mm	Hydrogel	50 pairs/box 4 boxes/case
2D72TG60	6						
2D72TG65	6½						
2D72TG70	7						
2D72TG75	7½						
2D72TG80	8						
2D72TG85	8½						
2D72TG90	9						

## Protexis™ Neoprene, Powder-Free

## Protexis™ Neoprene Surgical Gloves



Code	Size	Material	Color	Category	Fingertip thickness	Coatings	Sales Unit
2D73DP55	5½	Synthetic neoprene	Light brown	Standard	0.17 mm	Nitrile	50 pairs/box 4 boxes/case
2D73DP60	6						
2D73DP65	6½						
2D73DP70	7						
2D73DP75	7½						
2D73DP80	8						
2D73DP85	8½						
2D73DP90	9						

## Protexis™ Neoprene Essential Surgical Gloves



Code	Size	Material	Color	Category	Fingertip thickness	Coatings	Sales Unit
2D73DS55	5½	Synthetic neoprene	Light brown	Standard	0.17 mm	Nitrile	50 pairs/box 4 boxes/case
2D73DS60	6						
2D73DS65	6½						
2D73DS70	7						
2D73DS75	7½						
2D73DS80	8						
2D73DS85	8½						
2D73DS90	9						

## Protexis™ Latex, Powder-Free

### Protexis™ Latex Surgical Gloves




Code	Size	Material	Color	Category	Fingertip thickness	Coatings	Sales Unit
2D72NS55X	5½	Latex	Light brown	Standard	0.25 mm	Nitrile	50 pairs/box
2D72NS60X	6						4 boxes/case
2D72NS65X	6½						
2D72NS70X	7						
2D72NS75X	7½						
2D72NS80X	8						
2D72NS85X	8½						
2D72NS90X	9						

### Protexis™ Latex Essential Surgical Gloves



Code	Size	Material	Color	Category	Fingertip thickness	Coatings	Sales Unit
2D72LE55	5½	Latex	Cream	Standard	0.19 mm	Nitrile	50 pairs/box
2D72LE60	6						4 boxes/case
2D72LE65	6½						
2D72LE70	7						
2D72LE75	7½						
2D72LE80	8						
2D72LE85	8½						
2D72LE90	9						

### Protexis™ Latex Micro Surgical Gloves



Code	Size	Material	Color	Category	Fingertip thickness	Coatings	Sales Unit
2D72NT55X	5½	Latex	Light brown	Thin	0.17 mm	Nitrile	50 pairs/box
2D72NT60X	6						4 boxes/case
2D72NT65X	6½						
2D72NT70X	7						
2D72NT75X	7½						
2D72NT80X	8						
2D72NT85X	8½						
2D72NT90X	9						

### Protexis™ Latex with Neu-Thera™ Surgical Gloves



Code	Size	Material	Color	Category	Fingertip thickness	Coatings	Sales Unit
2D73TP55	5½	Latex	Light brown	Standard	0.25 mm	Nitrile and Neu-Thera™ Emollient Coating	50 pairs/box
2D73TP60	6						4 boxes/case
2D73TP65	6½						
2D73TP70	7						
2D73TP75	7½						
2D73TP80	8						
2D73TP85	8½						
2D73EB90	9						

Protexis™ Latex Blue with Neu-Thera™ Surgical Gloves



Code	Size	Material	Color	Category	Fingertip thickness	Coatings	Sales Unit
2D72LU55	5½	Latex	Blue	Underglove	0.19 mm	Nitrile and Neu-Thera™ Emollient Coating	50 pairs/box 4 boxes/case
2D72LU60	6						
2D72LU65	6½						
2D72LU70	7						
2D72LU75	7½						
2D72LU80	8						
2D72LU85	8½						
2D72LU90	9						

Protexis™ Latex Ortho Surgical Gloves



Code	Size	Material	Color	Category	Fingertip thickness	Coatings	Sales Unit
2D72LT60	6	Latex	Brown	Thick	0.34 mm	Hydrogel	40 pairs/box 4 boxes/case
2D72LT65	6½						
2D72LT70	7						
2D72LT75	7½						
2D72LT80	8						
2D72LT85	8½						
2D72LT90	9						

# Surgical Suction



Medi-Vac™ Suction and Fluid  
Collection Products

Argyle™ Surgical Suction Yankauers

Argyle™ Collection System

## Medi-Vac™ Flex Advantage™ Hardware

**Roll Stands** – Roll stands allow flexibility in the location and configuration of Medi-Vac™ Flex Advantage™ Canister installations.

**Canisters** – Reusable outer canisters used to hold disposable flexible liners

### Hardware / Roll Stands



Code	Description	Sales Unit
65652-586A	Universal Rollstand 76 cm	1 ea/cs
65652-596A	Regulator Mobile Suction Cart – an accessory will have to be added	1 ea/cs
65652-145A	Wall Plate W/Holes	60 ea/cs
65652-002A	Two Canister Manifold	10 ea/cs
65652-003A	Three Canister Manifold	10 ea/cs
65652-004A	Four Canister Manifold	10 ea/cs

### Canisters

Code	Description	Sales Unit
65652-611A	Medi-Vac™ 1000 ml Flex Canister with Internal On-Off Valve	10 ea/cs
65652-616A	Medi-Vac™ 1500 ml Flex Canister with Internal On-Off Valve	10 ea/cs
65652-631A	Medi-Vac™ 3000 ml Flex Canister with Internal On-Off Valve	8 ea/cs
65652-916A	Medi-Vac™ 1000 ml Flex Canister with T-Connection	10 ea/cs
65652-915A	Medi-Vac™ 1500 ml Flex Canister with T-Connection	10 ea/cs
65652-914A	Medi-Vac™ 3000 ml Flex Canister with T-Connection	8 ea/cs

### Other Hardware Options



Code	Description	Sales Unit
65652-911	Flex Advantage™ single canister base stand 1000 ml/1500 ml	1 ea/cs
65652-913A	Flex Advantage™ single canister base stand 3000 ml	10 ea/cs
65652-924A	Flex Advantage™ single canister roll stand 1000 ml/1500 ml	3 ea/cs
65652-928A	Flex Advantage™ double canister roll stand 1500 ml/3000 ml	10 ea/cs
65652-686A	Dual Head Rollstand	1 ea/cs

### Special Volume Options



Code	Description	Sales Unit
65651-120	Medi-Vac™ Guardian™ LVC Suction	4 ea/cs
65652-583A	Guardian™ LVC Rollstand – 96 cm	1 ea/cs
1-700CE	International Universal Critical Measurement – 300 ml	50 ea/cs



## Medi-Vac™ Flex Advantage™ Single Use Liners




**Single Use Liners** – Flex Advantage™ Flexible Liners are used to collect biological fluids by aspiration and have been designed and built for “high flow and high vacuum” applications (EN ISO 10079-3).

### Single Use Liners

Code	Description	Sales Unit
65651910	Flexible Liner, 1000 cc	50 ea/cs
65651920	Flexible Liner, 1500 cc	50 ea/cs
65651930	Flexible Liner, 3000 cc	50 ea/cs
65651910PG	Flexible Liner, 1000 cc + Gel Vac	45 ea/cs
65651920PG	Flexible Liner, 1500 cc + Gel Vac	45 ea/cs
65651930PG	Flexible Liner, 3000 cc + Gel Vac	45 ea/cs
65651916	Flexible Liner 1000 cc + 6 mm x 1.8 m tubing	50 ea/cs
65651926	Flexible Liner 1500 cc + 6 mm x 1.8 m tubing	50 ea/cs
65651929	Flexible Liner 1500 cc + 6 mm x 2.7 m tubing	50 ea/cs
65651936	Flexible Liner 3000 cc + 6 mm x 1.8 m tubing	50 ea/cs
65651939	Flexible Liner 3000 cc + 6 mm x 2.7 m tubing	50 ea/cs
65653916	Flexible Liner 1000 cc + 6 mm x 1.8 m tubing + Medistop Adult	50 ea/cs
65653926	Flexible Liner 1500 cc + 6 mm x 1.8 m tubing + Medistop Adult	50 ea/cs
65653929	Flexible Liner 1500 cc + 6 mm x 2.7 m tubing + Medistop Adult	50 ea/cs
65654916A	Flexible Liner 1000 cc + 6 mm x 1.8 m tubing + Medistop with tip	50 ea/cs
65654926A	Flexible Liner 1500 cc + 6 mm x 1.8 m tubing + Medistop with tip	50 ea/cs

## Medi-Vac™ Flex Advantage™ Accessories


**Tubing** – Flexible, connecting tubing to adapt your system for various setups.




### Wall Suction Tubing

Code	Description	Sales Unit
516500CE	Medi-Vac™ Universal Vacuum Connection Tubing, 8 mm I.D. x 30.5 m	4 ea/cs
<b>TANDEM</b>		
2-620CE	Flex Advantage™ Tandem Tube	30 ea/cs

### Specimen Socks



Code	Description	Sales Unit
2-122CE	Standard Specimen Sock	10 ea/cs
2-123CE	Short Specimen Sock	10 ea/cs
2-124CE	Large-port Specimen Sock	10 ea/cs
2-169CE	Female Specimen Sock	10 ea/cs





### Vacuum Control

Code	Description	Sales Unit
65655853	Medistop Vacuum Suction Regulator	200 ea/cs
65655854A	Medistop with Tip Vacuum Suction Regulator	200 ea/cs
NS366	Medistop Vacuum Suction Regulator Adult + 1.8 m tubing	30 ea/cs
NS369	Medistop Vacuum Suction Regulator Adult + 2.7 m tubing	30 ea/cs
NS466A	Medistop with Tip Vacuum Suction Regulator + 1.8 m tubing	30 ea/cs



### Gelling Agent


Code	Description	Sales Unit
GELVAC15	MediVac™ Gelling Agent – 15 g per sachet	1 CS containing 2 boxes of 100 sachets
SOLID500E	MediVac™ Solidifies up to 500 cc Bottle	96 ea/cs
SOLID1200E	MediVac™ Solidifies up to 1200 cc Bottle	96 ea/cs
SOLID1500E	MediVac™ Solidifies up to 1500 cc Bottle	96 ea/cs
SOLID2000E	MediVac™ Solidifies up to 2000 cc Bottle	96 ea/cs
SOLID3000E	MediVac™ Solidifies up to 3000 cc Bottle	96 ea/cs

This range of products are not medical devices and they are not CE marked.



### Accessories

Code	Description	Sales Unit
65652-109	Right angle connector For Flex Advantage™ patient port and CRD™/Guardian™	50 ea/cs
65652-119	Ortho port converter Accepts up to 0.95 cm I.D. tubing with male port	10 ea/cs



## Medi-Vac™ Guardian™

The Guardian™ Suction Canister System features single-use disposable hard locking lid attachments. Special brackets and stands are available for convenient mounting configurations and mobility. Designed for safety, economy and convenience, Guardian™ Suction Canister Systems are available in a variety of sizes to meet various suction and fluid collection needs.



### Guardian™ Suction Canisters

Code	Description	Sales Unit
65651-212	Rigid canister with mechanical filter – 1200 ml – 1200 cc	40/cs
65651-220	Rigid canister with mechanical filter – 2000 ml – 2000 cc	40/cs
65651-230	Rigid canister with mechanical filter – 3000 ml – 3000 cc	40/cs
65651-312	Guardian™ rigid canister with mechanical filter and tubing 1200 ml + 5 mm x 1.8 m	30 ea/cs
65651-896	Guardian™ rigid canister with mechanical filter and tubing 1200 ml + 6 mm x 1.8 m	30 ea/cs
65651-320	Guardian™ rigid canister with mechanical filter and tubing 2000 ml + 5 mm x 1.8 m	30 ea/cs
65651-395	Guardian™ rigid canister with mechanical filter and tubing 1200 ml + 5 mm x 1.8 m, 5 mm x 0.5 m	30 ea/cs
1-319CE	Guardian™ suction canister with tubing 1200 ml + mechanical inner diameter 6.4 mm x 2.7 m	30 ea/cs

Reusable polycarbonate bracket used to support the Guardian™ Canister in any slide mount such as the wall plate or mobile suction cart. Comes with on/off valve built into bracket. Has pre-attached tubing for connection to canister.



### Guardian™ Ring Bracket with Built-in On/Off Valve

Code	Description	Sales Unit
65652-159A	Guardian™ Ring Bracket with Internal On-Off Valve, 1200 ml	10 ea/cs
65652-139	Guardian™ Ring Bracket with Internal On-Off Valve, 2000 ml	1 ea/cs
2-542ACE	CRD™ and Guardian™ Tandem Tubes 46 cm x 8 mm I.D.	90 ea/cs





## Argyle™ Surgical Suction Yankauers

A comprehensive range of suction cannulas designed to meet various needs in the operating room.

### Argyle™ Yankauers


Complete range of rigid and flexible single use yankauers for suction of fluids

### Argyle™ Flexible Yankauers

- Sterile, single use
- Flexible yankauers can be shaped by the user to modify angle and adapt to specific operative sites



Code	Description	Vacuum Control	Fr/Ch	Length	Sales Unit (Box/packaging)
<b>ARGYLE™ MINI YANKAUER DEHP FREE 16 CM LENGTH</b>					
1180113001	No side holes	No	16	16 cm	100
1180113050	4 side holes	Yes	16	16 cm	100
1180113100	No side holes	Yes	16	16 cm	100
<b>ARGYLE™ FLEXIBLE YANKAUER DEHP FREE 25 CM LENGTH</b>					
8888504001	Fine tip	No	12	25 cm	50
8888504019	Fine tip	Yes	12	25 cm	50
8888501007	Medium tip	No	18	25 cm	50
8888501023	Medium tip	Yes	18	25 cm	50
8888501015	Medium tip with sump	No	18	25 cm	50
8888501031	Medium tip with sump	Yes	18	25 cm	50
1180501205	Medium tip with sure-grip	No	18	25 cm	50
8888502005	Large tip	No	22	25 cm	50
8888502013	Large tip with sump	Yes	22	25 cm	50
1180502039	Large tip with sure-grip	No	22	25 cm	50
1180504018	Flat tip	No	22	25 cm	50
1180503003	Extra large tip	No	30	25 cm	50
<b>ARGYLE™ FLEXIBLE YANKAUER SETS WITH NON-CONDUCTIVE FLEXIBLE TUBING</b>					
1180504613	Fine tip 5 mm dia tube	No	12	2.2 m**	25
1180504639	Fine tip 5 mm dia tube	No	12	3.0 m**	25
1180504621	Fine tip 5 mm dia tube	No	12	4.0 m**	25
1180501619	Medium tip 7 mm dia tube	No	18	2.2 m**	25
1180501635	Medium tip 7 mm dia tube	No	18	3.0 m**	25
1180501700	Medium tip 7 mm dia tube	Yes	18	3.0 m**	25
1180502617	Large tip 7 mm dia tube	No	22	2.2 m**	25
1180502625	Large tip 7 mm dia tube	No	22	3.0 m**	25
1180502618	Large tip 7 mm dia tube	No	22	2.2 m**	25
1180502658	Large tip 7 mm dia tube	No	22	3.0 m**	25
8888501502	Medium tip 7 mm dia tube	No	18	1.8 m**	50
8888501528	Medium tip with sump 7 mm dia tube	Yes	18	1.8 m**	50
8888502518	Large tip with sump 7 mm dia tube	Yes	22	1.8 m**	50



## Medivac™ Yankauers

Code	Description	Sales Unit
KEX60	Flexible Suction handle, Fine capacity	50 ea/cs
KEX61	Flexible Suction handle, Regular capacity	50 ea/cs
KEX61V	Flexible suction handle, regular capacity, vacuum control	50 ea/cs
KEX80	Yankauer, tapered bulbous tip – Ch. 26/29 cm	50 ea/cs
KEX82	Yankauer, tapered bulbous tip, vacuum control – Ch. 26/29 cm	50 ea/cs
KEX83A	Yankauer, tapered bulbous tip with 6 mm x 1.8 m sterile tubing	50 ea/cs
KEX84	Yankauer, tapered bulbous tip with 6 mm x 2.7 m sterile tubing	25 ea/cs
KEX85	Yankauer, tapered bulbous tip with 6 mm x 3.7 m sterile tubing	25 ea/cs



## Argyle™ Rigid Yankauers

- Sterile, single use

Code	Description	Vacuum Control	Fr/Ch	Length	Sales Unit (Box/packaging)
8888505115	Large tip	No	22	25 cm	50
8888505123	Large tip	Yes	22	25 cm	50
8888505016	Bulb tip	No	34	26 cm	50
8888505024	Bulb tip	Yes	34	26 cm	50



## Argyle™ Selec-Trol Yankauers

- Sterile, single use
- Available with large open tip, bulb tip or fine tip
- Open and bulb tips are made from transparent high quality PVC with lateral holes, rounded distal ends and bevelled edges for controlled atraumatic suction

Code	Description	Vacuum Control	Fr/Ch	Length	Sales Unit (Box/packaging)
1180505172	Large tip	On-off switch	24	27 cm	20
1180505164	Bulb tip	On-off switch	32	27 cm	20



## Argyle™ Poole Suction Instruments

- Sterile, single use
- Two-piece fixed construction: transparent inner and outer lumen with fingertip suction control

Code	Description	Fr/Ch	Length	Sales Unit (Box/packaging)
<b>ARGYLE™ POOLE INSTRUMENTS</b>				
8888508002	Flexible poole	34	27 cm	50
8888509703	Rigid poole	14-28	25 cm	50
<b>ARGYLE™ POOLE SETS</b>				
8888509729	Rigid poole with 5 mm dia cond tube	14-28	1.8 m**	50

## Argyle™ Connecting Tubing

### Argyle™ Suction Connecting Tubes

Sterile connecting tubes in medical grade PVC

Code	Diameter and Length	Sales Unit (Box/packaging)
<b>ARGYLE™ CONNECTING TUBE WITH MOLDED CONNECTORS</b>		
8888301507	5 mm x 50 mm	50
8888301515	5 mm x 1.8 m	50
8888301523	5 mm x 3.1 m	50
8888301531	5 mm x 3.7 m	20
8888301606	6 mm x 1.8 m	50
8888301614	6 mm x 3.1 m	50
8888301622	6 mm x 3.7 m	20
8888301705	7 mm x 1.8 m	50
8888301713	7 mm x 3.1 m	50
8888301721	7 mm x 3.7 m	20

## Argyle™ Connecting Tubing

### Argyle™ Suction Connecting Tubes

Sterile connecting tubes in medical grade PVC

Code	Diameter and Length	Sales Unit (Box/packaging)
<b>ARGYLE™ CONNECTING TUBES WITH FUNNEL TO FUNNEL CONNECTORS</b>		
8888284505	5 mm x 50 cm	50
8888284513	5 mm x 1.8 m	50
8888284539	5 mm x 3.1 m	50
8888284547	5 mm x 3.7 m	20
8888284554	5 mm x 6.1 m	10
8888284604	6 mm x 1.8 m	50
8888284612	6 mm x 3.1 m	50
8888284620	6 mm x 3.7 m	20
8888284638	6 mm x 6.1 m	10
8888284703	7 mm x 1.8 m	50
8888284737	7 mm x 3.7 m	20
8888284745	7 mm x 6.1 m	10
<b>ARGYLE™ CONNECTING TUBES WITH FUNNEL TO TAPERED END</b>		
8888284521	5 mm x 1.8 m	50
8888284711	7 mm x 1.8 m	50


## Argyle™ Surgical Suction Connecting Tubes

### Argyle™ Bubble Tubing

Can be easily cut to exact length and is available in conductive and non-conductive variations



Code	Diameter and Length	Sales Unit
8888280206	Argyle™ Bubble tubing 5 mm x 30 m – Non Conductive	1
8888280214	Argyle™ Bubble Tubing 5 mm x 30.5 m	1
1180839035	Argyle™ Bubble Tubing 6 mm x 50 m	1
8888280602	Argyle™ Bubble tubing 7 mm x 30 m – Non Conductive	1
8888280610	Argyle™ Bubble Tubing 7 mm x 30.5 m	1



## Medivac™ Suction Tubing

### Medivac™ Suction Tubing

Code	Diameter and Length	Sales Unit
NEX510	5 mm x 3.05 m Sterile Tubing	30 ea/cs
NEX56A	5 mm x 1.8 m Sterile Tubing	50 ea/cs
NEX610	6 mm x 3.05 m Sterile Tubing	25 ea/cs
NEX66A	6 mm x 1.8 m Sterile Tubing	45 ea/cs
NEX710	7 mm x 3.05 m Sterile Tubing	20 ea/cs
NEX712	7 mm x 3.60 m Sterile Tubing	20 ea/cs
NEX76A	7 mm x 1.8 m Sterile Tubing	40 ea/cs



## Argyle™ Specimen Trap

Latex free specimen trap for bronchial secretions

Code	Description	Capacity	Sales Unit (Box/package)
<b>ARGYLE™ SPECIMEN TRAP</b>			
8884724504	Specimen trap 70cc	70 ml	50
8884724500	Specimen trap 40cc	40 ml	50



## Argyle™ Lukens Mucus Collection

The Argyle™ Lukens specimen container comes with a whistle tip suction catheter with a tapered connector and a smooth Slide-Tex finish.

This product contains DEHP.

Tip

### Argyle™ Lukens Mucus Collection

Code	Description	Sales Unit (Box/package)
8888258608	Argyle™ Lukens Mucus Collection – 14Fr Catheter – 20 ml	50



Tapered Connector

## Argyle™ Tubing Connectors

### Argyle™ Connectors/Adaptors

Connectors allow customization of your suction system connecting sections of tubing together for various applications



Code	Description	Adaptable to Tube Diameters	Sales Unit (Box/package)
<b>ARGYLE™ ADAPTORS</b>			
8888270207	Vinyl for sterile tubes SIMS type	5-11 mm	100
8888271007	Vinyl for sterile tubes 5-in-1 type	5-11 mm	100
8888271411	Sterile polypropylene 5-in-1 type	8-14 mm	100
8888271502	Sterile nylon 5-in-1 type	8-14 mm	100
8888272013	Sterile polypropylene large 'Y' type	6-13 mm	100
8888270306	Vinyl for clean tubes SIMS type	5-11 mm	100
8888271403	Clean polypropylene 5-in-1 type	8-14 mm	100
8888271510	Clean nylon 5-in-1 type	8-14 mm	100
8888272005	Clean polypropylene large Y type	6-13 mm	50

# Thermometry



Genius™ 3 Tympanic  
Thermometer and Base

Filac™ Electronic Thermometer

## Genius™ 3 Tympanic Thermometer and Base



### Genius™ 3 Tympanic Thermometer and Base

Code	Description	Sales Unit (Box/packaging)
303013	Genius™ 3 Tympanic Thermometer and Base	1



### Genius™ 3 Tympanic Thermometer and Base – Accessories

Code	Description	Sales Unit (Box/packaging)
303030	Genius™ probe covers	2112 (22 x 96)
303097	Genius™ checker/calibrator	1
303058*	Locking wall mount	1
303059*	Cart with locking mount	1
F31950-SP*	Genius™ replacement cord	1
PT00058372*	Genius™ 3 Calibrator Cable	1

\* This range of products are not medical devices and they are not CE marked.

## Filac™ Electronic Thermometer



### Filac™ 3000 Electronic Thermometer

Code	Description	Sales Unit (Box/packaging)
502000	Filac™ 3000 probe covers	2000 (100 x 20)
500028*	Blue oral/axillary isolation chamber	1
500026*	Blue oral/axillary probe with 122 cm cord	1
500036*	Red rectal probe with 122 cm cord	1
500027*	Blue oral/axillary probe with 275 cm cord	1
500037*	Red rectal probe with 275 cm cord	1
500099*	Calibration plug	1
500046*	Mounting base for wall	1

\* This range of products are not medical devices and they are not CE marked.

# Urology



Silicone Foley Catheters  
Collection Devices  
Urine Meters  
Pre-connected Systems



## Silicone Foley Catheters

### Dover™ 100% Silicone Foley Catheter

- 100% silicone
- Standard/nelaton tip
- Sterile

Code	Fr/Ch	Length	Sales Unit
<b>2 WAY 1.5 ML BALLOON</b>			
8887603069	6	32 cm	10
<b>2 WAY 3 ML BALLOON</b>			
8887603085	8	32 cm	10
8887603101	10	32 cm	10
<b>2 WAY 5-10 ML BALLOON</b>			
8887605122	12	42 cm	10
8887605148	14	42 cm	10
8887605163	16	42 cm	10
8887605189	18	42 cm	10
8887605205	20	42 cm	10
8887605221	22	42 cm	10
8887605247	24	42 cm	10
8887605262	26	42 cm	10
<b>2 WAY 30 ML BALLOON</b>			
8887630187	18	42 cm	10
8887630203	20	42 cm	10
8887630229	22	42 cm	10
8887630245	24	42 cm	10
<b>3 WAY 30 ML BALLOON</b>			
8887665183	18	42 cm	10
8887665209	20	42 cm	10
8887665225	22	42 cm	10
<b>2 WAY 5-10 ML BALLOON – FEMALE CATHETER</b>			
8887615121	12	25 cm	10
8887615147	14	25 cm	10
8887615162	16	25 cm	10
8887615188	18	25 cm	10



## Silicone Foley Catheters

### Dover™ 100% Silicone Foley Catheter – Tiemann

- 100% silicone
- Tiemann tip
- Sterile

Code	Fr/Ch	Length	Sales Unit
<b>2 WAY 5-10 ML BALLOON</b>			
8887205121	12	44 cm	10
8887205147	14	44 cm	10
8887205162	16	44 cm	10
8887205188	18	44 cm	10
8887205204	20	44 cm	10
8887205220	22	44 cm	10
8887205246	24	44 cm	10



### Dover™ 100% Silicone Foley Catheter with Glyco-Block Pre-filled Syringe

- 100% silicone
- Standard/nelaton tip
- 10 ml sterile water/10% glycerine pre-filled syringe
- Sterile

Code	Fr/Ch	Length	Sales Unit
<b>2 WAY 5-10 ML BALLOON</b>			
8887405127	12	42 cm	5
8887405143	14	42 cm	5
8887405168	16	42 cm	5
8887405184	18	42 cm	5
8887405200	20	42 cm	5
8887405226	22	42 cm	5





## Silicone Foley Catheters

### Dover™ 100% Silicone Foley Catheter with Sterile Water Pre-filled Syringe

- 100% silicone
- Standard/nelaton tip
- 10 ml sterile water pre-filled syringe
- Sterile

Code	Fr/Ch	Length	Sales Unit
<b>2 WAY 3 ML BALLOON</b>			
8887103084	8	32 cm	5
8887103102	10	32 cm	5
<b>2 WAY 5-10 ML BALLOON</b>			
8887105123	12	42 cm	5
8887105149	14	42 cm	5
8887105164	16	42 cm	5
8887105180	18	42 cm	5
8887105206	20	42 cm	5
<b>2 WAY 5-10 ML BALLOON – FEMALE CATHETER</b>			
8887115122	12	25 cm	5
8887115148	14	25 cm	5
8887115163	16	25 cm	5




### Dover™ 100% Silicone Foley Catheter Community Catheter

- 100% silicone
- Standard/nelaton tip
- 10 ml sterile water pre-filled syringe
- 10 ml empty syringe
- Sterile

Code	Fr/Ch	Length	Sales Unit
<b>2 WAY 3 ML BALLOON</b>			
8887703083	8	32 cm	1
8887703109	10	32 cm	1
<b>2 WAY 5-10 ML BALLOON</b>			
8887705120	12	42 cm	1
8887705146	14	42 cm	1
8887705161	16	42 cm	1
8887705187	18	42 cm	1
8887705203	20	42 cm	1

## Collection Devices



### Dover™ Leg Straps

Code	Description	Sales Unit
31377	Self-adhesive adjustable leg straps	1 pair



### Dover™ Bedside Drainage Bags

- Integrated bedside hanger
- Anti-reflux device
- 2000 ml capacity
- Needle free sample port

Code	Description	Sales Unit
900010	120 cm inlet tubing, T Tap outlet valve	20



### MonoFlo™ Bedside Drainage Bags

- Integrated bedside hanger
- Mono-Flo™ anti-reflux device
- 2000 ml capacity
- Sample port
- Sterile

Code	Description	Sales Unit
6208-02	90 cm drainage inlet tube	30
6209-02	120 cm drainage inlet tube	30
6308	120 cm drainage inlet tube, Safeguard™ needle-free sample port and Splashguard™ II™ enclosed drain spout	20



## Collection Devices

### Dover™ Standard Bedside Drainage Bags

- Anti-reflux concept
- Allows for hanging the bag on a bed frame or sideboard
- Sample port
- Indicates that only the interior of the drainage bag and fluid pathway are guaranteed sterile

Code	Description	Sales Unit
3512	90 cm drainage inlet tube, 2000 ml	20



### Dover™ TURP Bag

- Anti-reflux concept
- Integrated bedside hanger
- Sample port
- Sterile

Code	Description	Sales Unit
10082EU	90 cm drainage inlet tube, 4000 ml	25



## Urine Meters

### Dover™ Urine Meter

- 400 ml meter
- 2000 ml drainage bag
- Sample port
- Integrated bedside hanger
- Sterile

Code	Description	Sales Unit
2000	150 cm drainage inlet tube, Splashguard™ II™ enclosed drain spout	10
7000	150 cm drainage inlet tube, Safeguard™ Needle-free sample port, Splashguard II™ enclosed drain spout	10
7000LL	150 cm drainage inlet tube, Safeguard™ Needle-free sample port, Splashguard II™ enclosed drain spout, with luer lock sample port	10



### Dover™ Urine Meter

- 200 ml meter
- 2000 ml drainage bag
- Sample port
- Integrated bedside hanger
- Sterile

Code	Description	Sales Unit
6500	120 cm drainage inlet tube	20
6502	150 cm drainage inlet tube	20
6500NF	120 cm drainage inlet tube, needle-free sample port	20
6502NF	150 cm drainage inlet tube, needle-free sample port	20



## Urine Meters

### Dover™ Urine Meter

- 500 ml meter
- 2600 ml drainage bag
- Anti-reflux valve
- Needle-free sample port
- Integrated bedside hanger
- Sterile

Code	Description	Sales Unit
900011	120 cm inlet tubing , T Tape outlet valve	1



### Dover™ Urine Meter

- 150 ml meter
- Sample port
- Anti-reflux device
- Integrated bedside hanger
- Sterile

Code	Description	Sales Unit
6600	120 cm drainage inlet tube	10



### Dover™ Urine Meter with Exchangeable Bag

- 500 ml meter
- 2000 ml drainage bag
- Needle-free sample port
- Integrated bedside hanger
- Sterile

Code	Description	Sales Unit
900012	120 cm inlet tubing, T tape outlet valve	10
900013	2000 ml replacement closed bag	100
900014	2000 ml replacement drainable bag	100

## Pre-connected Systems



### Mono-Flo™ Pre-connected Drainage Bags

- Dover™ 100% silicone foley catheter
- 2000 ml capacity
- Tamper evident seal
- Sample port
- Mono-Flo™ anti-reflux device
- Sterile

Code	Description	Sales Unit
6209121	12 Fr/Ch catheter	20
6209141	14 Fr/Ch catheter	20
6209161	16 Fr/Ch catheter	20
6209181	18 Fr/Ch catheter	20



### Dover™ Pre-connected Urine Meters

- 100% silicone foley catheter
- 200 ml urine meter
- Tamper evident seal
- Sample port

Code	Description	Sales Unit
6500121	12 Fr/Ch catheter	20
6500141	14 Fr/Ch catheter	20
6500161	16 Fr/Ch catheter	20
6500181	18 Fr/Ch catheter	20

# Wound Care



Antimicrobial Wound dressings

Foam Dressings

Alginate Dressings

Hydrogel Dressings

Film Dressings

Non-Adherent Dressings – Moist

Non-Adherent Dressings

Non Woven Dressings

Gauzes and Sponges

Bandages

Ancillary Products

## Antimicrobial Wound Dressings

### Kendall™ AMD Antimicrobial Foam Dressings



Kendall™ AMD Antimicrobial Foam Dressings provide a moist wound-healing environment for a wide variety of wounds. The dressing contains the antimicrobial agent PHMB that provides a barrier to bacteria penetration through the dressing and prevents colonization and proliferation of bacteria within the dressing for up to 7 days. Kendall™ AMD Antimicrobial Foam Dressings are indicated for use in management of postsurgical incisions, pressure sores, venous stasis ulcers, diabetic ulcers, donor sites, abrasions, lacerations, 1st and 2nd degree burns, dermatologic disorders, other wounds inflicted by trauma and, as a secondary dressing or cover dressing for packed wounds.

Code	Description	Qty Box	Qty Case
55522AMDX	5.1 cm x 5.1 cm	25	100
55544AMDX	10.2 cm x 10.2 cm	10	50
55548AMDX	10.2 cm x 20.3 cm	10	50
55566AMDX	15.2 cm x 15.2 cm	10	50
55588AMDX	20.3 cm x 20.3 cm	10	50
55535PAMD	8.8 cm x 7.5 cm Fenestrated, with Topsheet	10	50
55544PAMD	10 cm x 10 cm with Topsheet	10	50



### Kendall™ AMD Antimicrobial Foam – Bordered Dressing

Code	Description	Qty box	Qty case
55523BAMDX	4.4 cm x 8.3 cm, pad size 2.5 cm x 4.4 cm	10	50
55544BAMDX	8.9 cm x 8.9 cm, pad size 5.1 cm x 5.1 cm	10	50
55546BAMDX	8.9 cm x 14 cm, pad size 5.1 cm x 10.2 cm	10	50
55566BAMDX	14 cm x 14 cm, pad size 10.2 cm x 10.2 cm	10	50
55588BAMDX	19 cm x 19 cm, pad size 15.2 cm x 15.2 cm	10	50



### Kendall™ AMD Antimicrobial Foam – Disc

Code	Description	Qty Box	Qty Case
55511AMDX	2.54 cm diameter, 4 mm hole	10	40
55512AMDX	2.54 cm diameter, 7 mm hole	10	40



## Antimicrobial Wound Dressings



### Telfa™ AMD Antimicrobial Non-Adherent Dressing

- For use as primary dressings for surgical incisions, lacerations, abrasions, burns, donor sites, catheter sites, I.V. sites and central lines
- The antimicrobial activity of the PHMB in Telfa™ AMD dressing helps to resist bacterial colonization within the dressing and inhibit bacterial penetration through the dressing
- Sterile

Code	Description	Qty Box	Qty Case
<b>TELFA™ AMD PADS STERILE</b>			
7662	7.6 cm x 11.4 cm	50	900
7663	7.6 cm x 20.3 cm	50	600
<b>TELFA™ AMD ISLAND DRESSINGS</b>			
7665	10.2 cm x 12.7 cm (pad size 5 cm x 7.5 cm)	25	200
7666	10.2 cm x 20.3 cm (pad size 5 cm x 15 cm)	25	100
7667	10.2 cm x 25.4 cm (pad size 5 cm x 20 cm)	25	100
7668	10.2 cm x 35.6 cm (pad size 5 cm x 30 cm)	25	50



### Excilon™ AMD Antimicrobial Dressings

- For use as primary dressings for IV sites, tracheostomy tube sites, chest tube sites, catheter sites and drain sites
- The antimicrobial activity of the PHMB in Excilon AMD helps to resist bacterial colonization within the dressing and inhibit bacterial penetration through the dressing
- Packed in pairs for improved wound coverage and protection
- Sterile

Code	Description	Qty Box	Qty Case
7088	10.2 cm x 10.2 cm pad, drain sponge	25 Pairs	300 Pairs
7089	5.1 cm x 5.1 cm pad I.V. sponge	35 Pairs	700 Pairs



### Kerlix™ AMD

- For use as primary or secondary dressings for exuding wounds, surgical incisions, lacerations, abrasions, burns, wound packing, donor sites, catheter sites, I. V. sites and central lines. Also may be used for securement of primary dressings
- The antimicrobial activity of the PHMB in Kerlix AMD helps resist bacterial colonization within the dressing and inhibit bacterial penetration through the dressing
- Sterile

Code	Description	Qty Box	Qty Case
3331	Roll 11.4 cm x 3.7m	1	60
3332	Roll 11.4 cm x 3.7m soft pack	1	100
6662	Super sponge 15.2 cm x 17.1 cm	2	480
6665	Super sponge 15.2 cm x 17.1 cm	5	600
6660	Super sponge 15.2 cm x 17.1 cm	10	480



## Foam Dressings

### Kendall™ Hydrophilic Foam Dressings

- Highly absorbent hydrophilic foam wound dressing
- Soft, non-adherent pad
- Flexible nature allowing it to conform easily to all body contours
- Wide variety of convenient sizes

Code	Description	Qty Box	Qty Case
<b>KENDALL™ FOAM</b>			
55535	Fenestrated 8.9 cm x 7.6 cm	10	50
55522	5.1 cm x 5.1 cm	25	100
55533	7.6 cm x 7.6 cm	10	50
55544	10.2 cm x 10.2 cm	10	50
55548	10.2 cm x 20.3 cm	10	50
55555	12.7 cm x 12.7 cm	10	50
55566	15.2 cm x 15.2 cm	10	50
55588	20.3 cm x 20.3 cm	10	50
<b>KENDALL™ FOAM PLUS</b>			
55535P	Fenestrated 8.9 cm x 7.6 cm	10	50
55522P	5.1 cm x 5.1 cm	25	100
55533P	7.6 cm x 7.6 cm	10	50
55544P	10.2 cm x 10.2 cm	10	50
55548P	10.2 cm x 20.3 cm	10	50
55566P	15.2 cm x 15.2 cm	10	50
55588P	20.3 cm x 20.3 cm	10	50
<b>KENDALL™ FOAM ISLAND</b>			
55544B	5.1 cm x 5.1 cm pad	10	50
55566B	10.2 cm x 10.2 cm pad	10	50
55588B	15.2 cm x 15.2 cm pad	10	50



### Kendall™ Foam Gentle Border

Code	Description	Qty Box	Qty Case
55523BG	4.4 cm x 8.3 cm, pad size 2.5 cm x 4.4 cm	10	50
55544BG	8.9 cm x 8.9 cm, pad size 5.1 cm x 5.1 cm	10	50
55546BG	8.9 cm x 14.0 cm, pad size 5.1 cm x 10.2 cm	10	50
55566BG	14.0 cm x 14.0 cm, pad size 10.2 cm x 10.2 cm	10	50
55588BG	19.1 cm x 19.1 cm, pad size 15.2 cm x 15.2 cm	10	50



## Alginate Dressings

### Kendall™ Calcium Alginate Dressings

- Absorbs exudate from wound, control minor bleeding and provide a moist healing environment
- Sterile
- Plus version available for increased absorbency

Code	Description	Qty Box	Qty Case
<b>STERILE 1'S IN PEEL-BACK PACKAGE</b>			
9231	30.5 cm rope	5	20
9232	5.1 cm x 5.1 cm	10	100
9233	10.2 cm x 10.2 cm	10	50
9236	10.2 cm x 10.2 cm plus	10	50
9240	10.2 cm x 14 cm	10	50
9238	10.2 cm x 20.3 cm	5	50
9239	15.2 cm x 25.4 cm	10	50
9242	30.5 cm x 61 cm	5	5
9243	61 cm rope	5	20
9244	91.4 cm rope	5	20



## Hydrogel Dressings

### Kendall™ Hydrogel Breast Feeding Pad

- Indicated for nipple trauma associated with breastfeeding
- Distinctive disc shape conforms in a cup-like manner for nipple coverage

Code	Description	Qty Pack	Qty Box	Qty Case
9403	7.6 cm disc	1 pr.	5 pr.	25 pr.



## Film Dressings

### Kendall™ Transparent Film Dressings

- Dressing for I.V. sites, donor sites, burns, ulcers and sutures
- Simple, effective delivery system allows one-handed application
- Highly permeable to oxygen and moisture vapour

Code	Description	Qty Box	Qty Case
6651	3.8 cm x 3.8 cm	100	400
6640	5.1 cm x 7 cm	100	400
6641	10.2 cm x 12.1 cm	50	200
6647	10.2 cm x 20.3 cm	20	120
6648	20.3 cm x 25.4 cm	10	12



## Non-Adherent Dressings – Moist

### Xeroform™ Occlusive Gauze Dressing

- 3% Bismuth Tribromophenate in a special petrolatum blend on fine mesh gauze
- Non-adherent
- As a initial layer for wounds with light exudate where protection from contamination and/or deodorization is desire
- Sterile

Code	Size	Type	Qty Box	Qty Case
8884431302	2.5 cm x 20.3 cm	Strips in overwrap peelable foil	50	200
8884431605	12.7 cm x 22.9 cm	Strips in overwrap peelable foil	50	200
8884432000	10.2 cm x 2.7 m	Roll peelable foil packs	6	36
8884433301	2.5 cm x 20.3 cm	Patch peelable foil packs	50	200
8884433400	5.1 cm x 5.1 cm	Patch peelable foil packs	25	150
8884433500	10.2 cm x 10.2 cm	Patch peelable foil packs	25	150
8884433605	12.7 cm x 22.9 cm	Strip peelable foil packs	50	200



### Xeroform™ Non-Occlusive Gauze Dressing

- 3% Bismuth Tribromophenate in a water in-oil emulsion on fine mesh gauze
- For light exudate
- Clings and conforms
- Non-adherent construction
- Hydrophilic ointment
- As a initial layer for wounds with light exudate where protection from contamination and/or deodorization is desired
- Sterile

Code	Size	Type	Qty Box	Qty Case
8884436400	5.1 cm x 5.1 cm	Patch peelable foil pack	50	300
8884437100	12.7 cm x 22.9 cm	Strips in overwrap peelable foil	12	72
<b>STRIP OVERWRAP FOIL PACK</b>				
8884437201	2.5 cm x 20.3 cm	Strip peelable foil pack	50	200
8884437305	10.2 cm x 2.7 m	Roll	6	36



## Non-Adherent Dressings – Moist

### Vaseline Petrolatum Gauze

- Fine mesh, absorbent gauze impregnated with white petrolatum
- Non-adherent
- Petrolatum impregnant remains non-sensitizing, and non-irritating
- Sterile

Code	Size	Type	Qty Box	Qty Case
8884412600	2.5 cm x 91.4 cm	Strip in peelable foil packs	12	72
8884413605	7.6 cm x 22.9 cm	Strip in peelable foil packs	50	200
8884415600	7.6 cm x 91.4 cm	Strip in peelable foil packs	12	72
8884416600	15.2 cm x 91.4 cm	Strip in peelable foil packs	12	72
8884422600	2.5 cm x 91.4 cm	Strip in overwrap peelable foil packs	12	72
8884423600	7.6 cm x 22.9 cm	Strip in overwrap peelable foil packs	12	72
8884425600	7.6 cm x 91.4 cm	Strip in overwrap peelable foil packs	12	72
8884426600	15.2 cm x 91.4 cm	Strip in overwrap peelable foil packs	12	72



## Non-Adherent Dressings

### Telfa™ Non-Adherent Absorbent Dressing

- Non-adherent Telfa™ dressing won't disrupt the healing tissue
- Dressing for burns, surgical incisions, skin grafts, abrasions and other superficial wounds
- Sterile

Code	Description	Qty Box	Qty Case
<b>TELFA™ NON-ADHERENT DRESSINGS – STERILE 1'S IN PEEL-BACK PACKAGE</b>			
1055	5.1 cm x 7.6 cm	50	1200
1065	7.6 cm x 10.2 cm dressing	50	1200
1075	7.6 cm x 15.2 cm dressing	50	1200
1085	7.6 cm x 20.3 cm dressing	100	1200



## Non-Adherent Dressings

### Telfa™ Non-Adherent Clear Wound Dressings

- Precut sterile sizes
- Non-absorbent dressing
- Allows exudate to pass through apertures while limiting wound adherence to the secondary dressing
- Can be used with a variety of ointments

Code	Description	Qty Box	Qty Case
<b>STERILE 1'S IN PEEL-BACK PACKAGE</b>			
1109	7.6 cm x 7.6 cm	25	50
1111	10.2 cm x 12.7 cm	25	50
1113	30.5 cm x 30.5 cm	25	50
1114	30.5 cm x 61 cm	25	50

## Non Woven Dressings

### Excilon™ Drain Sponges

- 3-layer non-woven composition of drain sponge to promote absorption and wicking action
- Sterile

Code	Description	Qty Pack	Qty Box	Qty Case
7086	10.2 cm X 10.2 cm	2 pr	25 pr	600 (12 Box/25 pr each)

## Gauzes and Sponges

### Curity™ Gauze Sponges

#### 100% COTTON USP TYPE VII

Code	Description	Qty Tray	Qty Case
<b>1/PACK, STERILE</b>			
6132	7.6 cm x 7.6 cm – 12 ply	100	2400
6309	10.2 cm x 10.2 cm – 12 ply	100	1200
<b>NO STERILE, BULK</b>			
2346	7.6 cm x 7.6 cm – 12 ply	200	4000

## Bandages



### Kerlix™ Rolls – Sterile

- As primary dressing, open-weave design provides wicking action, aeration and absorbency

Code	Description	Qty Case
6725	Bandage roll, medium 8.6 cm x 3.3 m – 6 ply	96
6716	Bandage roll, large 11.4 cm x 2.8 m – 8 Ply	100
7880	Gauze roll, large 11.4 cm x 3.7 m – 8 Ply	100
7884	Super sponge – medium (15.2 cm x 17.1 cm)	480



### Webril™ Cotton Undercast Padding

- Crimped finish provides extra loft for cushioning

Code	Description	Qty Bag	Qty Case
1418	5 cm x 3.7 m roll	24	72
2059	7.6 cm x 3.7 m roll	12	72
3175	10.1 cm x 3.7 m	12	72
3489	15.2 cm x 3.7 m roll	6	36
2283	5.1 cm x 3.7 m roll – sterile 1 per pouch		50
2394	7.6 cm x 3.7 m roll – sterile 1 per pouch		50
2502	10.1 cm x 3.7 m – sterile 1 per pouch		50
2554	15.2 cm x 3.7 m roll – sterile 1 per pouch		24



### Webril™ Cotton Undercast Padding – Crimped Finish

- Crimped finish provides extra loft for cushioning

Code	Description	Qty Bag	Qty Case
2666	5.1 cm x 3.7 m roll – sterile 1 per pouch		50
2754	7.6 cm x 3.7 m roll – sterile 1 per pouch		50
2847	10.1 cm x 3.7 m – sterile 1 per pouch		50
2944	15.2 cm x 3.7 m roll – sterile 1 per pouch		24





## Ancillary Products

### Webcol™ Skin Cleansing Alcohol Preps

- Sterile, nonwoven sponges
- 70% isopropyl alcohol

Code	Description	Qty Box	Qty Case
6818-1	Medium 2 Ply	200	4000

# INDEX

## CARDIOTHORACIC

### SENTINEL SEAL™ CHEST DRAINAGE UNITS

#### Four Chamber, Dry Chest Drainage Unit Page 5

Code	Description	Collection Chamber	Sales Unit
8888571562	Sentinel Seal™ Chest Drainage Unit	2500 ml	5
8888571513	Sentinel Seal™ Chest Drainage Unit with dual collection chamber	1370 ml and 1030 ml	5
8888571489	Sentinel Seal™ Chest Drainage Unit with tube clamp & connector	2500 ml	5

### ALTITUDE™ CHEST DRAINAGE UNITS

#### Altitude™ Three Chamber, Dry Chest Drainage Unit Page 5

Code	Description	Collection Chamber	Sales Unit
8888571370	Altitude™ Chest Drainage Unit	2500 ml	5
8888571371	Altitude™ Chest Drainage Unit with dual collection chamber	1370 ml and 1030 ml	5

### AQUA-SEAL™ CHEST DRAINAGE UNITS

#### Three Chamber System Page 5

Code	Description	Collection Chamber	Sales Unit
8888571299	Aqua-Seal™ Chest Drainage Unit	2300 ml	5
8888571406	Aqua-Seal™ Chest Drainage Unit with dual collection chamber	1265 ml and 1050 ml	5
8888571315	Aqua-Seal™ Wet Suction Easy Change Connector	2300 ml	5

### THORA-SEAL™ CHEST DRAINAGE UNITS

#### Three Chamber System Page 6

Code	Description	Collection Chamber	Sales Unit
8884713308	Thora-Seal™ three chamber system	2500 ml	4
8884713900	Thora-Seal™ replacement collection chamber	2500 ml	6

#### Two Chamber System Page 6

Code	Description	Collection Chamber	Sales Unit
1180571570	Thora-Seal™ two chamber system	2600 ml	5

#### One Chamber System Page 6

Code	Description	Collection Chamber	Sales Unit
8884713100	Thora-Seal™ complete system	2000 ml	6
1814713105	Thora-Seal™ bottle only	2000 ml	6

### THORACIC CATHETERS AND TROCARS

#### Turkel™ Safety Thoracentesis and Paracentesis System Page 7

Code	Description	Sales Unit
8888566067	Turkel™ Safety Thoracentesis and Paracentesis System 8 Fr/Ch Safety Needle and Catheter – 12.5 cm 50 ml Monoject™ Luer Lock Syringe, 2000 ml Drainage Bag	10

#### Argyle™ Thoracic Catheters Page 7

Code	Description	Diameter Fr/Ch	N. of Bores	Sales Unit
<b>THORACIC CATHETERS</b>				
8888570507	Straight in PVC	12	4	10
8888570515	Straight in PVC	16	4	10
8888570523	Straight in PVC	20	5	10
8888570531	Straight in PVC	24	6	10
8888570549	Straight in PVC	28	6	10
8888570556	Straight in PVC	32	6	10
8888570564	Straight in PVC	36	6	10
8888570572	Straight in PVC	40	6	10
8888571018	Angled in PVC	16	4	10
8888571026	Angled in PVC	20	5	10
8888571034	Angled in PVC	24	6	10
8888571042	Angled in PVC	28	6	10
8888571059	Angled in PVC	32	6	10
8888571067	Angled in PVC	36	6	10
<b>Y CONNECTOR FOR CATHETERS</b>				
8884714200	"Y" extension for two thoracic catheters			10

#### Argyle™ Trocar Catheters Page 8

Code	Description	Diameter Fr/Ch	Length	Sales Unit
<b>TROCAR CATHETERS</b>				
8888560805	Sharp tip	8	23 cm	10
8888561019	Sharp tip	10	23 cm	10
8888561027	Sharp tip	12	23 cm	10
8888561035	Sharp tip	16	25 cm	10
8888561043	Sharp tip	20	40 cm	10
8888561050	Sharp tip	24	40 cm	10
8888561068	Sharp tip	28	40 cm	10
8888561076	Sharp tip	32	40 cm	10

## COMPRESSION

### KENDALL SCD™ COMFORT SLEEVES & 700 SERIES CONTROLLER

#### Kendall SCD™ Sequential Compression Comfort Sleeves – Thigh Length Page 10

Code	Description	Thigh Circumference	Size	Sales Unit
74010	Thigh length	Up to 40.6 cm	X-Small	5 pairs
74011	Thigh length	40.6 to 55.9 cm	Small	5 pairs
74012	Thigh length	55.9 to 71.1 cm	Medium	5 pairs
74013	Thigh length	71.1 to 91.4 cm	Large	3 pairs

#### Kendall SCD™ Sequential Compression Comfort Sleeves – Knee Length Page 10

Code	Description	Calf Circumference	Size	Sales Unit
74021	Knee length	Up to 35.5 cm	Small	5 pairs
74022	Knee length	35.5 to 53.3 cm	Medium	5 pairs
74023	Knee length	53.3 to 66.0 cm	Large	5 pairs

## Kendall SCD™ Sequential Compression Tear-Away – Thigh Length Page 10 ☹

Code	Description	Thigh Circumference	Size	Sales Unit
74041	Thigh length	40.6 to 55.9 cm	Small	5 pairs
74042	Thigh length	55.9 to 71.1 cm	Medium	5 pairs
74043	Thigh length	71.1 to 91.4 cm	Large	3 pairs

## Kendall SCD™ 700 Series Controller Page 10 ☹

Code	Description	Sales Unit
295250	SCD™ 700 series controller – EU	1
295251	SCD™ 700 series controller – UK	1

## KENDALL SCD™ EXPRESS SLEEVES

### Kendall SCD™ Express Sleeves Page 10 ☹

Code	Description	Calf/Thigh Circumference	Size	Sales Unit
73022	Knee length	34.3 to 53.3 cm	Medium	5 pairs
73023	Knee length	53.3 to 66.0 cm	Large	5 pairs
9790	Knee length	55.9 to 81.2 cm	X-Large	5 pairs
73011	Thigh length	31.1 to 55.9 cm	Small	5 pairs
73012	Thigh length	55.9 to 71.1 cm	Medium	5 pairs
73013	Thigh length	71.1 to 91.4 cm	Large	3 pairs

### Kendall SCD™ Express Tear-Away Sleeves Page 11 ☹

Code	Description	Thigh Circumference	Size	Sales Unit
73041	Thigh length	31.1 to 55.9 cm	Small	5 pairs
73042	Thigh length	55.9 to 71.1 cm	Medium	5 pairs
73043	Thigh length	71.1 to 91.4 cm	Large	3 pairs

### Kendall SCD™ Express Sterile Sleeve Page 11 ☹

Code	Description	Thigh Circumference	Size	Sales Unit
9736	Thigh length	55.9 to 71.1 cm	Medium	5 singles

### Kendall SCD™ Express Foot Cuff Page 11 ☹

Code	Description	Size	Sales Unit
73032	Foot cuff	Regular	10 singles
73033	Foot cuff	Large	10 singles

### Kendall SCD™ Express Tubing Page 11 ☹

Code	Description	Sales Unit
9528	Compression system tubing set 2.1 m	1
9595	Compression system tubing set (ambulatory extension) 1.2 m	1

## A-V IMPULSE™ SYSTEM

### A-V Impulse™ 6000 Series Controller Page 11 ☹

Code	Description	Sales Unit
AV6000-EUR	Dual channel controller	1 each
AV6000-UK	Dual Channel controller UK version	1 each

### ImPad™ Rigid Sole Foot Cover Page 11 ☹

Code	Description	Sales Unit
AV740-5	Medium	5 pairs
AV750-5	Large	5 pairs

### A-V Impulse™ Tubing Assembly Page 11 ☹

Code	Description	Sales Unit
AV810-01	Air hose blue, 3 m	1 each
AV820-01	Air hose red, 3 m	1 each
AV830-00	Air hose grey, 3 m	1 each

## T.E.D.™ ANTI-EMBOLISM STOCKINGS

### Thigh Length Style Page 12 ☹

Thigh Circumference	Calf Girth	Code	Leg Length	Size	Length Colour	Sales Unit
Less than 63.5 cm	Small	3071LF	Short	A	Yellow	6 pairs
		3130LF	Regular	B	White	6 pairs
		3222LF	Long	C	Yellow	6 pairs
	Medium	3310LF	Short	D	White	6 pairs
		3416LF	Regular	E	Yellow	6 pairs
		3549LF	Long	F	White	6 pairs
	Large	3634LF	Short	G	Yellow	6 pairs
		3728LF	Regular	H	White	6 pairs
		3856LF	Long	J	Yellow	6 pairs
63.5 to 81.3 cm	X-Large	4010LF	Short	K	White	6 pairs
		4114LF	Regular	L	Yellow	6 pairs
		4216LF	Long	M	White	6 pairs
		3180LF	Short	N	Yellow	6 pairs
		3181LF	Regular	P	White	6 pairs
		3182LF	Long	Q	Yellow	6 pairs
81.3 to 91.4 cm	XX-Large	3183LF	Short	R	White	6 pairs
		3184LF	Regular	S	Yellow	6 pairs
		3185LF	Long	T	White	6 pairs

### Thigh Length with Belt Style Page 12 ☹

Thigh Circumference	Calf Girth	Code	Leg Length	Size	Length Colour	Sales Unit	
Less than 63.5 cm	X-Small	3306	Regular	AA+	White	6 pairs	
		3320	Long	BB+	Blue	6 pairs	
	Small	3039	Regular	A+	White	6 pairs	
		3364	Long	B+	Blue	6 pairs	
	Medium	3144	Regular	C+	White	6 pairs	
		3449	Long	D+	Blue	6 pairs	
	Large	3221	Regular	E+	White	6 pairs	
		3523	Long	F+	Blue	6 pairs	
	63.5 to 81.3 cm	X-Large	3922	Regular	G+	White	6 pairs
			3995	Long	H+	Blue	6 pairs

## T.E.D.™ ANTI-EMBOLISM STOCKINGS

### Knee Length Style Page 12 ☹

Thigh Circumference	Calf Girth	Code	Leg Length	Size	Length Colour	Sales Unit
Less than 30.5 cm	Small	7071	Regular	A-	White	12 pairs
		7339	Long	B-	Blue	12 pairs
30.5 to 38 cm	Medium	7115	Regular	C-	White	12 pairs
		7480	Long	D-	Blue	12 pairs
38 to 44.5 cm	Large	7203	Regular	E-	White	12 pairs
		7594	Long	F-	Blue	12 pairs
44.5 to 51 cm	X-Large	7604	Regular	G-	White	12 pairs
		7802	Long	H-	Blue	12 pairs
51 to 58.4 cm	XX-Large	7470LF	Regular	J-	White	12 pairs
		7471LF	Long	K-	Blue	12 pairs
58.4 to 66 cm	XXX-Large	7472LF	Regular	L-	White	12 pairs
		7473LF	Long	M-	Blue	12 pairs

## T.E.D.™ ANTI-EMBOLISM STOCKINGS FOR CONTINUING CARE

### Knee Length Style Page 13 €€

Calf Circumference	Length	Description	Item Code			Sales Unit
			White	Beige	Black	
<30.5 cm	<41 cm	Small Regular	4278-	4265-	4434	1 pair
	>41 cm	Small Long	4282	4317-	4572-	1 pair
30.5 to 38 cm	<43 cm	Medium Regular	4279-	4271	4435	1 pair
	>43 cm	Medium Long	4283	4323-	4573	1 pair
38 to 44.5 cm	<46 cm	Large Regular	4280	4289	4436	1 pair
	>46 cm	Large Long	4284	4335	4574	1 pair
44.5 to 51 cm	<46 cm	X-Large Regular	4281	4296-	4437-	1 pair
	>46 cm	X-Large Long	4285	4344	4575	1 pair

### Thigh Length Style Page 13 €€

Thigh Circumference	Calf Circumference		Length			Sales Unit
	Small	Medium	Short <74 cm	Regular 74-84 cm	Long >84 cm	
<63.5 cm	Small	30.5 cm	4297-	4300-	4303	1 pair
	Medium	30.5 to 38 cm	4298-	4301-	4304	1 pair
	Large	38 to 44.5 cm	4299	4302	4305-	1 pair

Available in white only

## DENTAL

### ENDODONTIC NEEDLES

#### Monoject™ Endodontic Needles – Sterile Page 15 €€<sub>0123</sub>

Code	Needle Gauge and Length	Colour	Sales Unit
8881471232	23G x 32 mm	Orange	4 x (25 x 1 pc)
8881471273	27G x 32 mm	Yellow	4 x (25 x 1 pc)

#### Monoject™ Endodontic Irrigation Syringes with Needle – Sterile

Page 15 €€<sub>0123</sub>

Code	Size Syringe	Needle Gauge and Length	Sales Unit
1180323114E	3 ml	23G x 32 mm (1.25 inches)	800 (8x100)
1180327114E	3 ml	27G x 32 mm (1.25 inches)	800 (8x100)

## ELECTRODES

### KENDALL™ ECG ELECTRODES FOR ADULTS

#### General monitoring Electrodes - Long Term monitoring Page 17 €€

Code	Match Code	Description	Sales Unit
31.1946.21	H345G	Radiolucent, suitable for X-ray, CT and MRI with solid gel, strong adhesive (72 hour wearable, shower and wash resistant), with integrated abradar for skin preparation, 50 x 45 mm Ø, with stud adaptor	30/300
31.1925.21	H925G	With solid gel for long-term use, 57 x 34 mm, with stud adaptor	50/500
31.1355.21	H1355G	With solid gel for long-term use, also suitable for neurology and paediatric applications, 43 x 35 mm Ø, stud adaptor	50/500
31.1355.821	H1355TSG	As H1355G, but radiolucent, suitable for CT and MRI	50/500
31.1245.21	H1245G	With special hydrogel for long-term use, also suitable for neurology and paediatric applications, 30 x 24 mm Ø, stud adaptor	50/500
31.1917.21	H9155G	As H925G, but equipped with strong adhesive (recommended for long-term ECG monitoring and stress test), 57 x 34 mm, with stud adaptor	50/500
31.1915.821	H91TSG	As H9155G, but radiolucent, suitable for X-ray, CT and MRI, 57 x 34 mm, with stud adaptor	50/500

31.1935.21	H935G	As H925G, however only 42 x 24 mm	50/500
31439725	Medi-Trace™ 133	With conductive adhesive hydrogel, packed in perforated strips of 3, 30 mm Ø with stud adaptor	3/600
31050522	Medi-Trace™ 200	With conductive adhesive hydrogel, packed in perforated strips of 5, 36,5 mm tear drop with stud adaptor	100/1000
31118733	Medi-Trace™ 100	With conductive adhesive hydrogel, Packed in perforated strips of 5, 30 mm Ø with stud adaptor	100/1000
31013926	Medi-Trace™ 530	With conductive adhesive hydrogel, pack of 30 loose electrodes, 45 mm tear drop with stud adaptor	30/600
31.1914.21	H9145G	With solid gel for long term use, 57 x 34 mm, with 4 mm banana adaptor	30/300
31.1914.821	H914TSG	As H9145G but radiolucent, suitable for CT and MRI, with stud adaptor	30/300
RR3305		With conductive adhesive hydrogel, radiolucent, perforated strips of 5, 31.8x36.8 mm with stud adaptor, repositionable 2 times within 72 hours	5/600
ES40076	Medi-Trace™ H59P	Cloth backing with conductive adhesive hydrogel, repositionable, 38 x 22 mm with stud adaptor	30/450
ES40281	Medi-Trace™ H69P	Cloth backing with conductive adhesive hydrogel, repositionable, 22 mm squared with stud adaptor	30/450
22750	Medi-Trace™ 750	Clear Tape Electrodes, Conductive Adhesive Hydrogel, Radiolucent	50/1000

#### General monitoring Electrodes - Short Term monitoring Page 17 €€

Code	Match Code	Description	Sales Unit
31.1986.21	H985G	Microporous non-woven backing, with solid gel, 55 mm Ø, with stud adaptor	30/300
31.1996.21	H995G	Non-woven backing, with solid gel, 55 mm Ø, with stud adaptor	30/300

### KENDALL™ ECG ELECTRODES FOR ADULTS

#### Diagnostic Monitoring Electrodes Page 18 €€

Code	Match Code	Description	Sales Unit
31.1985.21	H98LG	Microporous non-woven backing, with liquid gel, 55 mm Ø, with stud adaptor	30/300
31.1345.21	H34LG	With liquid gel, wearable for 72 hours, shower and wash resistant, particularly suitable for diagnostic and stress test, 50 x45 mm Ø, with stud adaptor	30/300
31.1663.21	H66LG	With liquid gel, 55 mm Ø, with stud adaptor	30/300
31.1995.21	H99LG	Non-woven backing, with liquid gel, 55 mm Ø, with stud adaptor	30/300

#### Disposable ECG TAB Electrodes Page 18 €€

Code	Match Code	Description	Sales Unit
31433538	Kendall™ 5400	Tab Electrodes 23 x 33 mm	100/4000
30807732	Kendall™ 5500	Tab Electrodes 23 x 33 mm	100/2000

### ACCESSORIES Page 18 €€

Code	Match Code	Description	Sales Unit
22660A	A34	Universal adaptor for stud and TAB electrodes, adapts to 4 mm banana plug	10

### KENDALL™ NEONATAL ECG ELECTRODES

#### Disposable ECG Electrodes – Transparent backing Page 19 €€

Code	Match Code/ Lead Length	Colour	Connector Type	Shape	Size (mm)	Special Features	Sales Unit
31424780	1041PTS/ 60 cm	AAMI (white, black, green)	1.5 mm DIN safety plug	Round	25 Ø	Suitable for incubator, X-ray radiolucent, repositionable	3/300
31424785	1042PTS/ 60 cm	AAMI (white, black, green)	1.5 mm DIN safety plug	Rectangular	25 x 30	Suitable for incubator, X-ray radiolucent, repositionable	3/300

## Disposable ECG Electrodes – Cloth backing Page 19 <<

Code	Match Code/ Lead Length	Colour Code	Connector Type	Shape	Size (mm)	Special Features	Sales Unit
31424768	1050NPSM/ 60 cm	AAMI (white, black, green)	1.5 mm, DIN safety plug	Round	30 Ø	Suitable for incubator, repositionable	3/300
31424743	1051NPSM/ 60 cm	AAMI (white, black, green)	1.5 mm, DIN safety plug	Round	25 Ø	Suitable for incubator, repositionable, radiolucent	3/300
31424735	1052NPSM/ 60 cm	AAMI (white, black, green)	1.5 mm, DIN safety plug	Round	28 Ø	Suitable for incubator, repositionable	3/300
31424719	1053NPSM/ 60 cm	AAMI (white, black, green)	1.5 mm, DIN safety plug	Round	28 Ø	Suitable for incubator, repositionable, radiolucent	3/300
31424768	1050NPSM/ 60 cm	AAMI (white, black, green)	1.5 mm, DIN safety plug	Round	0 Ø	Suitable for incubator, repositionable	3/300

## Disposable ECG Electrodes – Foam backing Page 19 <<

Code	Match Code/ Lead Length	Colour Code	Connector Type	Shape	Size (mm)	Special Features	Sales Unit
30710001	4103/ 45.7 cm	AAMI (white, black, green)	1.5 mm, DIN safety plug	Round	23 Ø	Suitable for incubator, repositionable	3/300
31368072	4203/ 45.7 cm	IEC (red, yellow, black)	1.5 mm, DIN safety plug	Round	23 Ø	Suitable for incubator, repositionable	3/300
31086575	4105/ 4.7 cm	AAMI (white, black, green)	1.5 mm, DIN safety plug	Oval	20 x 30	Suitable for incubator, repositionable	3/300

## KENDALL™ NEONATAL ECG ELECTRODES

### Disposable ECG Electrodes – non-woven tape backing Page 20 <<

Code	Match Code/ Lead Length	Colour Code	Connector Type	Shape	Size (mm)	Special Features	Sales Unit
31.5872.04	H87PG/F/ 60 cm	IEC (red, yellow, black)	1.5 mm DIN safety plug	Round	30 Ø	Suitable for incubator, repositionable	3/120
31.5072.03	H207PG/F/ 60 cm	IEC (red, yellow, black)	1.5 mm DIN safety plug	Round	22 x 30	Suitable for incubator, repositionable	3/120
31.5032.26	H203PG/ 60 cm	White	4 mm banana socket	Oval	22 x 30	Suitable for incubator, repositionable	15/120
31.5033.26	H203PGK/ 60 cm	White	4 mm banana, socket, extra short lead	Oval	22 x 30	Suitable for incubator, repositionable	15/120
31.5072824	H207PGT/ 90 cm	Carbon leads, transparent leads jacketing	1.5 mm DIN safety plug	Oval	22 x 30	Suitable for incubator, radiolucent, repositionable	3/90
31.5039124	H203PGT/ 90 cm	Carbon leads, transparent leads jacketing	4 mm banana socket	Oval	22 x 30	Suitable for incubator, radiolucent, repositionable	3/90
31.3876.04	H87H/F/ 60 cm	IEC (red, yellow, black)	1.5 mm DIN safety plug	Round	30 Ø	Suitable for standard care, repositionable	3/120
31.3076.26	H207H/ 60 cm	White	1.5 mm DIN safety plug	Oval	22 x 30	Suitable for standard care, repositionable	15/120
31.3276.03	H27HF/ 60 cm	IEC (red, yellow, black)	1.5 mm DIN safety plug	Round	24 Ø	Suitable for standard care, repositionable	3/120

## KENDALL™ INFANT ECG ELECTRODES

### Infant ECG Electrode for Prematures & Speciality Application Page 20 <<

Code	Match Code/ Lead Length	Colour Code	Connector Type	Shape	Size (mm)	Special Features	Sales Unit
EP30016	LB36-LB/ 45.7 cm	AAMI (white, black, green)	1.5 mm, DIN safety plug	Rectangular	17 x 22	Suitable for incubator, with integrated fastening band, repositionable	3/300

## ELECTRODES FOR POLYSOMNOGRAPHY

### Page 20 <<

Code	Match Code	Description	Sales Unit
31.6071.26	H207LG/F60	Electrode for PSG applications with integrated 60 cm lead wire red, yellow, black	12/180
31.6072.26	H207LGN	Electrode for PSG applications with integrated 150 cm lead wire white	10/140
31.6073.26	H207LG/F150	Electrode for PSG applications with integrated 150 cm lead wire red, yellow, black	12/180
31.6074.26	H207LG/F200	Electrode for PSG applications with integrated 150 cm lead wire red, yellow, black	6/90

## MULTIFUNCTION DEFIBRILLATION ELECTRODES

### Medi-Trace™ Classic Page 21 <<

Code	Match Code	Can be used with compatible defibrillator modes	Quantity/ Pouch	Pairs/ Case
31319281	1310P	LIFEPAK (LP 9, LP 12, LP 15, LP 20, LP 20E, LP 500, LP 1000), LIFEPAK CR Plus, LIFEPAK EXPRESS	1 Pair	10
31177721	1210H	CodeMaster XL+ (M1722A, M1722B, M1723A, M1723B), Heartstream XLT M3500B <sup>1</sup> , HeartStart MRx (M3535A <sup>1</sup> , M3536A <sup>1</sup> ), HeartStart XL M4735A <sup>1</sup> , HeartStart 4000 M5500B <sup>1</sup>	1 Pair	10
31469219	1410Z	Zoll PD1400, PD2000, R Series BLS, R Series Plus, R Series ALS, X Series, Propaq MD	1 Pair	10
40000006	1710H	Heartstream XLT M3500B <sup>2</sup> , HeartStart MRx (M3535A <sup>2</sup> , M3536A <sup>2</sup> ), HeartStart XL M4735A <sup>2</sup> , HeartStart 4000 M5500B <sup>2</sup> , HeartStart FR2+ (M3840A, M3841A, M3860A, M3861A <sup>2</sup> )	1 Pair	10

1 Requires Philips M3508A hands-free cable

2 AED self-test will fail if electrodes are pre-connected

3 Requires Philips M3507A (barrel style) hands-free cable

## MULTIFUNCTION DEFIBRILLATION ELECTRODES

### Medi-Trace™ Cadence PC (pre-connect) Page 21 <<

Code	Can be used with compatible defibrillator modes	Quantity/ Pouch	Pairs/ Case
22550PC	LIFEPAK (LP 9, LP 12, LP 15, LP 20, LP 20E, LP 500, LP 1000), LIFEPAK CR Plus, LIFEPAK EXPRESS	1 Pair	10
22660PC	Heartstream XLT M3500B <sup>1</sup> , HeartStart MRx (M3535A, M3536A <sup>1</sup> ), HeartStart XL M4735A, HeartStart 4000 M5500B, HeartStart FR2+ (M3840A, M3841A, M3860A, M3861A <sup>2</sup> )	1 Pair	10
22661PC	CodeMaster XL+ (M1722A, M1722B, M1723A, M1723B), Heartstream XLT M3500B <sup>3</sup> , HeartStart MRx (M3535A, M3536A <sup>3</sup> ), HeartStart XL M4735A <sup>3</sup> , HeartStart 4000 M5500B <sup>3</sup>	1 Pair	5
22770PC	Zoll PD1400, PD2000, R Series BLS, R Series Plus, R Series ALS, X Series, Propaq MD	1 Pair	10

All above codes are suitable for patients ≥ 8 years or ≥ 25 kg

1 Requires Philips M3508A hands-free cable

2 AED self-test will fail if electrodes are pre-connected

3 Requires Philips M3507A (barrel style) hands-free cable

### Medi-Trace™ Cadence RTS (radio-transparent) Page 21 <<

Code	Can be used with compatible defibrillator modes	Quantity/ Pouch	Pairs/ Case
22550A	LIFEPAK (LP 9, LP 12, LP 15, LP 20, LP 20E, LP 500, LP 1000), LIFEPAK CR Plus, LIFEPAK EXPRESS	1 Pair	10
22550R	LIFEPAK (LP 9, LP 12, LP 15, LP 20, LP 20E, LP 500, LP 1000), LIFEPAK CR Plus, LIFEPAK EXPRESS	1 Pair	10
22550P*	LIFEPAK (LP 9, LP 12, LP 15, LP 20, LP 20E, LP 500, LP 1000)	1 Pair	5
22770R	Zoll PD1400, PD2000, R Series, R Series Plus, R Series ALS, Propaq MD	1 Pair	10
22770P*	Zoll PD1400, PD2000, R Series BLS, R Series Plus, R Series ALS, X Series, Propaq MD	1 Pair 1 Pair	5 5
22660R	Heartstream XLT M3500B <sup>1</sup> , HeartStart MRx (M3535A, M3536A), HeartStart XL M4735A <sup>1</sup> , HeartStart 4000 M5500B <sup>1</sup> , HeartStart FR2+ (M3840A, M3841A, M3860A, M3861A <sup>2</sup> )	1 Pair	10

22660P*	Heartstream XLT M3500B, HeartStart MRx (M3535A, M3536A, M4735A), HeartStart 4000 M5500B	1 Pair	10
22661R	CodeMaster XL+ (M1722A, M1722B, M1723A, M1723B), Heartstream XLT M3500B <sup>3</sup> , HeartStart MRx (M3535A, M3536A <sup>3</sup> ), HeartStart XL M4735A <sup>3</sup> , HeartStart 4000 M5500B <sup>3</sup>	1 Pair	10
22661P	CodeMaster XL+ (M1722A, M1722B, M1723A, M1723B), Heartstream XLT M3500B <sup>3</sup> , HeartStart MRx (M3535A, M3536A <sup>3</sup> ), HeartStart XL M4735A <sup>3</sup> , HeartStart 4000 M5500B <sup>3</sup>	1 Pair	5

\* These pads are designed for patients < 8 years or < 25 kg.

1 Requires Philips M3508A hands-free cable

2 AED self-test will fail if electrodes are pre-connected

3 Requires Philips M3507A (barrel style) hands-free cable

## KENDALL™ DL – DISPOSABLE LEAD-WIRES

### Kendall™ DL Bedside Lead Wire Syst Page 22 €

IEC	AAMI	Description	Quantity (each)
33103EG	–	3 Lead bedside red, green, yellow	10
33103EB	–	3 Lead bedside black, red, yellow	10
33105E	33105	5 Lead bedside	10
33110E	33110	10 Lead bedside	5
33111E	33111	Chest Leads (V-leads)	10

### Kendall™ DL Telemetry Lead Sets Direct Connect Page 22-23 €

IEC	AAMI	Description	Quantity (each)
33112E	33112	For use with 5 Lead Telemetry System, manufactured by GE Medical Systems Information Technologies, Inc., under its trademark ApexPro™	10
33113E	33113	For use with 5 Lead Telemetry System manufactured by Koninklijke Philips Electronics N.V. under its trademark Veridia™	10
33114E	33114	For use with 5 Lead Telemetry Systems manufactured by Koninklijke Philips Electronics N.V. under its trademark IntelliVue™	10
–	33125	For use with 6 Lead CH Telemetry System manufactured by GE Medical Systems Information Technologies, Inc., under its trademark ApexPro™	10
–	33126	For use with 6 lead FH Telemetry Systems manufactured by GE Medical Systems Information Technologies, Inc., under its trademark ApexPro™	10
33132E	–	For use with 5 Lead Draeger Medical, Inc., Dual Pin Cables	10
–	33133	For use with 5 Lead Spacelabs Healthcare Telemetry Systems	10
–	33134	For use with 6 Lead Nihon Kohden Telemetry Systems	10
–	33147	5 lead non-SPO2. For use with Philips IntelliVue™ MX40 Telemetry System	100
–	33148	6 lead non-SPO2. For use with Philips IntelliVue™ MX40 Telemetry System – Direct Connect	100
–	33150	6 lead SPO2. For use with Philips IntelliVue™ MX40 Telemetry System, compatible with Nellcor SpO <sub>2</sub> Sensors – Direct Connect	100
–	33151	Disposable 5 Lead. For use with Philips IntelliVue™ MX40 Telemetry System, compatible with Philips SpO <sub>2</sub> Sensors – Direct Connect	100
–	33152	Disposable 6 Lead. For use with Philips IntelliVue™ MX40 Telemetry System, compatible with Philips SpO <sub>2</sub> Sensors – Direct Connect	100

### Kendall™ DL Dual Connect System Page 23 €

IEC	AAMI	Description	Quantity (each)
33135E	–	Dual connect Lead Set + Cable - 5 lead	5
33136E	–	Dual connect Lead Set only - 5 lead	25
–	33145	Dual connect Cable only	5

### Radiolucent Disposable Lead Sets Page 23 €

IEC	AAMI	Description	Quantity (each)
–	33136R36	Radiolucent Lead Set, 5 Lead 90 cm	10
–	33136R72	Radiolucent Lead Set, 5 Lead 180 cm	10
–	33136RV36	Radiolucent Chest Lead Set, 5 Lead 90 cm	10

## NUTRITIONAL INSUFFICIENCY

### KANGAROO™ EPUMP ENTERAL FEEDING PUMP Page 26 €

IEC	Description	Sales Unit
482400J	Kangaroo™ ePump including pole clamp, power adaptor and operating manual	1

### Kangaroo™ ePump Enteral Pump Feeding Sets, DEHP-Free PVC, Sterile Page 26 €

#### Feeding sets with spikes

IEC	Description	Sales Unit
<b>FEED SETS ONLY</b>		
777401	ENPlus spike set	30
777402	Dual ENPlus spike set	30
777403	ENPlus 3-in-1 set	30
<b>FEED AND WATER SETS</b>		
777405	ENPlus spike & 1000ml flush set	30
777406	ENPlus 3-in-1 & 1000ml flush set	36
777407	ENPlus spike feed & vented spike flush set	30
777408	ENPlus 3-in-1 feed & 3-in-1 flush set	30

### Feeding sets with integrated bags Page 26 €

IEC	Description	Sales Unit
<b>FEED SETS ONLY</b>		
777064	500 ml set	30
777106	1000 ml set	30
777163	1600 ml set	30
<b>FEED AND WATER SETS</b>		
778104	1000 ml feed & 1000 ml flush set	30
778161	1600 ml feed & 1000 ml flush set	30
673662	1000 ml feed & 1000 ml flush set, nonsterile	30

### Adaptors Page 26 €

Code	Description	Sales Unit
777222	ENFit™ to ENLock Oral-luer medication port adaptor	30
777444	ENFit™ to funnel feeding tube adaptor	30
1814717995	3in1 adaptor	50

## KANGAROO™ EPUMP ACCESSORY PRODUCTS Page 27

Code	Description	Sales Unit
382491	Power adaptor €	1
382492	Pole clamp €	1
1054715*	Battery pack	1
382493	Plug attachment set €	1
382499	RS232 Serial Port communications cable €	1

\* This product is not a medical device and is not CE marked.

## KANGAROO™ JOEY ENTERAL FEEDING PUMP Page 27 €

Code	Description	Sales Unit
383400	Kangaroo™ Joey Enteral Feeding Pump including pole clamp, power adaptor and operating manual	1

### Kangaroo™ Joey Enteral Pump Feeding Sets, DEHP-Free PVC, Sterile Page 27 €

#### Feeding sets with spikes

Code	Description	Sales Unit
<b>FEED SETS ONLY</b>		
777501	ENPlus spike set	30
777502	Dual ENPlus spike set	30
777503	ENPlus 3-in-1 set	30
<b>FEED AND WATER SETS</b>		
777505	ENPlus spike & 1000ml flush set	30
777506	ENPlus 3-in-1 & 1000ml flush set	36



## Feeding sets with integrated bags

Code	Description	Sales Unit
<b>FEED ONLY SETS</b>		
666064	500 ml set	30
666106	1000 ml set	30
<b>FEED AND WATER SETS</b>		
668104	1000 ml feed & 1000 ml flush set	30

## Adaptors Page 28 <sup>CE 0123</sup>

Code	Description	Sales Unit
777222	ENFit™ to ENLock Oral-luer medication port adaptor	30
777444	ENFit™ to funnel feeding tube adaptor	30
1814717995	3in1 adaptor	50

## KANGAROO™ JOEY ENTERAL PUMP ACCESSORIES AND BACKPACKS Page 28

Code	Description	Sales Unit
383491	Power adaptor <sup>CE 0123</sup>	1
383493	Pole clamp <sup>CE</sup>	1
F010506	Battery pack*	1
382493	Plug attachment set <sup>CE 0123</sup>	1
770026	Mini backpack – blue <sup>CE</sup>	1
770028	Super mini backpack – blue <sup>CE</sup>	1
770034	Mini backpack – pink <sup>CE</sup>	1
770029	Super mini backpack – pink <sup>CE</sup>	1
770036	Large backpack – black <sup>CE</sup>	1

\*This product is not a medical device and is not CE marked.

## ENTERAL NUTRITION GRAVITY FEEDING SETS

### Enteral Nutrition Delivery Gravity Sets, DEHP-Free PVC, Sterile

Page 28 <sup>CE 0123</sup>

Code	Description	Sales Unit
1180262154	Paediatric 60ml feeding set	36
777901	ENPlus gravity feeding set	30

## PAEDIATRIC FEEDING TUBES

### Kangaroo™ Neonatal and Paediatric Feeding Tubes –

Polyurethane (PUR), Sterile Page 29 <sup>CE 0123</sup>

Code	Fr/Ch	Length	Sales Unit
461355E	3.5	31 cm with 1 cm depth markings	10
461412E	5	51 cm with 1 cm depth markings	10
461420E	6.5	51 cm with 1 cm depth markings	10
461438E	6.5	91 cm with 1 cm depth markings	10
461701E	8	51 cm with 1 cm depth markings	10
461503E	5	91 cm with 1 cm depth markings	10
461800E	8	107 cm with 1 cm depth markings	10
461042E	10	107 cm with 1 cm depth markings	10

### Argyle™ Neonatal and Paediatric Feeding Tubes –

DEHP-Free Polyvinylchloride (PVC), Sterile Page 29 <sup>CE 0123</sup>

Code	Fr/Ch	Length	Sales Unit
461206E	3.5	31 cm	10
460802E	5	41 cm	50
461008E	5	91 cm	50
460711E	6.5	41 cm	50
460604E	8	41 cm	50
460406E	8	107 cm	50
460208E	10	107 cm	50

## KANGAROO™ NASOGASTRIC FEEDING TUBES WITH ENFIT™ CONNECTION Page 29 <sup>CE 0123</sup>

Code	Description	Length	Weight	Stylet	Sales Unit
8884720817E	8 Fr/Ch	91 cm	3 g	NO	10
8884720825E	8 Fr/Ch	91 cm	3 g	YES	10
8884720841E	8 Fr/Ch	109 cm	3 g	NO	10
8884720858E	8 Fr/Ch	109 cm	3 g	YES	10
8884710859E	8 Fr/Ch	109 cm	7 g	YES	10
8884710867E	8 Fr/Ch	109 cm	5 g	YES	10
8884710826E	8 Fr/Ch	140 cm	7 g	YES	10
8884721088E	10 Fr/Ch	109 cm	5 g	YES	10
8884721096E	10 Fr/Ch	109 cm	5 g	NO	10
8884711006E	10 Fr/Ch	109 cm	7 g	YES	10
8884721055E	10 Fr/Ch	140 cm	5 g	YES	10
8884721211E	12 Fr/Ch	91 cm	5 g	NO	10
8884721237E	12 Fr/Ch	109 cm	5 g	NO	10
8884721252E	12 Fr/Ch	109 cm	5 g	NO	10
8884711246E	12 Fr/Ch	109 cm	7 g	NO	10
8884711253E	12 Fr/Ch	109 cm	7 g	YES	10
8884721255E	12 Fr/Ch	140 cm	5 g	YES	10

## SHORT TERM ACCESS – NASO-JEJUNAL FEEDING TUBES

### Kangaroo™ Naso-Jejunal Feeding Tube – Polyurethane (PUR),

Non-Sterile Page 30 <sup>CE 0123</sup>

Code	Fr/Ch	Length	Weight	Guide/wire	Sales Unit
8884752505E	12	152 cm	7 g	YES	2

### Kangaroo™ Naso-Jejunal Feeding and Gastric Decompression Tube –

Polyurethane (PUR), Non-Sterile Page 30 <sup>CE 0123</sup>

Code	Fr/Ch	Length	Weight	Guide/wire	Sales Unit
8884752513E	16	Gastric tube section	-	YES	2

## GASTROINTESTINAL ACCESS – DUAL LUMEN TUBES

### Salem Sump Stomach Tube – Silicone, Sterile Page 30 <sup>CE 0123</sup>

Code	Fr/Ch	Length	Sales Unit
8888265108	10	90 cm	10
8888265116	12	120 cm	10
8888265124	14	120 cm	10
8888265132	16	120 cm	10
8888265140	18	120 cm	10

## GASTROINTESTINAL ACCESS – DUAL LUMEN TUBES

### Salem Sump Stomach Tube with Anti-Reflux Valve – Silicone,

Sterile Page 31 <sup>CE 0123</sup>

Code	Fr/Ch	Length	Sales Unit
8888266221	12	120 cm	10
8888266247	14	120 cm	10
8888266262	16	120 cm	10
8888266288	18	120 cm	10



**Salem Sump Gastroduodenal Tube – Polurethane (PUR), Sterile**

Page 31 <sup>€</sup><sub>0123</sub>

Code	Fr/Ch	Length	Sales Unit
1180264523	14	120 cm	25
1180264531	16	120 cm	25
1180264549	18	120 cm	25
1180264556	20	120 cm	25

**Salem Sump Gastroduodenal Tube – DEHP-Free Polyvinylchloride (PVC), Sterile**

Page 31 <sup>€</sup><sub>0123</sub>

Code	Fr/Ch	Length	Sales Unit
8888264804	10	90 cm	25
8888264820	12	108 cm	25
8888264846	14	120 cm	25
8888264861	16	120 cm	25
8888264887	18	120 cm	25
1180264895	20	120 cm	25

**Salem Sump Gastroduodenal Tube with Moulded Connector – DEHP-Free Polyvinylchloride (PVC), Sterile**

Page 31 <sup>€</sup><sub>0123</sub>

Code	Fr/Ch	Length	Sales Unit
1180264408	10	90 cm	25
1180264416	12	120 cm	25
1180264424	14	120 cm	25
1180264432	16	120 cm	25
1180264440	18	120 cm	25
1180264457	20	120 cm	25

**GASTROINTESTINAL ACCESS – DUAL LUMEN TUBES**

**Salem Sump Stomach Tube with Anti-Reflux Valve – DEHP-Free Polyvinylchloride (PVC), Sterile**

Page 32 <sup>€</sup><sub>0123</sub>

Code	Fr/Ch	Length	Sales Unit
266122	14	120 cm	50
266130	16	120 cm	50
266148	18	120 cm	50

**Salem Sump Dual Lumen Stomach Tube With ENFit Connection**

<sup>€</sup><sub>0123</sub>

Code	Fr/Ch	Length	Sales Unit
8888268060E	6	61 cm	10
8888268086E	8	61 cm	10
8888264911E	10	91 cm	50
8888264929E	12	122 cm	50
8888264945E	14	122 cm	50
8888264960E	16	122 cm	50
8888264986E	18	122 cm	50
266197	Anti-Reflux Valve		10

**Salem Sump Anti-Reflux Valve, Sterile**

Page 33 <sup>€</sup><sub>0123</sub>

Code	Sales Unit
266197	10

**Kangaroo™ P.E.G. Accessory Products, Non-Sterile**

Page 33 <sup>€</sup><sub>0123</sub>

Code	Fr/Ch	Length	Sales Unit
8884751622E	16	PEG adaptor – ENFit™	20
8884752042E	20	PEG adaptor – ENFit™	20

**REPLACEMENT LONG TERM ACCESS – SKIN LEVEL DEVICES**

**Kangaroo™ Skin Level Balloon Gastrostomy Kit – Silicone, Sterile**

Page 34 <sup>€</sup><sub>0123</sub>

Code	Fr/Ch	Length	Sales Unit
712080	12	0.8 cm	1
712100	12	1.0 cm	1
712120	12	1.2 cm	1
712150	12	1.5 cm	1
712170	12	1.7 cm	1
712200	12	2.0 cm	1

712230	12	2.3 cm	1
712250	12	2.5 cm	1
712270	12	2.7 cm	1
712300	12	3.0 cm	1
712350	12	3.5 cm	1
712400	12	4.0 cm	1
712450	12	4.5 cm	1
712500	12	5.0 cm	1
714080	14	0.8 cm	1
714100	14	1.0 cm	1
714120	14	1.2 cm	1
714150	14	1.5 cm	1
714170	14	1.7 cm	1
714200	14	2.0 cm	1
714230	14	2.3 cm	1
714250	14	2.5 cm	1
714270	14	2.7 cm	1
714300	14	3.0 cm	1
714350	14	3.5 cm	1
714400	14	4.0 cm	1
714450	14	4.5 cm	1
714500	14	5.0 cm	1
716080	16	0.8 cm	1
716100	16	1.0 cm	1
716120	16	1.2 cm	1
716150	16	1.5 cm	1
716170	16	1.7 cm	1
716200	16	2.0 cm	1
716230	16	2.3 cm	1
716250	16	2.5 cm	1
716270	16	2.7 cm	1
716300	16	3.0 cm	1
716350	16	3.5 cm	1
716400	16	4.0 cm	1
716450	16	4.5 cm	1
716500	16	5.0 cm	1
718080	18	0.8 cm	1
718100	18	1.0 cm	1
718120	18	1.2 cm	1
718150	18	1.5 cm	1

**REPLACEMENT LONG TERM ACCESS – SKIN LEVEL DEVICES**

**Kangaroo™ Skin Level Balloon Gastrostomy Kit – Silicone, Sterile**

Page 35 <sup>€</sup><sub>0123</sub>

Code	Fr/Ch	Length	Sales Unit
718170	18	1.7 cm	1
718200	18	2.0 cm	1
718230	18	2.3 cm	1
718250	18	2.5 cm	1
718270	18	2.7 cm	1
718300	18	3.0 cm	1
718350	18	3.5 cm	1
718400	18	4.0 cm	1
718450	18	4.5 cm	1
718500	18	5.0 cm	1
720080	20	0.8 cm	1
720100	20	1.0 cm	1
720120	20	1.2 cm	1
720150	20	1.5 cm	1
720170	20	1.7 cm	1
720200	20	2.0 cm	1
720230	20	2.3 cm	1
720250	20	2.5 cm	1
720270	20	2.7 cm	1
720300	20	3.0 cm	1
720350	20	3.5 cm	1
720400	20	4.0 cm	1
720450	20	4.5 cm	1
720500	20	5.0 cm	1
724080	24	0.8 cm	1
724100	24	1.0 cm	1
724120	24	1.2 cm	1
724150	24	1.5 cm	1
724170	24	1.7 cm	1
724200	24	2.0 cm	1
724230	24	2.3 cm	1
724250	24	2.5 cm	1
724270	24	2.7 cm	1

724300	24	3.0 cm	1
724350	24	3.5 cm	1
724400	24	4.0 cm	1
724450	24	4.5 cm	1
724500	24	5.0 cm	1

### Skin Level Gastrostomy Accessory Products Page 36 <sup>€</sup><sub>0123</sub>

Code	Fr/Ch	Length	Sales Unit
8884741813		Bolus feeding set, non-sterile	10
8884761126		Bolus feeding set, sterile	5
8884741821		Continuous feeding set, non-sterile	10
8884761100		Continuous feeding set, sterile	5

## KANGAROO™ GASTROSTOMY FEEDING TUBE WITH ENFIT™ Y-PORT

### KANGAROO™ gastrostomy feeding tube Page 36 <sup>€</sup><sub>0123</sub>

Code	Description	Sales Unit
8884715122E	12 Fr/Ch diameter, 5 ml balloon	5
8884715148E	14 Fr/Ch diameter, 5 ml balloon	5
8884720163E	16 Fr/Ch diameter, 20 ml balloon	5
8884720189E	18 Fr/Ch diameter, 20 ml balloon	5
8884720205E	20 Fr/Ch diameter, 20 ml balloon	5
8884720221E	22 Fr/Ch diameter, 20 ml balloon	5
8884720247E	24 Fr/Ch diameter, 20 ml balloon	5
8884720265E	26 Fr/Ch diameter, 20 ml balloon	5
8884720285E	28 Fr/Ch diameter, 20 ml balloon	5

### Accessory Page 36 <sup>€</sup><sub>0123</sub>

Code	Description	Sales Unit
8884752042E	Replacement ENFit™ Y-port	20

## OBSTETRICS AND PERINATAL

### FETAL SPIRAL ELECTRODES

#### Kendall™ Fetal Spiral Electrode System Page 38

Code	Description	Sales Unit
<b>KENDALL™ FETAL SPIRAL ELECTRODES</b>		
31479549	Single helix with safety connector <sup>€</sup>	50
<b>KENDALL™ FETAL SPIRAL ELECTRODE ATTACHMENT PADS</b>		
31481396	Small attachment pad 4.8 cm <sup>€</sup>	50
50000095	Large attachment pad 9.5 cm x 4.8 cm <sup>€</sup>	50
<b>KENDALL™ FETAL SPIRAL ELECTRODE CABLES <sup>€</sup></b>		
The system has cables for all major fetal monitors		

### BELTS

#### Kendall™ Abdominal Belts Page 38 <sup>€</sup>

Code	Size	frequency	Hole ends	Finished Colour	Sales Unit
<b>KENDALL™ BUTTONHOLE BELTS</b>					
31410270	6 cm x 1.2 m	3.2 cm	Yes	Pink and blue	50 pairs
40000007	6 cm x 1.2 m	3.2 cm	Yes	Pink and blue	100 pairs

## ACCESSORIES FOR NEONATOLOGY

#### Argyle™ IV Support Boards Page 38 <sup>€</sup>

Code	Description	Dimension	Sales Unit
MI00880	Argyle™ I.V. Support Board, Preemie	2.5 cm x 7.6 cm	50
MI00775	Argyle™ I.V. Support Board, Small	4.1 cm x 9.8 cm	50
MI00850	Argyle™ I.V. Support Board, Large	4.8 cm x 12.4 cm	50

#### Argyle™ Replogle Suction Catheters Page 39 <sup>€</sup><sub>0123</sub>

Code	Fr size	Length	Sales Unit
8888256545	6	61 cm	10/case
8888256529	8	61 cm	10/case
8888256503	10	61 cm	10/case

## OPERATING ROOM

### DEVON™ INSTRUMENT HANDLING PRODUCTS

#### Devon™ Sponge Counting System Page 41

Code	Description	Sales Unit
23250-50C	Sponge Counter Bag 5-10 count Devon™ Clear Non-Sterile	1
23250-50W	Sponge Counter Bag 5-10 count Devon™ White Non-Sterile	1
31144119	Dispenser box, clear back*	50 box/250 case
31144127	Dispenser box, opaque back*	50 box/250 case
31144168	Dispenser rack for I.V. pole*	1
31144176	Dispenser rack for wall*	1
31144184	Dispenser rack for I.V. pole, triple*	1

These codes are not medical devices and are not CE marked.

### DEVON™ HANDS-FREE-TRANSFER PRODUCTS

#### Devon™ Magnetic Drapes Page 41 <sup>€</sup><sub>0123</sub>

Code	Description	Size	Sales Unit
31140570	Magnetic drape	25.4 cm x 40.6 cm	30
31140562	Magnetic drape	40.6 cm x 50.8 cm	20
31140596	Magnetic drape with stays	25.4 cm x 40.6 cm	30
31140547	Magnetic drape with stays	40.6 cm x 50.8 cm	20

#### Devon™ Magnetic Drapes with Neutral Zone Page 42 <sup>€</sup><sub>0123</sub>

Code	Description	Neutral Zone	Size	Sales Unit
31140612	Magnetic drape with stay	1	25.4 cm x 40.6 cm	30
31140539	Magnetic drape with stays	1	40.6 cm x 50.8 cm	20
31140505	Magnetic drape	2	25.4 cm x 40.6 cm	26
31140430	Magnetic drape	2	40.6 cm x 50.8 cm	16

## DEVON™ SURGICAL MARKERS

#### Devon™ Surgical Markers with Ruler Cap Page 42 <sup>€</sup><sub>0123</sub>

Code	Description	Sales Unit
31145926	Skin marker ruler cap-reg	25 box/100 case
31145918	Skin marker ruler cap-reg-ruler/label	25 box/100 case
31145942	Skin marker ruler cap-fine	25 box/100 case
31145892	Skin marker ruler cap-fine-ruler	25 box/100 case
31145843	Skin marker ruler cap-fine-ruler/label	25 box/100 case

#### Devon™ Dual Tip Surgical Markers with Ruler Cap Page 43 <sup>€</sup><sub>0123</sub>

Code	Description	Sales Unit
31145900	Dual Tip	25 box/100 case
31145868	Dual Tip with 9 labels	25 box/100 case
31160560	Dual tip with Flexi-ruler	25 box/100 case
31145934	Dual tip with 9 labels and flexi ruler	25 box/100 case

#### Devon™ Surgical Markers Page 43 <sup>€</sup><sub>0123</sub>

Code	Description	Sales Unit
19005-01CE	Surgical marking pen traditional regular tip	25 box/100 case
19005-02CE	Surgical marking pen traditional regular tip with ruler	25 box/100 case
19005-03CE	Surgical marking pen traditional regular tip with ruler and labels	25 box/100 case
19004-01CE	Surgical marking pen traditional fine tip	25 box/100 case
19004-02CE	Surgical marking pen traditional fine tip with ruler	25 box/100 case
19004-03CE	Surgical marking pen traditional fine tip with ruler and labels	25 box/100 case

**Devon™ Skin Marker – Non-Sterile** Page 43 CE 0123

Code	Description	Sales Unit
31146020	Skin marker, non-sterile	200/case

**DEVON™ UTILITY MARKERS** Page 43 CE 0123

Code	Description	Sales Unit
31145884	Black ink, labels, sterile	25 box/100 case

**DEVON™ LIGHT HANDLE SYSTEMS**

**Devon™ Light Glove** Page 44 CE 2787

Code	Description	Sales Unit
LT-F01B	Light glove	1 pack/25 box/100 case
LT-F02	Light glove	2 pack/40 box/160 case
LT-F03	Light glove	3 pack/48 box/192 case

**Devon™ Light Handle** Page 44 CE 0123

Code	Description	Sales Unit
31140133	Light handle	1 pack/16 box/128 case
31140125	Light handle	2 pack/28 box/112 case
31140117	Light handle	3 pack/21 box/84 case

**Devon™ Operating Room, Camera Light Cover** Page 45 CE 0123

Code	Description	Sales Unit
15010-1CE	OR Camera Light Cover	10 (10 covers per box)

**Devon™ Light Sleeve** Page 45 CE 0123

Code	Description	Sales Unit
31140273	Light sleeve	1 pack/70 case
31140240	Light sleeve	2 pack/120 case

**DEVON™ NEEDLE COUNTERS**

**Devon™ Foam Block Needle Counters with Boxlocks** Page 46

Code	Description	Count Capacity	Sales Unit
31142352	Blade removal, double black magnet*	100/100	8 box/64 case
31142360	Blade removal, adhesive*	20/40	12 box/96 case
31142386	Blade removal, single black magnet*	20/40	12 box/96 case
31142394	Blade removal, single black magnet*	10/20	12 box/96 case
31142469	Double foam*	80/80	8 box/64 case
31142402	Adhesive*	40/70	8 box/64 case
31142311	Single black magnet*	40/70	8 box/64 case
31142501	Adhesive*	20/40	12 box/96 case
31142493	Single black magnet*	20/40	12 box/96 case

\*This range of products are not medical devices and they are not CE marked.

**DEVON™ NEEDLE COUNTERS**

**Devon™ Double Blade Removers** Page 46

Code	Description	Count Capacity	Sales Unit
31181459	Single foam block with magnet*	30/30	12 box/96 case
31181467	Double magnet*	25/25	12 box/96 case
31313425	Single foam strips with single magnet*	25/25	12 box/96 case

\*This range of products are not medical devices and they are not CE marked.

**Devon™ Magnetic Needle Counters** Page 47

Code	Description	Count Capacity	Sales Unit
31142238	Boxlocks, blade removal, double white magnet*	30/30	12 box/96 case
31142113	Boxlocks, blade removal, double black magnet*	30/30	12 box/96 case
31142154	Boxlocks, double black magnet*	60/60	8 box/64 case
31142147	Boxlocks, double black magnet*	30/30	12 box/96 case
31142196	Double black magnet*	20/20	21 box/168 case
31142139	Clear top, boxlocks, single black magnet*	30/30	8 box/64 case
31142410	Clear top, single black magnet*	10/10	21 box/168 case

\*This range of products are not medical devices and they are not CE marked.

**DEVON™ NEEDLE COUNTERS**

**Devon™ Foam Strip Needle Counters with Boxlocks** Page 47

Code	Description	Count Capacity	Sales Unit
31142451	Blade removal, double foam*	30/30	12 box/96 case
31142295	Blade removal, black magnet*	15/30	12 box/96 case
31142303	Double foam*	30/30	12 box/96 case
31142188	Black magnet*	30/60	8 box/64 case
31142212	Clear top*	15/15	12 box/96 case

\*This range of products are not medical devices and they are not CE marked.

**Devon™ Needle Counters – Non-Sterile** Page 47

Code	Description	Count Capacity	Sales Unit
31155370	Needle counter, plastic box, single magnet, clear top*	10	200
31155362	Needle counter, plastic box, single foam strips, clear top*	15	200
31155487	Needle counter, plastic box, double magnet*	20	200
31155263	BladeGuard, plastic box, double white magnet*	30	200
31155552	Pop'n Count plastic box, foam block and magnet, medium*	10	200
31155503	Pop'n Count plastic box, foam block and adhesive, medium*	10	200
31155545	Pop'n Count BladeGuard plastic box, foam block and magnet, medium*	10	200
31155651	Pop'n Count BladeGuard plastic box, foam box and magnet, extra large*	20	200
31155107	Pop'n Count BladeGuard plastic box, foam block and double magnet, extra large*	80	100

\*This range of products are not medical devices and they are not CE marked.

**DEVON™ OR ESSENTIALS**

**Devon™ Cautery Tipolisher** Page 48 CE 0123

Code	Description	Sales Unit
300-50CE	Cleaner cautery tip sterile	25 box/100 case

**DEVON™ POSITIONING PRODUCTS**

**Head Positioning** Page 48 CE

Code	Description	Sales Unit
FP-HDCCE	Child head positioner	1 pack/24 case
FP-HDCRCE	Adult head positioner	1 pack/24 case
FP-HDSLCE	Slotted adult head positioner	1 pack/24 case
FP-HEAD9CE	Head Donut 9in head positioner	1 pack/36 case
FP-HEAD7CE	Head Donut 7in Head Positioner	1 pack/72 case
FP-HDMCE	Foam Multi Ring Head Donut	1 pack/36 case
FP-HDFRCE	Frame Head Positioner	1 pack/28 case
FP-HDSFCE	Soft touch head positioner	1 pack/28 case

Code	Description	Size	Sales Unit
FP-UNICE	Convuluted ulnar nerve protector	41 x 15 x 5 cm	2 pack/72 case
FP-LARMCCE	Foam Laminectomy Arm Cradle	61 x 13 x 8 cm	2 pack/24 case
FP-ARTKNCE	Arthroscopic knee holder	-	1 pack/10 case
FP-ARTLECE	Well leg holder	-	1 pack/18 case
FP-HEELICE	Convuluted foot & heel protector	-	2 pack/72 case
FP-ORTB2CE	Convuluted OR table pad	183 x 51 x 5 cm	1 pack/12 case
FP-CHESTCE	Chest roll	43 x 13 cm	2 pack/12 case
FP-OVERICE	Utility pad	51 x 30 x 30 cm	1 pack/12 case
FP-ARMBICE	Convuluted arm board pad	51 x 20 x 5 cm	12 pair
FP-ABDSCE	Small abduction pillow	46 x 30 x 15 cm	6 case
FP-ABDMCE	Medium abduction pillow	56 x 38 x 15 cm	6 case
FP-ABDLCE	Large abduction pillow	63 x 46 x 15 cm	6 case

Body Positioners Page 49 <sup>CE</sup>

Code	Description	Size	Sales Unit
FP-BAW2CE	Foam Body Alignment Wedge 18 inch	46 x 20 x 20 cm	8
FP-BAW1CE	Foam Body Alignment Wedge 22 inch	55 x 30 x 22 cm	8

DEVON™ VESSEL LOOPS Page 49 <sup>CE</sup>

Made of Silicone. Radiopaque bands that find practical application in situation where arteries veins, nerves or tendons must be looped for identification, retraction or occlusion. Available in four different colors and 2 sizes.

Code	Description	Sales Unit
2400-1RCE	Silicone Vessel Loop, Red, Mini, Radiopaque, 2/PK, Sterile	20 pack/200 per case
2400-1BCE	Silicone Vessel Loop, Blue, Mini, Radiopaque, 2/PK, Sterile	20 pack/200 per case
2400-1YCE	Silicone Vessel Loop, Yellow, Mini, Radiopaque, 2/PK, Sterile	20 pack/200 per case
2400-1WCE	Silicone Vessel Loop, White, Mini, Radiopaque, 2/PK, Sterile	20 pack/200 per case
2400-2RCE	Silicone Vessel Loop, Red, Maxi, Radiopaque, 2/PK, Sterile	20 pack/200 per case
2400-2BCE	Silicone Vessel Loop, Blue, Maxi, Radiopaque, 2/PK, Sterile	20 pack/200 per case
2400-2YCE	Silicone Vessel Loop, Yellow, Maxi, Radiopaque, 2/PK, Sterile	20 pack/200 per case
2400-2WCE	Silicone Vessel Loop, White, Maxi, Radiopaque, 2/PK, Sterile	20 pack/200 per case

# SAFETY NEEDLES & SYRINGES AND BLOOD COLLECTION

## SAFETY NEEDLES

Magellan™ Safety Hypodermic Needles Page 51 <sup>CE</sup>

Code	Needle Gauge and Length	Sales Unit
8881850558	25G x 5/8" (16 mm)	10 x 50 pcs
8881850510	25G x 1" (25 mm)	10 x 50 pcs
8881850358	23G x 5/8" (16 mm)	10 x 50 pcs
8881850310	23G x 1" (25 mm)	10 x 50 pcs
8881850210	22G x 1" (25 mm)	10 x 50 pcs
8881850215	22G x 1-1/2" (38 mm)	10 x 50 pcs
8881850158	21G x 5/8" (16 mm)	10 x 50 pcs
8881850110	21G x 1" (25 mm)	10 x 50 pcs
8881850115	21G x 1-1/2" (38 mm)	10 x 50 pcs
8881850010	20G x 1" (25 mm)	10 x 50 pcs
8881850015	20G x 1-1/2" (38 mm)	10 x 50 pcs
8881850910	19G x 1" (25 mm)	10 x 50 pcs
8881850915	19G x 1-1/2" (38 mm)	10 x 50 pcs
8881850810	18G x 1" (25 mm)	10 x 50 pcs
8881850815	18G x 1-1/2" (38 mm)	10 x 50 pcs

Magellan™ Insulin and Tuberculin Safety Syringe Page 51 <sup>CE</sup>

Code	Size	Needle Gauge and Length	Colour	Cannula	Sales Unit
<b>INSULIN</b>					
8881892910	1 ml	29G x 1/2" (13 mm)	Orange	Fixed	50 box/500 case
8881892950	0.5 ml	29G x 1/2" (13 mm)	Orange	Fixed	50 box/500 case
8881892930	0.3 ml	29G x 1/2" (13 mm)	Orange	Fixed	50 box/500 case
8881893010	1 ml	30G x 5/16" (8 mm)	Orange	Fixed	50 box/500 case
8881893050	0.5 ml	30G x 5/16" (8 mm)	Orange	Fixed	50 box/500 case
8881893030	0.3 ml	30G x 5/16" (8 mm)	Orange	Fixed	50 box/500 case
<b>TUBERCULIN</b>					
8881882812	1 ml	28G x 1/2" (13 mm)	Blue	Fixed	50 box/500 case

## SAFETY NEEDLES

Monoject™ Safety Insulin Syringes – Sterile Page 52 <sup>CE</sup>

Code	Size	Graduation	Needle Gauge and Length	Cannula	Sales Unit
8881511344	0.3 ml	U-100	30G x 11/16" (8 mm)	Fixed	5 x 100 pcs
8881511336	0.5 ml	U-100	30G x 11/16" (8 mm)	Fixed	5 x 100 pcs

Monoject™ Hypodermic Safety Needle Page 52 <sup>CE</sup>

Code	Pack Colour	Needle Gauge and Length	Box/Case Qty
1181810		18G x 1" (1.27 mm x 2.5 cm)	100/800
1181815		18G x 1-1/2" (1.27 mm x 3.8 cm)	100/800
1181910		19G x 1" (1.067 mm x 2.5 cm)	100/800
1181915		19G x 1-1/2" (1.067 mm x 3.8 cm)	100/800
1182010		20G x 1" (0.902 mm x 2.5 cm)	100/800
1182015		20G x 1-1/2" (0.902 mm x 3.8 cm)	100/800
1182110		21G x 1" (0.813 mm x 2.5 cm)	100/800
1182115		21G x 1-1/2" (0.813 mm x 3.8 cm)	100/800
1182210		22G x 1" (0.711 mm x 2.5 cm)	100/800
1182215		22G x 1-1/2" (0.711 mm x 3.8 cm)	100/800
1182310		23G x 1" (0.635 mm x 2.5 cm)	100/800
1182558		25G x 5/8" (0.508 mm x 1.6 cm)	100/800
1182510		25G x 1" (0.508 mm x 2.5 cm)	100/800
1182705		27G x 1/2" (0.406 mm x 1.3 cm)	100/800
1182734		27G x 3/4" (0.406 mm x 1.9 cm)	100/800
1183005		30G x 1/2" (0.305 mm x 1.3 cm)	100/800

Monoject™ Bluntfill Needle Page 52 <sup>CE</sup>

Code	Needle Gauge and Length	Box/Case Qty
11811022F	Bluntfill with filter 18G x 1-1/2" (1.27 mm x 3.8 cm)	100/1000

## SHARPS DISPOSAL




### PHLEBOTOMY AND MULTI-PURPOSE CONTAINERS

SharpSafety Phlebotomy Sharps Containers Page 54

Code	Description	Litres Capacity	Dimensions HxDxW	Brackets & Holders	Colour	Sales Unit
89005A	SharpSafety phlebotomy sharps container adapts to present blood drawing trays and techniques*	1	15.9 x 10.8 x 11.4 cm	1		100
8906	Small Monoject™ sharps container, with clear top for easier assessment of fill level*	1.5	13.3 x 20 x 13.3 cm	-		40
15225A	SharpSafety phlebotomy sharps container with rotor opening*	2	18.4 x 12 x 16.5 cm	2		60

\*This range of products are not medical devices and they are not CE marked.






## MONOJECT™ MULTI-PURPOSE SHARPS CONTAINERS Page 54

Code	Description	Litres Capacity	Dimensions HxDxW	Brackets & Holders	Colour	Sales Unit
1100676376 1100676517	Monoject™ sharps container with chimney drop top lid. Clear top on translucent container for easier assessment of fill level. Can be decontaminated in an autoclave*	4	13.9 x 17.3 x 28.6 cm	11		40
1100676384 1100676533	Monoject™ sharps container with chimney drop top lid. Clear top on translucent container for easier assessment of fill level. Can be decontaminated in an autoclave*	8	23.7 x 17.3 x 28.6 cm	11		20
1100676392 1100676558	Monoject™ sharps container with chimney drop top lid. Clear top on translucent container for easier assessment of fill level. Can be decontaminated in an autoclave*	11	40.8 x 17.3 x 28.6 cm	11		10

\* This range of products are not medical devices and they are not CE marked.

## MULTI-PURPOSE CONTAINERS




### SharpSafety Multi-Purpose Sharps Containers Page 55

Code	Description	Litres Capacity	Dimensions HxDxW	Brackets & Holders	Colour	Sales Unit
8970 8970Y	Multi-purpose container for all sharps waste with rotor opening*	7.4	25.4 x 18.4 x 26.6 cm	2 7 8		20
8990SA	Multi-purpose container with clear hinged lid*	7.4	25.4 x 18.4 x 26.6 cm	2 7 8		20
8990Y	Multi-purpose container with special gasket hinged lid*	7.4	25.4 x 18.4 x 26.6 cm	2 7 8		20
8967Y	Multi-purpose container with horizontal drop opening lid*	7.4	32.4 x 18.4 x 26.6 cm	2 7 8		20
8527R	Multi-purpose container with rotor opening and hinged lid to accommodate a variety of sharps sizes and provide temporary and permanent closure lid*	11	35 x 15.2 x 35 cm	4 9		10

\* This range of products are not medical devices and they are not CE marked.



## IN-ROOM SYSTEM

### SharpSafety Counterbalanced Lid Containers Page 55

Code	Description	Litres Capacity	Dimensions HxDxW	Brackets & Holders	Colour	Sales Unit
8506Y 8507SA	Sharps container with counterbalanced lid and counter-balanced door*	4.6	31.7 x 14 x 27.3 cm	3 6 15		20
8534SA	Sharps container with counterbalanced lid and counter-balanced door*	7.4	30 x 16.5 x 35 cm	9 16		10
8536SA 8537SA 8537Y	Sharps container with counterbalanced lid and counter-balanced door*	11	42 x 15.2 x 35 cm	4 9 16		10

\* This range of products are not medical devices and they are not CE marked.



## SharpSafety Always-Open Lid Containers Page 56

Code	Description	Litres Capacity	Dimensions HxDxW	Brackets & Holders	Colour	Sales Unit
85321R	Sharps container with always-open lid provides convenient, quick and easy point-of-use sharps disposal*	7.4	29.2 x 15.2 x 35 cm	9 14		10
85221R	Sharps container with always-open lid provides convenient, quick and easy point-of-use sharps disposal*	11	41.3 x 15.2 x 35 cm	4 9 15		10

\* This range of products are not medical devices and they are not CE marked.

## LARGE VOLUME CONTAINERS




### SharpSafety Large Volume Containers Page 56

Code	Description	Litres Capacity	Dimensions HxDxW	Brackets & Holders	Colour	Sales Unit
8980S, 8935, 8935Y, 8938	Large volume sharps container with sliding lid*	30 45 67	45 x 28 x 39.4 cm 47.6 x 32.4 x 46.3 cm 66 x 32.4 x 46.3 cm	18, 19 17, 18, 22 17, 18, 22		10 10 5
8980, 8980Y 8991	Large volume sharps container hinged lid and separate rotor opening for smaller sharps waste*	30 45 67	45 x 28 x 39.4 cm 47.6 x 32.4 x 46.3 cm 66 x 32.4 x 46.3 cm	18, 20 18, 21 18, 21		10 10 5

\* This range of products are not medical devices and they are not CE marked.

## CHEMOTHERAPY CONTAINERS


### SharpSafety Large Volume Chemotherapy Containers Page 57 <sup>CE 0123</sup>

Code	Description	Litres Capacity	Dimensions HxDxW	Brackets & Holders	Colour	Sales Unit
8982	Chemotherapy sharps container with hinged lid	7.4	25.4 x 18.4 x 26.7 cm	2, 7, 8		20
8985 8931 8989	Chemotherapy sharps container with hinged lid	30 45 67	45 x 28 x 39.4 cm 47.6 x 32.4 x 46.3 cm 66 x 32.4 x 46.3 cm	18, 20 18, 21 18, 21		10 10 5
8934 8939	Chemotherapy sharps container with sliding lid	45 65	47.6 x 32.4 x 46.3 cm 66 x 32.4 x 46.3 cm	17, 18, 22 17, 18, 22		10 5

\* This range of products are not medical devices and they are not CE marked.


## PGII LARGE VOLUME SHARPS CONTAINERS

### SharpSafety Large Volume PGII Containers Page 58

Code	Description	Litres Capacity	Dimensions HxDxW	Brackets & Holders	Colour	Sales Unit
8998PG2	PGII container large volume, with hinged and sealing gasket lid	30 45 67	45 x 28 x 39.4 cm 47.6 x 32.4 x 46.3 cm 66 x 32.4 x 46.3 cm	18, 20 18, 21 18, 21		10 10 5

\* This range of products are not medical devices and they are not CE marked.

### SharpSafety Chemotherapy PGII Containers Page 58

Code	Description	Litres Capacity	Dimensions HxDxW	Brackets & Holders	Colour	Sales Unit
8985PG2 8989PG2	Chemotherapy sharps container with hinged lid	30 45 67	45 x 28 x 39.4 cm 47.6 x 32.4 x 46.3 cm 66 x 32.4 x 46.3 cm	18, 19 18, 21 18, 21		10 10 5

\* This range of products are not medical devices and they are not CE marked.



## BRACKETS AND HOLDERS

### Brackets and Holders Page 58

Code	Description	Dimensions HxDxW	Sales Unit
8974	Holder for 7.4 litre phlebotomy, multi-purpose and ChemoSafety containers. For containers code 1525SA, 8996, 8970, 8970Y, 89671, 8967Y, 8982, 8990SA, 8990Y	8.2 x 19 x 28 cm	5
8518X	Locking brackets for 2 & 4.6 litre In-Room™ containers. For containers code 85121, 850201, 850301, 851201, 85120Y	21 x 16.5 x 28.6 cm	5
8975	Wall/cart bracket for 7.4 litre multi-purpose and ChemoSafety containers. For containers code 8970, 8970Y, 89671, 8967Y, 8990Y, 8996, 8982	19 x 16.5 x 25.4 cm	5
8912SA	Wall Bracket for 1 litre phlebotomy container. For containers code 8900SA, 8906	–	5

\* This range of products are not medical devices and they are not CE marked.

## WALL ENCLOSURES

### Wall Enclosures Page 59

Code	Description	Dimensions HxDxW	Sales Unit
8556H	Wall enclosure for 2 & 4.6 litre sharps container with glove dispenser and incl. a 4.6 litre SharpSafety™ horizontal lid container	47 x 14 x 33 cm	1
8550LG	Glove dispenser for 85161H and 85301H wall enclosure	17.1 x 11.4 x 30.5 cm	10
8550B	Separate glove dispenser with attached wall mounting bracket	18.4 x 11.4 x 31.1 cm	5

\* This range of products are not medical devices and they are not CE marked.

## FOOT PEDAL CARTS AND FLOOR BRACKETS

### Foot Pedal Carts and Floor Brackets Page 59

Code	Description	Dimensions HxDxW	Sales Unit
8938FP	Sharps cart foot pedal cart for 45 & 67 litre sliding lid large volume and ChemoSafety containers (container not included)	78.7 x 57.8 x 38.7 cm	1
8992H	Floor bracket with wheels for 30, 45 & 67 litre large volume and ChemoSafety containers (container not included)	55.9 x 32.4 x 57.1 cm	1
8981FP	Sharps cart foot pedal cart for 30 litre sliding lid large volume and ChemoSafety containers (container not included)	55.9 x 54.6 x 40.6 cm	1
8980FP	Sharps cart foot pedal cart for 30 litre hinged lid large volume and ChemoSafety containers (container not included)	55.9 x 54.6 x 40.6 cm	1
8991FP	Sharps cart foot pedal cart for 45 & 67 litre hinged lid large volume and ChemoSafety containers (container not included)	78.7 x 57.8 x 38.7 cm	1

\* This range of products are not medical devices and they are not CE marked.

## SURGICAL GLOVES

### PROTEXIS™ POLYISOPRENE, POWDER-FREE

#### Protexis™ PI Surgical Gloves Page 61 <sup>CE</sup><sub>2797</sub>

Code	Size	Material	Color	Category	Fingertip thickness	Coatings	Sales Unit
2D72PT55X	5½	Synthetic polyisoprene (PI)	Cream	Standard	0.23 mm		50 pairs/box
2D72PT60X	6						4 boxes/case
2D72PT65X	6½						
2D72PT70X	7						
2D72PT75X	7½						
2D72PT80X	8						
2D72PT85X	8½						
2D72PT90X	9						

#### Protexis™ PI Micro Surgical Gloves Page 61 <sup>CE</sup><sub>2797</sub>

Code	Size	Material	Color	Category	Fingertip thickness	Coatings	Sales Unit
2D73PM55	5½	Synthetic polyisoprene (PI)	Cream	Thin	0.20 mm		50 pairs/box
2D73PM60	6						4 boxes/case
2D73PM65	6½						
2D73PM70	7						
2D73PM75	7½						
2D73PM80	8						
2D73PM85	8½						
2D73PM90	9						

#### Protexis™ PI with Neu-Thera™ Surgical Gloves Page 61 <sup>CE</sup><sub>2797</sub>

Code	Size	Material	Color	Category	Fingertip thickness	Coatings	Sales Unit
2D73TE55	5½	Synthetic polyisoprene (PI)	Cream	Standard	0.23 mm	Neu-Thera™ Emollient Coating	50 pairs/box
2D73TE60	6						4 boxes/case
2D73TE65	6½						
2D73TE70	7						
2D73TE75	7½						
2D73TE80	8						
2D73TE85	8½						
2D73TE90	9						

#### Protexis™ PI Blue with Neu-Thera™ Surgical Gloves Page 61 <sup>CE</sup><sub>2797</sub>

Code	Size	Material	Color	Category	Fingertip thickness	Coatings	Sales Unit
2D73EB55	5½	Synthetic polyisoprene (PI)	Royal blue	Underglove	0.20 mm	Neu-Thera™ Emollient Coating	50 pairs/box
2D73EB60	6						4 boxes/case
2D73EB65	6½						
2D73EB70	7						
2D73EB75	7½						
2D73EB80	8						
2D73EB85	8½						
2D73EB90	9						

#### Protexis™ PI Ortho Surgical Gloves Page 62 <sup>CE</sup><sub>2797</sub>

Code	Size	Material	Color	Category	Fingertip thickness	Coatings	Sales Unit
2D73ET60	6	Synthetic polyisoprene (PI)	Cream	Thick	0.30 mm		40 pairs/box
2D73ET65	6½						4 boxes/case
2D73ET70	7						
2D73ET75	7½						
2D73ET80	8						
2D73ET85	8½						
2D73ET90	9						

**Protexis™ PI Orthopaedic Surgical Gloves** Page 62 CE 2797

Code	Size	Material	Color	Category	Fingertip thickness	Coatings	Sales Unit
2D73HT60	6	Synthetic polyisoprene (PI)	Brown	Thick	0.34 mm	Hydrogel	40 pairs/box
2D73HT65	6½						4 boxes/case
2D73HT70	7						
2D73HT75	7½						
2D73HT80	8						
2D73HT85	8½						
2D73HT90	9						

**Protexis™ PI Textured Surgical Gloves** Page 62 CE 2797

Code	Size	Material	Color	Category	Fingertip thickness	Coatings	Sales Unit
2D72TG55	5½	Synthetic polyisoprene (PI)	Cream	Medium	0.27 mm	Hydrogel	50 pairs/box
2D72TG60	6						4 boxes/case
2D72TG65	6½						
2D72TG70	7						
2D72TG75	7½						
2D72TG80	8						
2D72TG85	8½						
2D72TG90	9						

**PROTEXIS™ NEOPRENE, POWDER-FREE**

**Protexis™ Neoprene Surgical Gloves** Page 62 CE 2797

Code	Size	Material	Color	Category	Fingertip thickness	Coatings	Sales Unit
2D73DP55	5½	Synthetic neoprene	Light brown	Standard	0.17 mm	Nitrile	50 pairs/box
2D73DP60	6						4 boxes/case
2D73DP65	6½						
2D73DP70	7						
2D73DP75	7½						
2D73DP80	8						
2D73DP85	8½						
2D73DP90	9						

**Protexis™ Neoprene Essential Surgical Gloves** Page 62 CE 2797

Code	Size	Material	Color	Category	Fingertip thickness	Coatings	Sales Unit
2D73DS55	5½	Synthetic neoprene	Light brown	Standard	0.17 mm	Nitrile	50 pairs/box
2D73DS60	6						4 boxes/case
2D73DS65	6½						
2D73DS70	7						
2D73DS75	7½						
2D73DS80	8						
2D73DS85	8½						
2D73DS90	9						

**PROTEXIS™ LATEX, POWDER-FREE**

**Protexis™ Latex Surgical Gloves** Page 63 CE 2797

Code	Size	Material	Color	Category	Fingertip thickness	Coatings	Sales Unit
2D72NS55X	5½	Latex	Light brown	Standard	0.25 mm	Nitrile	50 pairs/box
2D72NS60X	6						4 boxes/case
2D72NS65X	6½						
2D72NS70X	7						
2D72NS75X	7½						
2D72NS80X	8						
2D72NS85X	8½						
2D72NS90X	9						

**Protexis™ Latex Essential Surgical Gloves** Page 63 CE 2797

Code	Size	Material	Color	Category	Fingertip thickness	Coatings	Sales Unit
2D72LE55	5½	Latex	Cream	Standard	0.19 mm	Nitrile	50 pairs/box
2D72LE60	6						4 boxes/case
2D72LE65	6½						
2D72LE70	7						
2D72LE75	7½						
2D72LE80	8						
2D72LE85	8½						
2D72LE90	9						

**Protexis™ Latex Micro Surgical Gloves** Page 63 CE 2797

Code	Size	Material	Color	Category	Fingertip thickness	Coatings	Sales Unit
2D72NT55X	5½	Latex	Light brown	Thin	0.17 mm	Nitrile	50 pairs/box
2D72NT60X	6						4 boxes/case
2D72NT65X	6½						
2D72NT70X	7						
2D72NT75X	7½						
2D72NT80X	8						
2D72NT85X	8½						
2D72NT90X	9						

**Protexis™ Latex with Neu-Thera™ Surgical Gloves** Page 63 CE 2797

Code	Size	Material	Color	Category	Fingertip thickness	Coatings	Sales Unit
2D73TP55	5½	Latex	Light brown	Standard	0.25 mm	Nitrile and Neu-Thera™ Emollient Coating	50 pairs/box
2D73TP60	6						4 boxes/case
2D73TP65	6½						
2D73TP70	7						
2D73TP75	7½						
2D73TP80	8						
2D73TP85	8½						
2D73EB90	9						

**Protexis™ Latex Blue with Neu-Thera™ Surgical Gloves** Page 64 CE 2797

Code	Size	Material	Color	Category	Fingertip thickness	Coatings	Sales Unit
2D72LU55	5½	Latex	Blue	Underglove	0.19 mm	Nitrile and Neu-Thera™ Emollient Coating	50 pairs/box
2D72LU60	6						4 boxes/case
2D72LU65	6½						
2D72LU70	7						
2D72LU75	7½						
2D72LU80	8						
2D72LU85	8½						
2D72LU90	9						

**Protexis™ Latex Ortho Surgical Gloves** Page 64 CE 2797

Code	Size	Material	Color	Category	Fingertip thickness	Coatings	Sales Unit
2D72LT60	6	Latex	Brown	Thick	0.34 mm	Hydrogel	40 pairs/box
2D72LT65	6½						4 boxes/case
2D72LT70	7						
2D72LT75	7½						
2D72LT80	8						
2D72LT85	8½						
2D72LT90	9						



# SURGICAL SUCTION

## MEDI-VAC™ FLEX ADVANTAGE™ HARDWARE

### Hardware / Roll Stands Page 66

Code	Description	Sales Unit
65652-586A	Universal Rollstand 76 cm <b>CE</b>	1 ea/cs
65652-596A	Regulator Mobile Suction Cart – an accessory will have to be added <b>CE</b>	1 ea/cs
65652-145A	Wall Plate W/Holes <b>CE</b>	60 ea/cs
65652-002A	Two Canister Manifold <b>CE</b> 2797	10 ea/cs
65652-003A	Three Canister Manifold <b>CE</b> 2797	10 ea/cs
65652-004A	Four Canister Manifold <b>CE</b> 2797	10 ea/cs

### Canisters Page 66 **CE** 2797

Code	Description	Sales Unit
65652-611A	Medi-Vac™ 1000 ml Flex Canister with Internal On-Off Valve	10 ea/cs
65652-616A	Medi-Vac™ 1500 ml Flex Canister with Internal On-Off Valve	10 ea/cs
65652-631A	Medi-Vac™ 3000 ml Flex Canister with Internal On-Off Valve	8 ea/cs
65652-916A	Medi-Vac™ 1000 ml Flex Canister with T-Connection	10 ea/cs
65652-915A	Medi-Vac™ 1500 ml Flex Canister with T-Connection	10 ea/cs
65652-914A	Medi-Vac™ 3000 ml Flex Canister with T-Connection	8 ea/cs

### Other Hardware Options Page 66 **CE**

Code	Description	Sales Unit
65652-911	Flex Advantage™ single canister base stand 1000 ml/1500 ml	1 ea/cs
65652-913A	Flex Advantage™ single canister base stand 3000 ml	10 ea/cs
65652-924A	Flex Advantage™ single canister roll stand 1000 ml/1500 ml	3 ea/cs
65652-928A	Flex Advantage™ double canister roll stand 1500 ml/3000 ml	10 ea/cs
65652-686A	Dual Head Rollstand	1 ea/cs

### Special Volume Options Page 66

Code	Description	Sales Unit
65651-120	Medi-Vac™ Guardian™ LVC Suction <b>CE</b> 2797	4 ea/cs
65652-583A	Guardian™ LVC Rollstand – 96 cm <b>CE</b>	1 ea/cs
1-700CE	International Universal Critical Measurement – 300 ml <b>CE</b>	50 ea/cs

## MEDI-VAC™ FLEX ADVANTAGE™ SINGLE USE LINERS

### Single Use Liners Page 67 **CE** 2797

Code	Description	Sales Unit
65651910	Flexible Liner, 1000 cc	50 ea/cs
65651920	Flexible Liner, 1500 cc	50 ea/cs
65651930	Flexible Liner, 3000 cc	50 ea/cs
65651910PG	Flexible Liner, 1000 cc + Gel Vac	45 ea/cs
65651920PG	Flexible Liner, 1500 cc + Gel Vac	45 ea/cs
65651930PG	Flexible Liner, 3000 cc + Gel Vac	45 ea/cs
65651916	Flexible Liner 1000 cc + 6 mm x 1.8 m tubing	50 ea/cs
65651926	Flexible Liner 1500 cc + 6 mm x 1.8 m tubing	50 ea/cs
65651929	Flexible Liner 1500 cc + 6 mm x 2.7 m tubing	50 ea/cs
65651936	Flexible Liner 3000 cc + 6 mm x 1.8 m tubing	50 ea/cs
65651939	Flexible Liner 3000 cc + 6 mm x 2.7 m tubing	50 ea/cs

65653916	Flexible Liner 1000 cc + 6 mm x 1.8 m tubing + Medistop Adult	50 ea/cs
65653926	Flexible Liner 1500 cc + 6 mm x 1.8 m tubing + Medistop Adult	50 ea/cs
65653929	Flexible Liner 1500 cc + 6 mm x 2.7 m tubing + Medistop Adult	50 ea/cs
65654916A	Flexible Liner 1000 cc + 6 mm x 1.8 m tubing + Medistop with tip	50 ea/cs
65654926A	Flexible Liner 1500 cc + 6 mm x 1.8 m tubing + Medistop with tip	50 ea/cs

## MEDI-VAC™ FLEX ADVANTAGE™ ACCESSORIES

### Wall Suction Tubing Page 67 **CE**

Code	Description	Sales Unit
516500CE	Medi-Vac™ Universal Vacuum Connection Tubing, 8 mm I.D. x 30.5 m	4 ea/cs
<b>TANDEM</b>		
2-620CE	Flex Advantage™ Tandem Tube	30 ea/cs

### Specimen Socks Page 67 **CE**

Code	Description	Sales Unit
2-122CE	Standard Specimen Sock	10 ea/cs
2-123CE	Short Specimen Sock	10 ea/cs
2-124CE	Large-port Specimen Sock	10 ea/cs
2-169CE	Female Specimen Sock	10 ea/cs

### Vacuum Control Page 68 **CE**

Code	Description	Sales Unit
65655853	Medistop Vacuum Suction Regulator	200 ea/cs
65655854A	Medistop with Tip Vacuum Suction Regulator	200 ea/cs
NS366	Medistop Vacuum Suction Regulator Adult + 1.8 m tubing	30 ea/cs
NS369	Medistop Vacuum Suction Regulator Adult + 2.7 m tubing	30 ea/cs
NS466A	Medistop with Tip Vacuum Suction Regulator + 1.8 m tubing	30 ea/cs

### Gelling Agent Page 68

Code	Description	Sales Unit
GELVAC15	Medivac™ Gelling Agent – 15 g per sachet	1 CS containing 2 boxes of 100 sachets
SOLID500E	MediVac™ Solidifies up to 500 cc Bottle	96 ea/cs
SOLID1200E	MediVac™ Solidifies up to 1200 cc Bottle	96 ea/cs
SOLID1500E	MediVac™ Solidifies up to 1500 cc Bottle	96 ea/cs
SOLID2000E	MediVac™ Solidifies up to 2000 cc Bottle	96 ea/cs
SOLID3000E	MediVac™ Solidifies up to 3000 cc Bottle	96 ea/cs

This range of products are not medical devices and they are not CE marked.

### Accessories Page 68 **CE**

Code	Description	Sales Unit
65652-109	Right angle connector For Flex Advantage™ patient port and CRD™/Guardian™	50 ea/cs
65652-119	Ortho port converter Accepts up to 0.95 cm I.D. tubing with male port	10 ea/cs

## MEDI-VAC™ GUARDIAN™

### Guardian™ Suction Canisters Page 69 **CE** 2797

Code	Description	Sales Unit
65651-212	Rigid canister with mechanical filter – 1200 ml – 1200 cc	40/cs

65651-220	Rigid canister with mechanical filter – 2000 ml – 2000 cc	40/cs
65651-230	Rigid canister with mechanical filter – 3000 ml – 3000 cc	40/cs
65651-312	Guardian™ rigid canister with mechanical filter and tubing 1200 ml + 5 mm x 1.8 m	30 ea/cs
65651-896	Guardian™ rigid canister with mechanical filter and tubing 1200 ml + 6 mm x 1.8 m	30 ea/cs
65651-320	Guardian™ rigid canister with mechanical filter and tubing 2000 ml + 5 mm x 1.8 m	30 ea/cs
65651-395	Guardian™ rigid canister with mechanical filter and tubing 1200 ml + 5 mm x 1.8 m, 5 mm x 0.5 m	30 ea/cs
1-319CE	Guardian™ suction canister with tubing 1200 ml + mechanical inner diameter 6.4 mm x 2.7 m	30 ea/cs

**Guardian™ Ring Bracket with Built-in On/Off Valve** Page 69

Code	Description	Sales Unit
65652-159A	Guardian™ Ring Bracket with Internal On-Off Valve, 1200 ml <sup>CE 2797</sup>	10 ea/cs
65652-139	Guardian™ Ring Bracket with Internal On-Off Valve, 2000 ml <sup>CE 2797</sup>	1 ea/cs
2-542ACE	CRD™ and Guardian™ Tandem Tubes 46 cm x 8 mm I.D. <sup>CE</sup>	90 ea/cs

**ARGYLE™ SURGICAL SUCTION YANKAUERS**

**Argyle™ Flexible Yankauers** Page 70 <sup>CE 0123</sup>

Code	Description	Vacuum Control	Fr/Ch	Length	Sales Unit (Box/package)
<b>ARGYLE™ MINI YANKAUER DEHP FREE 16 CM LENGTH</b>					
1180113001	No side holes	No	16	16 cm	100
1180113050	4 side holes	Yes	16	16 cm	100
1180113100	No side holes	Yes	16	16 cm	100
<b>ARGYLE™ FLEXIBLE YANKAUER DEHP FREE 25 CM LENGTH</b>					
8888504001	Fine tip	No	12	25 cm	50
8888504019	Fine tip	Yes	12	25 cm	50
8888501007	Medium tip	No	18	25 cm	50
8888501023	Medium tip	Yes	18	25 cm	50
8888501015	Medium tip with sump	No	18	25 cm	50
8888501031	Medium tip with sump	Yes	18	25 cm	50
1180501205	Medium tip with sure-grip	No	18	25 cm	50
8888502005	Large tip	No	22	25 cm	50
8888502013	Large tip with sump	Yes	22	25 cm	50
1180502039	Large tip with sure-grip	No	22	25 cm	50
1180504018	Flat tip	No	22	25 cm	50
1180503003	Extra large tip	No	30	25 cm	50
<b>ARGYLE™ FLEXIBLE YANKAUER SETS WITH NON-CONDUCTIVE FLEXIBLE TUBING</b>					
1180504613	Fine tip 5 mm dia tube	No	12	2.2 m**	25
1180504639	Fine tip 5 mm dia tube	No	12	3.0 m**	25
1180504621	Fine tip 5 mm dia tube	No	12	4.0 m**	25
1180501619	Medium tip 7 mm dia tube	No	18	2.2 m**	25
1180501635	Medium tip 7 mm dia tube	No	18	3.0 m**	25
1180501700	Medium tip 7 mm dia tube	Yes	18	3.0 m**	25
1180502617	Large tip 7 mm dia tube	No	22	2.2 m**	25
1180502625	Large tip 7 mm dia tube	No	22	3.0 m**	25
1180502618	Large tip 7 mm dia tube	No	22	2.2 m**	25
1180502658	Large tip 7 mm dia tube	No	22	3.0 m**	25
8888501502	Medium tip 7 mm dia tube	No	18	1.8 m**	50
8888501528	Medium tip with sump 7 mm dia tube	Yes	18	1.8 m**	50
8888502518	Large tip with sump 7 mm dia tube	Yes	22	1.8 m**	50

**Medivac™ Yankauers** Page 71 <sup>CE 2797</sup>

Code	Description	Sales Unit
KEX60	Flexible Suction handle, Fine capacity	50 ea/cs
KEX61	Flexible Suction handle, Regular capacity	50 ea/cs
KEX61V	Flexible suction handle, regular capacity, vacuum control	50 ea/cs
KEX80	Yankauer, tapered bulbous tip – Ch. 26/29 cm	50 ea/cs
KEX82	Yankauer, tapered bulbous tip, vacuum control – Ch. 26/29 cm	50 ea/cs
KEX83A	Yankauer, tapered bulbous tip with 6 mm x 1.8 m sterile tubing	50 ea/cs
KEX84	Yankauer, tapered bulbous tip with 6 mm x 2.7 m sterile tubing	25 ea/cs
KEX85	Yankauer, tapered bulbous tip with 6 mm x 3.7 m sterile tubing	25 ea/cs

**Argyle™ Rigid Yankauers** Page 71 <sup>CE 0123</sup>

Code	Description	Vacuum Control	Fr/Ch	Length	Sales Unit (Box/package)
8888505115	Large tip	No	22	25 cm	50
8888505123	Large tip	Yes	22	25 cm	50
8888505016	Bulb tip	No	34	26 cm	50
8888505024	Bulb tip	Yes	34	26 cm	50

**Argyle™ Selec-Trol Yankauers** Page 71 <sup>CE 0123</sup>

Code	Description	Vacuum Control	Fr/Ch	Length	Sales Unit (Box/package)
1180505172	Large tip	On-off switch	24	27 cm	20
1180505164	Bulb tip	On-off switch	32	27 cm	20

**Argyle™ Poole Suction Instruments** Page 71 <sup>CE 0123</sup>

Code	Description	Fr/Ch	Length	Sales Unit (Box/package)
<b>ARGYLE™ POOLE INSTRUMENTS</b>				
8888508002	Flexible poole	34	27 cm	50
8888509703	Rigid poole	14-28	25 cm	50
<b>ARGYLE™ POOLE SETS</b>				
8888509729	Rigid poole with 5 mm dia cond tube	14-28	1.8 m**	50

**ARGYLE™ CONNECTING TUBING**

**Argyle™ Suction Connecting Tubes** Page 72 <sup>CE 0123</sup>

Code	Diameter and Length	Sales Unit (Box/package)
<b>ARGYLE™ CONNECTING TUBE WITH MOLDED CONNECTORS</b>		
8888301507	5 mm x 50 mm	50
8888301515	5 mm x 1.8 m	50
8888301523	5 mm x 3.1 m	50
8888301531	5 mm x 3.7 m	20
8888301606	6 mm x 1.8 m	50
8888301614	6 mm x 3.1 m	50
8888301622	6 mm x 3.7 m	20
8888301705	7 mm x 1.8 m	50
8888301713	7 mm x 3.1 m	50
8888301721	7 mm x 3.7 m	20

## ARGYLE™ CONNECTING TUBING

### Argyle™ Suction Connecting Tubes Page 72 <sup>CE</sup><sub>0123</sub>

Code	Diameter and Length	Sales Unit (Box/packaging)
<b>ARGYLE™ CONNECTING TUBES WITH FUNNEL TO FUNNEL CONNECTORS</b>		
8888284505	5 mm x 50 cm	50
8888284513	5 mm x 1.8 m	50
8888284539	5 mm x 3.1 m	50
8888284547	5 mm x 3.7 m	20
8888284554	5 mm x 6.1 m	10
8888284604	6 mm x 1.8 m	50
8888284612	6 mm x 3.1 m	50
8888284620	6 mm x 3.7 m	20
8888284638	6 mm x 6.1 m	10
8888284703	7 mm x 1.8 m	50
8888284737	7 mm x 3.7 m	20
8888284745	7 mm x 6.1 m	10
<b>ARGYLE™ CONNECTING TUBES WITH FUNNEL TO TAPERED END</b>		
8888284521	5 mm x 1.8 m	50
8888284711	7 mm x 1.8 m	50

## ARGYLE™ SURGICAL SUCTION CONNECTING TUBES

### Argyle™ Bubble Tubing Page 73 <sup>CE</sup>

Code	Diameter and Length	Sales Unit
8888280214	Argyle™ Bubble Tubing 5 mm x 30.5 m	1
8888280206	Argyle™ Bubble tubing 5 mm x 30 m – Non Conductive	1
1180839035	Argyle™ Bubble Tubing 6 mm x 50 m	1
8888280610	Argyle™ Bubble Tubing 7 mm x 30.5 m	1
8888280602	Argyle™ Bubble tubing 7 mm x 30 m – Non Conductive	1

## MEDIVAC™ SUCTION TUBING

### Medivac™ Suction Tubing Page 73 <sup>CE</sup><sub>2797</sub>

Code	Diameter and Length	Sales Unit
NEX510	5 mm x 3.05 m Sterile Tubing	30 ea/cs
NEX56A	5 mm x 1.8 m Sterile Tubing	50 ea/cs
NEX610	6 mm x 3.05 m Sterile Tubing	25 ea/cs
NEX66A	6 mm x 1.8 m Sterile Tubing	45 ea/cs
NEX710	7 mm x 3.05 m Sterile Tubing	20 ea/cs
NEX712	7 mm x 3.60 m Sterile Tubing	20 ea/cs
NEX76A	7 mm x 1.8 m Sterile Tubing	40 ea/cs

### Argyle™ Specimen Trap Page 74 <sup>CE</sup><sub>0123</sub>

Code	Description	Capacity	Sales Unit (Box/packaging)
<b>ARGYLE™ SPECIMEN TRAP</b>			
8884724504	Specimen trap 70cc	70 ml	50
8884724500	Specimen trap 40cc	40 ml	50

## ARGYLE™ LUKENS MUCUS COLLECTION

### Argyle™ Lukens Mucus Collection Page 74 <sup>CE</sup><sub>0123</sub>

Code	Description	Sales Unit (Box/packaging)
8888258608	Argyle™ Lukens Mucus Collection – 14Fr Catheter – 20 ml	50

## ARGYLE™ TUBING CONNECTORS

### Argyle™ Connectors/Adaptors Page 75 <sup>CE</sup><sub>0123</sub>

Code	Description	Adaptable to Tube Diameters	Sales Unit (Box/packaging)
<b>ARGYLE™ ADAPTORS</b>			
8888270207	Vinyl for sterile tubes SIMS type	5-11 mm	100
8888271007	Vinyl for sterile tubes 5-in-1 type	5-11 mm	100
8888271411	Sterile polypropylene 5-in-1 type	8-14 mm	100
8888271502	Sterile nylon 5-in-1 type	8-14 mm	100
8888272013	Sterile polypropylene large 'Y' type	6-13 mm	100
8888270306	Vinyl for clean tubes SIMS type	5-11 mm	100
8888271403	Clean polypropylene 5-in-1 type	8-14 mm	100
8888271510	Clean nylon 5-in-1 type	8-14 mm	100
8888272005	Clean polypropylene large 'Y' type	6-13 mm	50

## THERMOMETRY

### GENIUS™ 3 TYMPANIC THERMOMETER & BASE

#### Genius™ 3 Tympanic Thermometer and Base Page 77 <sup>CE</sup><sub>0123</sub>

Code	Description	Sales Unit (Box/packaging)
303013	Genius™ 3 Tympanic Thermometer and Base	1

#### Genius™ 3 Tympanic Thermometer and Base – Accessories Page 77

Code	Description	Sales Unit (Box/packaging)
303030	Genius™ probe covers <sup>CE</sup>	2112 (22 x 96)
303097	Genius™ checker/calibrator <sup>CE</sup>	1
303058*	Locking wall mount	1
303059*	Cart with locking mount	1
F31950-SP*	Genius™ replacement cord	1
PT00058372*	Genius™ 3 Calibrator Cable	1

\* This range of products are not medical devices and they are not CE marked.

#### Filac™ 3000 Electronic Thermometer Page 77

Code	Description	Sales Unit (Box/packaging)
502000	Filac™ 3000 probe covers <sup>CE</sup>	2000 (100 x 20)
500028*	Blue oral/axillary isolation chamber	1
500026*	Blue oral/axillary probe with 122 cm cord	1
500036*	Red rectal probe with 122 cm cord	1
500027*	Blue oral/axillary probe with 275 cm cord	1
500037*	Red rectal probe with 275 cm cord	1
500099*	Calibration plug	1
500046*	Mounting base for wall	1

\* This range of products are not medical devices and they are not CE marked.

## UROLOGY

## SILICONE FOLEY CATHETERS

**Dover™ 100% Silicone Foley Catheter** Page 99 

Code	Fr/Ch	Length	Sales Unit
<b>2 WAY 1.5 ML BALLOON</b>			
8887603069	6	32 cm	10
<b>2 WAY 3 ML BALLOON</b>			
8887603085	8	32 cm	10
8887603101	10	32 cm	10
<b>2 WAY 5-10 ML BALLOON</b>			
8887605122	12	42 cm	10
8887605148	14	42 cm	10
8887605163	16	42 cm	10
8887605189	18	42 cm	10
8887605205	20	42 cm	10
8887605221	22	42 cm	10
8887605247	24	42 cm	10
8887605262	26	42 cm	10
<b>2 WAY 30 ML BALLOON</b>			
8887630187	18	42 cm	10
8887630203	20	42 cm	10
8887630229	22	42 cm	10
8887630245	24	42 cm	10
<b>3 WAY 30 ML BALLOON</b>			
8887665183	18	42 cm	10
8887665209	20	42 cm	10
8887665225	22	42 cm	10
<b>2 WAY 5-10 ML BALLOON – FEMALE CATHETER</b>			
8887615121	12	25 cm	10
8887615147	14	25 cm	10
8887615162	16	25 cm	10
8887615188	18	25 cm	10

## SILICONE FOLEY CATHETERS

**Dover™ 100% Silicone Foley Catheter – Tiemann** Page 80 

Code	Fr/Ch	Length	Sales Unit
<b>2 WAY 5-10 ML BALLOON</b>			
8887205121	12	44 cm	10
8887205147	14	44 cm	10
8887205162	16	44 cm	10
8887205188	18	44 cm	10
8887205204	20	44 cm	10
8887205220	22	44 cm	10
8887205246	24	44 cm	10

**Dover™ 100% Silicone Foley Catheter with Glyco-Block Pre-filled Syringe**Page 80 

Code	Fr/Ch	Length	Sales Unit
<b>2 WAY 5-10 ML BALLOON</b>			
8887405127	12	42 cm	5
8887405143	14	42 cm	5
8887405168	16	42 cm	5
8887405184	18	42 cm	5
8887405200	20	42 cm	5
8887405226	22	42 cm	5

## SILICONE FOLEY CATHETERS

**Dover™ 100% Silicone Foley Catheter with Sterile Water Pre-filled**Syringe Page 81 

Code	Fr/Ch	Length	Sales Unit
<b>2 WAY 3 ML BALLOON</b>			
8887103084	8	32 cm	5
8887103102	10	32 cm	5
<b>2 WAY 5-10 ML BALLOON</b>			
8887105123	12	42 cm	5
8887105149	14	42 cm	5
8887105164	16	42 cm	5
8887105180	18	42 cm	5
8887105206	20	42 cm	5
<b>2 WAY 5-10 ML BALLOON – FEMALE CATHETER</b>			
8887115122	12	25 cm	5
8887115148	14	25 cm	5
8887115163	16	25 cm	5

**Dover™ 100% Silicone Foley Catheter Community Catheter** Page 81 

Code	Fr/Ch	Length	Sales Unit
<b>2 WAY 3 ML BALLOON</b>			
8887703083	8	32 cm	1
8887703109	10	32 cm	1
<b>2 WAY 5-10 ML BALLOON</b>			
8887705120	12	42 cm	1
8887705146	14	42 cm	1
8887705161	16	42 cm	1
8887705187	18	42 cm	1
8887705203	20	42 cm	1

## COLLECTION DEVICES

**Dover™ Leg Straps** Page 82 

Code	Description	Sales Unit
31377	Self-adhesive adjustable leg straps	1 pair

**Dover™ Bedside Drainage Bags** Page 82

Code	Description	Sales Unit
90010	120 cm inlet tubing, T Tap outlet valve	20

**MonoFlo™ Bedside Drainage Bags** Page 82 

Code	Description	Sales Unit
6208-02	90 cm drainage inlet tube	30
6209-02	120 cm drainage inlet tube	30
6308	120 cm drainage inlet tube, Safeguard™ needle-free sample port and Splashguard™ II™ enclosed drain spout	20

**Dover™ Standard Bedside Drainage Bags** Page 83 

Code	Description	Sales Unit
3512	90 cm drainage inlet tube, 2000 ml	20

**Dover™ TURP Bag** Page 83 

Code	Description	Sales Unit
10082EU	90 cm drainage inlet tube, 4000 ml	25

## URINE METERS

### Dover™ Urine Meter Page 83 <sup>CE</sup><sub>0123</sub>

Code	Description	Sales Unit
2000	150 cm drainage inlet tube, Splashguard™ II™ enclosed drain spout	10
7000	150 cm drainage inlet tube, Safeguard™ Needle-free sample port, Splashguard II™ enclosed drain spout	10
7000LL	150 cm drainage inlet tube, Safeguard™ Needle-free sample port, Splashguard II™ enclosed drain spout, with luer lock sample port	10

### Dover™ Urine Meter Page 83 <sup>CE</sup><sub>0123</sub>

Code	Description	Sales Unit
6500	120 cm drainage inlet tube	20
6502	150 cm drainage inlet tube	20
6500NF	120 cm drainage inlet tube, needle-free sample port	20
6502NF	150 cm drainage inlet tube, needle-free sample port	20

### Dover™ Urine Meter Page 84 <sup>CE</sup><sub>0123</sub>

Code	Description	Sales Unit
900011	120 cm inlet tubing , T Tape outlet valve	1

### Dover™ Urine Meter Page 84

Code	Description	Sales Unit
6600	120 cm drainage inlet tube	10

### Dover™ Urine Meter with Exchangeable Bag Page 84 <sup>CE</sup><sub>0123</sub>

Code	Description	Sales Unit
900012	120 cm inlet tubing, T tape outlet valve	10
900013	2000 ml replacement closed bag	100
900014	2000 ml replacement drainable bag	100

## PRE-CONNECTED SYSTEMS

### Mono-Flo™ Pre-connected Drainage Bags Page 85 <sup>CE</sup><sub>0123</sub>

Code	Description	Sales Unit
6209121	12 Fr/Ch catheter	20
6209141	14 Fr/Ch catheter	20
6209161	16 Fr/Ch catheter	20
6209181	18 Fr/Ch catheter	20

## PRE-CONNECTED SYSTEMS

### Dover™ Pre-connected Urine Meters Page 85 <sup>CE</sup><sub>0123</sub>

Code	Description	Sales Unit
6500121	12 Fr/Ch catheter	20
6500141	14 Fr/Ch catheter	20
6500161	16 Fr/Ch catheter	20
6500181	18 Fr/Ch catheter	20

## WOUND CARE

### ANTIMICROBIAL WOUND DRESSINGS

#### Kendall™ AMD Antimicrobial Foam Dressings Page 87 <sup>CE</sup><sub>0123</sub>

Code	Description	Qty Box	Qty Case
55522AMD	5.1 cm x 5.1 cm	25	100
55544AMD	10.2 cm x 10.2 cm	10	50

55548AMD	10.2 cm x 20.3 cm	10	50
55566AMD	15.2 cm x 15.2 cm	10	50
55588AMD	20.3 cm x 20.3 cm	10	50
55535PAMD	8.8 cm x 7.5 cm Fenestrated, with Topsheet	10	50
55544PAMD	10 cm x 10 cm with Topsheet	10	50

#### Kendall™ AMD Antimicrobial Foam – Bordered Dressing Page 87 <sup>CE</sup><sub>0123</sub>

Code	Description	Qty box	Qty case
55523BAMD	4.4 cm x 8.3 cm, pad size 2.5 cm x 4.4 cm	10	50
55544BAMD	8.9 cm x 8.9 cm, pad size 5.1 cm x 5.1 cm	10	50
55546BAMD	8.9 cm x 14 cm, pad size 5.1 cm x 10.2 cm	10	50
55566BAMD	14 cm x 14 cm, pad size 10.2 cm x 10.2 cm	10	50
55588BAMD	19 cm x 19 cm, pad size 15.2 cm x 15.2 cm	10	50

#### Kendall™ AMD Antimicrobial Foam – Disc Page 87 <sup>CE</sup><sub>0123</sub>

Code	Description	Qty Box	Qty Case
55511AMD	2.54 cm diameter, 4 mm hole	10	40
55512AMD	2.54 cm diameter, 7 mm hole	10	40

## ANTIMICROBIAL WOUND DRESSINGS

#### Telfa™ AMD Antimicrobial Non-Adherent Dressing Page 88 <sup>CE</sup><sub>0123</sub>

Code	Description	Qty Box	Qty Case
<b>TELFA™ AMD PADS STERILE</b>			
7662	7.6 cm x 11.4 cm	50	900
7663	7.6 cm x 20.3 cm	50	600
<b>TELFA™ AMD ISLAND DRESSINGS</b>			
7665	10.2 cm x 12.7 cm (pad size 5 cm x 7.5 cm)	25	200
7666	10.2 cm x 20.3 cm (pad size 5 cm x 15 cm)	25	100
7667	10.2 cm x 25.4 cm (pad size 5 cm x 20 cm)	25	100
7668	10.2 cm x 35.6 cm (pad size 5 cm x 30 cm)	25	50

#### Excilon™ AMD Antimicrobial Dressings Page 88 <sup>CE</sup><sub>0123</sub>

Code	Description	Qty Box	Qty Case
7088	10.2 cm x 10.2 cm pad, drain sponge	25 Pairs	300 Pairs
7089	5.1 cm x 5.1 cm pad I.V. sponge	35 Pairs	700 Pairs

#### Kerlix™ AMD Page 88 <sup>CE</sup><sub>0123</sub>

Code	Description	Qty Box	Qty Case
3331	Roll 11.4 cm x 3.7m	1	60
3332	Roll 11.4 cm x 3.7m soft pack	1	100
6662	Super sponge 15.2 cm x 17.1 cm	2	480
6665	Super sponge 15.2 cm x 17.1 cm	5	600
6660	Super sponge 15.2 cm x 17.1 cm	10	480

## FOAM DRESSINGS

#### Kendall™ Hydrophilic Foam Dressings Page 89 <sup>CE</sup><sub>0123</sub>

Code	Description	Qty Box	Qty Case
<b>KENDALL™ FOAM</b>			
55535	Fenestrated 8.9 cm x 7.6 cm	10	50
55522	5.1 cm x 5.1 cm	25	100
55533	7.6 cm x 7.6 cm	10	50
55544	10.2 cm x 10.2 cm	10	50
55548	10.2 cm x 20.3 cm	10	50
55555	12.7 cm x 12.7 cm	10	50
55566	15.2 cm x 15.2 cm	10	50
55588	20.3 cm x 20.3 cm	10	50



KENDALL™ FOAM PLUS				
Code	Description	Qty Box	Qty Case	
55535P	Fenestrated 8.9 cm x 7.6 cm	10	50	
55522P	5.1 cm x 5.1 cm	25	100	
55533P	7.6 cm x 7.6 cm	10	50	
55544P	10.2 cm x 10.2 cm	10	50	
55548P	10.2 cm x 20.3 cm	10	50	
55566P	15.2 cm x 15.2 cm	10	50	
55588P	20.3 cm x 20.3 cm	10	50	

KENDALL™ FOAM ISLAND				
Code	Description	Qty Box	Qty Case	
55544B	5.1 cm x 5.1 cm pad	10	50	
55566B	10.2 cm x 10.2 cm pad	10	50	
55588B	15.2 cm x 15.2 cm pad	10	50	

### Kendall™ Foam Gentle Border Page 89 CE 0123

Code	Description	Qty Box	Qty Case
55523BG	4.4 cm x 8.3 cm, pad size 2.5 cm x 4.4 cm	10	50
55544BG	8.9 cm x 8.9 cm, pad size 5.1 cm x 5.1 cm	10	50
55546BG	8.9 cm x 14.0 cm, pad size 5.1 cm x 10.2 cm	10	50
55566BG	14.0 cm x 14.0 cm, pad size 10.2 cm x 10.2 cm	10	50
55588BG	19.1 cm x 19.1 cm, pad size 15.2 cm x 15.2 cm	10	50

## ALGINATE DRESSINGS

### Kendall™ Calcium Alginate Dressings Page 90 CE 0123

Code	Description	Qty Box	Qty Case
STERILE 1'S IN PEEL-BACK PACKAGE			
9231	30.5 cm rope	5	20
9232	5.1 cm x 5.1 cm	10	100
9233	10.2 cm x 10.2 cm	10	50
9236	10.2 cm x 10.2 cm plus	10	50
9240	10.2 cm x 14 cm	10	50
9238	10.2 cm x 20.3 cm	5	50
9239	15.2 cm x 25.4 cm	10	50
9242	30.5 cm x 61 cm	5	5
9243	61 cm rope	5	20
9244	91.4 cm rope	5	20

## HYDROGEL DRESSINGS

### Kendall™ Hydrogel Breast Feeding Pad Page 90 CE 0123

Code	Description	Qty Pack	Qty Box	Qty Case
9403	7.6 cm disc	1 pr.	5 pr.	25 pr.

## FILM DRESSINGS

### Kendall™ Transparent Film Dressings Page 90 CE 0123

Code	Description	Qty Box	Qty Case
6651	3.8 cm x 3.8 cm	100	400
6640	5.1 cm x 7 cm	100	400
6641	10.2 cm x 12.1 cm	50	200
6647	10.2 cm x 20.3 cm	20	120
6648	20.3 cm x 25.4 cm	10	12

## NON-ADHERENT DRESSINGS – MOIST

### Xeroform™ Occlusive Gauze Dressing Page 91 CE 0123

Code	Size	Type	Qty Box	Qty Case
8884431302	2.5 cm x 20.3 cm	Strips in overwrap peelable foil	50	200
8884431605	12.7 cm x 22.9 cm	Strips in overwrap peelable foil	50	200
8884432000	10.2 cm x 2.7 m	Roll peelable foil packs	6	36
8884433301	2.5 cm x 20.3 cm	Patch peelable foil packs	50	200
8884433400	5.1 cm x 5.1 cm	Patch peelable foil packs	25	150
8884433500	10.2 cm x 10.2 cm	Patch peelable foil packs	25	150
8884433605	12.7 cm x 22.9 cm	Strip peelable foil packs	50	200

### Xeroform™ Non-Occlusive Gauze Dressing Page 91 CE 0123

Code	Size	Type	Qty Box	Qty Case
8884436400	5.1 cm x 5.1 cm	Patch peelable foil pack	50	300
8884437100	12.7 cm x 22.9 cm	Strips in overwrap peelable foil	12	72
STRIP OVERWRAP FOIL PACK				
8884437201	2.5 cm x 20.3 cm	Strip peelable foil pack	50	200
8884437305	10.2 cm x 2.7 m	Roll	6	36

## NON-ADHERENT DRESSINGS – MOIST

### Vaseline Petrolatum Gauze Page 92 CE 0123

Code	Size	Type	Qty Box	Qty Case
8884412600	2.5 cm x 91.4 cm	Strip in peelable foil packs	12	72
8884413605	7.6 cm x 22.9 cm	Strip in peelable foil packs	50	200
8884415600	7.6 cm x 91.4 cm	Strip in peelable foil packs	12	72
8884416600	15.2 cm x 91.4 cm	Strip in peelable foil packs	12	72
8884422600	2.5 cm x 91.4 cm	Strip in overwrap peelable foil packs	12	72
8884423600	7.6 cm x 22.9 cm	Strip in overwrap peelable foil packs	12	72
8884424600	7.6 cm x 45.7 cm	Strip in overwrap peelable foil packs	12	72
8884425600	7.6 cm x 91.4 cm	Strip in overwrap peelable foil packs	12	72
8884426600	15.2 cm x 91.4 cm	Strip in overwrap peelable foil packs	12	72

## NON-ADHERENT DRESSINGS

### Telfa™ Non-Adherent Absorbent Dressing Page 92 CE 0123

Code	Description	Qty Box	Qty Case
TELFA™ NON-ADHERENT DRESSINGS – STERILE 1'S IN PEEL-BACK PACKAGE			
1055	5.1 cm x 7.6 cm	50	1200
1065	7.6 cm x 10.2 cm dressing	50	1200
1075	7.6 cm x 15.2 cm dressing	50	1200
1085	7.6 cm x 20.3 cm dressing	100	1200

## NON-ADHERENT DRESSINGS

### Telfa™ Non-Adherent Clear Wound Dressings Page 93 CE 0123

Code	Description	Qty Box	Qty Case
STERILE 1'S IN PEEL-BACK PACKAGE			
1109	7.6 cm x 7.6 cm	25	50
1111	10.2 cm x 12.7 cm	25	50
1113	30.5 cm x 30.5 cm	25	50
1114	30.5 cm x 61 cm	25	50

## NON WOVEN DRESSINGS

### Excilon™ Drain Sponges Page 93 CE 0123

Code	Description	Qty Pack	Qty Box	Qty Case
7086	10.2 cm X 10.2 cm	2 pr	25 pr	600 (12 Box/25 pr each)

## GAUZES AND SPONGES

### Curity™ Gauze Sponges Page 93

Code	Description	Qty Tray	Qty Case
1/PACK, STERILE			
6132	7.6 cm x 7.6 cm – 12 ply <small>CE 0123</small>	100	2400
6309	10.2 cm x 10.2 cm – 12 ply <small>CE 0123</small>	100	1200
NO STERILE, BULK			
2346	7.6 cm x 7.6 cm – 12 ply <small>CE</small>	200	4000

## BANDAGES

### Kerlix™ Rolls – Sterile Page 94 <sup>CE</sup><sub>0123</sub>

Code	Description	Qty Case
6725	Bandage roll, medium 8.6 cm x 3.3 m – 6 ply	96
6716	Bandage roll, large 11.4 cm x 2.8 m – 8 Ply	100
7880	Gauze roll, large 11.4 cm x 3.7 m – 8 Ply	100
7884	Kerlix Super Sponge – Medium (15.2 cm x 17.1 cm)	480

### Webriil™ Cotton Undercast Padding Page 94

Code	Description	Qty Bag	Qty Case
1418	5 cm x 3.7 m roll <sup>CE</sup>	24	72
2059	7.6 cm x 3.7 m roll <sup>CE</sup>	12	72
3175	10.1 cm x 3.7 m <sup>CE</sup>	12	72
3489	15.2 cm x 3.7 m roll <sup>CE</sup>	6	36
2283	5.1 cm x 3.7 m roll – sterile 1 per pouch <sup>CE</sup> <sub>0123</sub>		50
2394	7.6 cm x 3.7 m roll – sterile 1 per pouch <sup>CE</sup> <sub>0123</sub>		50
2502	10.1 cm x 3.7 m – sterile 1 per pouch <sup>CE</sup> <sub>0123</sub>		50
2554	15.2 cm x 3.7 m roll – sterile 1 per pouch <sup>CE</sup> <sub>0123</sub>		24

### Webriil™ Cotton Undercast Padding – Crimped Finish Page 94 <sup>CE</sup><sub>0123</sub>

Code	Description	Qty Bag	Qty Case
2666	5.1 cm x 3.7 m roll – sterile 1 per pouch		50
2754	7.6 cm x 3.7 m roll – sterile 1 per pouch		50
2847	10.1 cm x 3.7 m – sterile 1 per pouch		50
2944	15.2 cm x 3.7 m roll – sterile 1 per pouch		24

## ANCILLIARY PRODUCTS

### Webcol™ Skin Cleansing Alcohol Preps Page 95 <sup>CE</sup><sub>0123</sub>

Code	Description	Qty Box	Qty Case
6818-1	Medium 2 Ply	200	4000



## The Essential Medical Solutions Product Guide



This catalogue contains information on products for a wide range of countries and is not limited to products or information that are available and accessible in your specific country. You may contact our customer service or a Cardinal Health representative for more specific information on the products available in your location.

For healthcare professionals only

CARDINAL HEALTH, Cardinal Health LOGO, ESSENTIAL TO CARE, ALTITUDE, ARGYLE, A-V IMPULSE, DEVON, DOVER, EXCILON, FILAC, GENIUS, IMPAD, KANGAROO, KENDALL, KENDALL SCD, KERLIX, MAGELLAN, MEDI-VAC, MONO-FLO, MONOJECT, MEDI-TRACE, SENTINEL SEAL, TELFA, THORA-SEAL, T.E.D., TURKEL, WEBCOL, WEBRIL and XEROFORM are trademarks of Cardinal Health and may be registered in the US and/or in other countries. All other marks are the property of their respective owners.

**Important information:** Before using any medical device, review all relevant package information, including the label and/or the Instructions for Use.

© 2021 Cardinal Health. All Rights Reserved.

2SRV19-1035151-4 (01-2022)